# ImmuGlo™ Anti-Heart Antibody IFA

IVD

REF

1101H

48 Determinations

## INTENDED USE

An indirect immunofluorescence antibody test for the detection and semi-quantitation of anti-heart antibodies in serum of patients with *myocarditis* and *cardiomyopathy*.

## SUMMARY AND EXPLANATION

Anti-heart antibodies have been reported in patients with biopsy proven *myocarditis*, *cardiomyopathy* and to a lesser extent in patients with *ischemic* heart disease. Results in serum of healthy individuals show negative reactions. In addition, other studies suggest an association of heart antibodies to bacterial and viral infections.

Well standardized, indirect immunofluorescence (IF) methods on primate, human or rodent heart tissue provide a simple, non-invasive test to identify individuals who may have immune-mediated cardiac dysfunction<sup>1-6</sup>.

Cardiac specific autoantibodies have been reported in 25% of patients with dilated *cardiomyopathy* at diagnosis<sup>4</sup>. Anti-heart antibodies levels decline over time with disease progression. This suggests that anti-heart antibodies may be predictors of the disease.

## PRINCIPLES OF PROCEDURE

In the indirect immunofluorescence method used in this kit, patients' sera are incubated on optimized preparations of tissue sections to allow binding of antibodies to the substrate. Any antibodies not bound are removed by rinsing the slide. Bound antibodies of the IgG class are detected by incubation of the substrate with fluorescein-labeled, anti-human IgG conjugate. Reactions are observed under a fluorescence microscope equipped with appropriate filters. The presence of anti-heart antibodies is characterized by fibrillar, sacrolemmal and diffuse cytoplasmic reactions.

## PRODUCT INFORMATION

### Storage and preparation

Store all reagents at 2-8°C. Reagents are ready for use after equilibration to room temperature.

### Materials provided

IMMCO Anti-Heart Antibody IFA

REF

1101H

Kits contain sufficient reagents to perform 48 determinations each.

8 x

SORB | SLD | 6

6 well Rat Heart Substrate Slides

1 x 0.5 ml

CONTROL + | CMA \*

**Positive Control.** Contains human serum with BSA.

1 x 0.5 ml

CONTROL - | \*

**Negative Control.** Contains human serum with BSA.

1 x 5 ml

IgG-CONJ | FITC \*

**Anti-human IgG FITC Conjugate.** Contains BSA. **Protect from light.**

- 1 x 60 ml BUF\*
- 2 vials BUF WASH
- 1 x 5.0 ml MOUNTING MEDIUM\*
- 1 x 1.0 ml EVANS
- 1 x 12 COVER SLD

- Sample Diluent.** Contains BSA.
- Phosphate Buffered Saline (PBS).** Dissolve each vial to 1 liter.
- Mounting Medium.** Do not freeze.
- Evan's Blue Counterstain.**
- Coverslips.** Reconstitute to one liter each.

\* Contains < 0.1% NaN<sub>3</sub>

**Material required but not provided**

- Fluorescence microscope
- Micropipette or Pasteur pipette
- Serological pipettes
- Staining dish (e.g. Coplin jar)
- Small test tubes (e.g. 13 x 75 mm) and test tube rack
- Distilled or deionized water
- 1 liter container
- Wash bottle
- Paper towels
- Incubation chamber

**WARNINGS AND PRECAUTIONS**

For *in vitro* Diagnostic Use. All human derived components used have been tested for HbsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. All human serum specimens and human derived products should be treated as potentially hazardous, regardless of their origin. Follow good laboratory practices in storing, dispensing and disposing of these materials<sup>14</sup>.

**WARNING - Sodium azide (NaN<sub>3</sub>)** may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

Instructions should be followed exactly as they appear in this insert to ensure valid results. Do not interchange kit components with those from other sources other than the same catalog number from IMMCO. Do not use beyond expiration date.

**SPECIMEN COLLECTION AND PREPARATION**

Only serum specimens should be used for this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of this test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum should be frozen at -20°C. Avoid repeated freezing and thawing of samples.

**PROCEDURE**

**Test Method**

**A. Screening**

1. Dilute each patient serum 1:10 with the Sample Diluent provided (0.1 ml serum + 0.9 ml Diluent). Do not dilute Positive or Negative Controls. Save the undiluted sera to determine antibody titers if screening tests are found to be positive.
2. Allow pouches containing substrate slides to equilibrate to room temperature for 10-15 minutes. Carefully remove the slides without touching the substrate.
3. Label the slides and place them in an incubation chamber lined with paper towels moistened with water to prevent drying.
4. Invert dropper vial and gently squeeze to apply 1 drop (approximately 50 µl of the Negative Control to well #1. Similarly apply 1 drop of Positive Control to well #2. Avoid overfilling the wells.

5. Using a micropipette or Pasteur Pipette, apply 1 drop of patient's diluted serum (approximately 50 µl) to the other wells. Avoid overfilling the wells.
6. Place the lid on the incubation chamber and incubate slides 30 minutes at room temperature.
7. Remove a slide from the incubation chamber. Hold slide at tab end and rinse gently with approximately 10 ml PBS using a pipette, or rinse slide in beaker filled with PBS. Do not use wash bottle. Transfer slide immediately into Coplin jar and wash 10 minutes. Repeat process with all remaining slides.
8. Remove slide(s) from Coplin jar. Blot the edge of the slide on a paper towel to remove excess PBS. Place the slide in the incubation chamber. Immediately invert the Conjugate dropper vial and gently squeeze to apply 1 drop (approximately 50 µl) to each well.
9. Repeat steps **7 and 8** for each slide.
10. Replace the lid on the incubation chamber. Incubate 30 minutes at room temperature.
11. Remove a slide from incubator. Hold the slide at the tab end and dip the slide in a beaker containing PBS to remove excess conjugate. Place slide(s) in a staining dish filled with PBS for 10 minutes. If desired, 2-3 drops of Evans blue counterstain may be added to the final wash. Repeat for the remaining slides. NOTE: Improper washing may lead to increased background fluorescence.
12. Remove a slide from the staining dish. Blot the edge of the slide on a paper towel to remove excess PBS. **To prevent slide from drying, proceed immediately with next step while slide is still wet.**
13. Mount the coverslip by applying **3 drops** of Mounting Medium evenly on the coverslip and place coverslip over slide. Avoid applying undue pressure and prevent lateral movement of the coverslip.
14. Repeat steps 12 and 13 for each slide.
15. Examine for specific fluorescence under a fluorescence microscope at a magnification of 200x or greater.

Slides may be read as soon as prepared. However, because of the presence of antifading agent in the mounting medium, no significant loss of staining intensity occurs if reading is delayed for up to 48 hours. Slides should be stored in the dark at 2-8°C.

**B. Endpoint Determination (titration)**

A serum positive in the screening test may be further tested following steps 5 through 13 to determine the titer. Each test run should include the Positive and Negative Controls. Make serial two-fold dilutions starting at 1:10. The reciprocal of the highest dilution producing a positive reaction is the titer.

**Preparation of Serial Dilutions**

Number four tubes 1 through 4. Add 0.9 ml of Sample Diluent to tube 1 and 0.2 ml to tubes 2 through 4. Pipette 0.1 ml of undiluted serum to tube 1 and mix thoroughly. Transfer 0.2 ml from tube 1 to tube 2 and mix thoroughly. Continue transferring 0.2 ml from one tube to the next after mixing to yield the dilutions depicted in the following table:

<b>Tubes</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>
<b>Serum</b>	<b>0.1 ml</b>			
	<b>+</b>			
<b>Buffered Diluent</b>	<b>0.9 ml</b>	<b>0.1 ml</b>	<b>0.1 ml</b>	<b>0.1 ml</b>
				
<b>Transfer</b>		<b>0.2 ml</b>	<b>0.2 ml</b>	<b>0.2 ml</b>
<b>Final dilution</b>	<b>1:10</b>	<b>1:20</b>	<b>1:40</b>	<b>1:80 etc.</b>

**QUALITY CONTROL**

Both a Positive and Negative Control should be included with each test run. The Negative Control should show no specific fluorescence, whereas the Positive Control should have

2+ or greater staining intensity.

If expected results are not obtained, the run should be repeated. If inadequate results continue to occur with the controls, these may be due to:

- Turbidity. Discard and use another control
- Problems with the optical system of the fluorescence microscope. These may include : improper alignment, bulb beyond useful life expectancy, etc.
- Allowing the slide to dry during the procedure.

### **INTERPRETATION OF RESULTS**

Test results for anti-heart antibodies should be read negative (<10) or positive with titer. Read for specific staining of the heart tissue. Organ specific heart antibodies give fine striational IF reactions on cardiac tissue but are either negative or weakly positive on skeletal muscle tissue<sup>9</sup>.

The following two main staining patterns have been associated with myocarditis:

**1) fibrillar**

**2) sacrolemmal**

See images at the end of this document.

Please note that, low titer anti-heart antibodies could occur in individuals with no heart disease. Various other tissue antibodies such as anti-nuclear antibodies (ANA), and anti-mitochondrial antibodies (AMA) may also be observed on heart tissue. Sera exhibiting nuclear reactions may be tested on HEp-2 cell and Liver substrates. Any sera giving smooth muscle or mitochondrial staining reactions should be tested on Mouse Kidney/ Stomach substrate.

### **LIMITATIONS OF THE PROCEDURE**

In some cases, sera positive for cardiac antibodies may either be very weak or negative at the initial screening dilution (prozone phenomenon). In such doubtful cases the sera should be screened at higher dilutions and, if positive, antibody titers determined.

Two or more types of autoantibodies in the serum may react with heart substrate causing interference with their detection and correct identification by immunofluorescence. This results either in the failure to detect anti-cardiac antibodies or suppression of their titer, if the interfering antibody has a higher titer.

### **EXPECTED VALUES**

The incidence of anti-heart antibodies appears at the end of this document.

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
## REFERENCES•REFERENCIAS•LITERATUR•RIFERIMENTI


1. Neumann DA, Burek CL, Baughman KL, Rose NR, Herskowitz A. Circulating heart antibodies in patients with myocarditis or cardiomyopathy. *J Am Coll Cardiol* 16:839-846, 1990
2. Maisch B. Autoreactivity to the cardiac myocyte, connective tissue and the extracellular matrix in heart disease and postcardiac injury. *Springer Semin Immunopathol* 11:368-395, 1989.
3. Majeed HA, Yousof AM, Pokorny J, Bicova R, Bahr G, Behbahani K, Rotta J. Human heart sacrolemmal sheath antibodies in children with non-suppurative sequelae of group A streptococcal infections: a follow up study. *Ann Rheum Dis* 50:752-754, 1991.
4. Caforio ALP, Goldman JH, Baig MK et al. Cardiac autoantibodies in dilated cardiomyopathy become undetectable with disease progression. *Heart* 77:62-67, 1997
5. Bartels C, Honig R, Burger G et al. The significance of anticardiolipin antibodies and anti-heart muscle antibodies for the diagnosis of postpericardiotomy syndrome. *Europ Heart Journal* 15:1494-1499, 1994.
6. Rose NR, Neumann DA, Burek CL and Herskowitz A. Autoimmune Heart Diseases. *The Autoimmune Diseases II*; 12:303-316, 1992.
7. Beutner EH, Kumar V, Krasny SA and Chorzelski TP. Defined immunofluorescence in immunodermatology. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds. John Wiley and Sons, New York, 3rd Ed., 3-40, 1987.
8. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health, 1993 [HHS Pub. No. (CDC) 93-8395].
9. Nicholson GC, Dawkins RL, McDonald BL and Wetherall JD. A Classification of Anti-Heart Antibodies: Differentiation between Heart-Specific and Heterophile Antibodies. *Clin Immunol Immunopath* 7:349-363, 1977.
10. Maisch B. Autoreactivity to the Cardiac Myocyte, Connective Tissue and the Extracellular Matrix of Heart Disease and Postcardiac Injury. *Spring Semin Immunopathol* 11:369-395, 1989.

## **Types of anti-heart antibody reactions:**

**1) fibrillar**

**2) sarcolemmal**

### ***Incidence of Anti-Heart Autoantibodies***

<b>Condition</b>	<b>% positive</b>
Myocarditis/perimyocarditis	60-100
Pericarditis	30-100
Heart transplantation	65-100
Dilated cardiomyopathy	35-80
Normals	10-30

Adapted from Maisch et al<sup>10</sup>.

*For technical assistance please contact:*

**IMMCO**<sup>®</sup>  
DIAGNOSTICS

**IMMCO Diagnostics, Inc.**

**60 Pineview Drive**

**Buffalo, NY 14228-2120**

**Telephone: (716) 691-0091**

**Fax: (716) 691-0466**

**Toll Free USA/Canada: 1-800-537-TEST**

**E-Mail: [info@immcodiagnostics.com](mailto:info@immcodiagnostics.com)**

*or your local product distributor*

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante  
Autorizzato/Representante Autorizado/Représentant Autorisé

EMERGO Group, Inc.

Molenstraat 15, 2513 BH, The Hague,

The Netherlands

Tel (+31) 345 8570, Fax (+31) 346 7299

[www.emergogroup.com](http://www.emergogroup.com)