AUTOANTIBODY TEST SYSTEM

IVD PRODUCT INSERT

REF 1102-60 60 Determinations

REF 1102 100 Determinations

REF 1102-120 120 Determinations

REF 1103 200 Determinations

REF 1103-240 240 Determinations

REF 1103-480 480 Determinations

REF 1103-512 512 Determinations

REF 1125 HEp-2/Mouse Kidney COMVI I 100 Determinations

REF 1134 HEp-2/Mouse Kidney/Stomach COMVI II 100 Determinations

INTENDED USE

Indirect immunofluorescence F antibody tests for the detection and quantitation of anti nuclear antibodies ANA **REF** 1102, 1102-120, 1103, 1103-480, 1103-512, 1102-60, 1103-240, 1125, 1134, anti mitochondrial antibodies AMA anti smooth muscle antibodies ASMA **REF** 1125, 1134 and anti gastric parietal cell antibodies AGPA **REF** 1134 in human serum.

SUMMARY AND EXPLANATION

Antinuclear antibodies (ANA), detected by indirect immunofluorescence aid in the diagnosis of connective tissue disorders including systemic lupus erythematosus SLE mixed connective tissue disease Sjögren's syndrome and scleroderma ^{1 5}. ANA occur in about 95% of SLE patients as well as patients with other connective tissue diseases. ANA may also occur in other disorders such as chronic active hepatitis and primary biliary cirrhosis ^{6 8}.

Anti-mitochondrial antibodies (AMA) occur in over 90% of primary biliary cirrhosis cases 3 11% of chronic active hepatitis patients and are absent in patients with extra hepatic biliary obstruction and in other liver diseases. The universal presence of anti mitochondrial antibodies in primary biliary cirrhosis and their virtual absence in extra hepatic jaundice makes their detection of considerable value in the differential diagnosis ^{6 12}.

Anti-smooth muscle antibodies (ASMA) in high titer >160 occur in the majority of cases of chronic active hepatitis and in intermediate titers 40 80 in acute viral hepatitis. Occasionally they may occur in cases of primary biliary cirrhosis where they are also found in intermediate titers. The significance of titers of 20 40 is doubtful since these titers may occur in normal individuals ^{13 14}.

Anti-gastric parietal cell antibodies (AGPA) are commonly associated with pernicious anemia and chronic atrophic gastritis where they occur in about 90% and 50% of cases respectively. However they are not disease specific as they may occur in low frequency in other disorders. Although healthy individuals may have gastric parietal cell antibodies this finding may reflect asymptomatic atrophic gastritis. Negative findings for gastric parietal cell antibodies provide strong evidence for excluding pernicious anemia ^{15 17}.

PRINCIPLES OF PROCEDURE

In the indirect F method used in this kit patients' sera are incubated on a variety of substrates HEp 2 cells or HEp 2 cells and mouse kidney/stomach sections to allow binding of antibodies. Any antibodies not bound are removed by rinsing. Bound antibodies of the gG class are detected by incubation of the substrate with fluorescein labeled anti human gG conjugate. Reactions are observed under a fluorescence microscope equipped with appropriate filters. The presence of ANA ASMA AMA and AGPA is demonstrated by an apple green fluorescence of specific histologic structures in the tissue. The titers the reciprocal of the highest dilution giving a positive reaction are then determined by testing serial dilutions ¹⁸.

PRODUCT INFORMATION

Storage and preparation

Store all reagents at 2-8°C. Reagents are ready for use after equilibration to room temperature.

EN

Materials provided

HEp-2 Cells **REF** 1103-240, 1102-60, 1102, 1102-120, 1103, 1103-480, 1103-512

HEp-2/Mouse Kidney COMVI I **REF** 1125

HEp-2/Mouse Kidney-Stomach COMVI II **REF** 1134

20 x	SORB SLD 6	6 well Substrate Slides (1102-120)
6 x	SORB SLD 10	10 well Substrate Slides (1102-60)
10 x	SORB SLD 10	10 well Substrate Slides (1102, 1125, 1134)
20 x	SORB SLD 10	10 well Substrate Slides (1103)
20 x	SORB SLD 12	12 well Substrate Slides (1103-240)
40 x	SORB SLD 12	12 well Substrate Slides (1103-480)
32 x	SORB SLD 16	16 well Substrate Slides (1103-512)
1 x 0.5 ml	CONTROL + ANA *	ANA Positive Control. Contains human serum.
1 x 0.5 ml	CONTROL + AMA *	AMA Positive Control. Contains human serum. (1125, 1134)
1 x 0.5 ml	CONTROL - *	Negative Control. Contains human serum.
1 x 5 ml	IgG-CONJ FITC *	Anti-human gG F TC Conjugate with or without Evan's Blue. Protect from light. 1 vial (1102-60), 2 vials (1102, 1102-120, 1125, 1134), 3 vials (1103, 1103-240), 6 vials (1103-480, 1103-512)
1 x 60 ml	BUF *	Buffered Diluent. 1 vial (1102-60, 1102, 1102-120, 1125, 1134) 2 vials (1103, 1103-240, 1103-480, 1103-512)
2 x	BUF WASH	Phosphate Buffered Saline (PBS). Dissolve each vial to 1 liter. 2 vials (1102-60, 1102, 1102-120, 1125, 1134) 3 vials (1103, 1103-240, 1103-480, 1103-512)
1 x 5.0 ml	MOUNTING MEDIUM *	Mounting Medium. Do not freeze. 1 vial (1102-60, 1102, 1125, 1134), 2 vials (1102-120, 1103), 3 vials (1103-240), 4 vials (1103-480, 1103-512)
1 x 12	COVER SLD	Coverslips. 1 box (1102-60, 1102), 2 boxes (1102-120, 1103), 3 boxes (1103-240, 1103-480)
1 x 12	COVER SLD LONG	Long Coverslips (1125, 1134), 3 boxes (1103-512)

Optional components

1 x 1.0 ml **EVANS** Evan's Blue Counterstain.

Contains < 0.1% NaN₃

EN

Symbols used on labels:

-  Lot number
-  Catalog number
-  Use by
-  Storage temperature
-  Read instructions for use
-  *n vitro* diagnostic use
-  Manufacturer
-  Number of Tests

Material required but not provided

- Fluorescence microscope
- Micropipette or Pasteur pipette
- Serological pipettes
- Staining dish (e.g. Coplin jar)
- Small test tubes (e.g. 13 x 75 mm) and test tube rack
- Distilled or deionized water
- 1 liter container
- Wash bottle
- Paper towels
- Incubation chamber

WARNINGS AND PRECAUTIONS

For *in vitro* Diagnostic Use. All human derived components used have been tested for HbsAg, HCV, H V-1 and 2 and HTLV- and found negative by FDA required tests. All human serum specimens and human derived products should be treated as potentially hazardous, regard-less of their origin. Follow good laboratory practices in storing, dispensing and disposing of these materials¹⁹.

WARN NG - Sodium azide (NaN₃) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

Instructions should be followed exactly as they appear in this insert to ensure valid results. Do not interchange kit components with those from other sources other than the same catalog number from Immco Diagnostics Inc. Do not use beyond expiration date.

SPECIMEN COLLECTION AND PREPARATION

Only serum specimens should be used for this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of this test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum should be frozen at -20°C. Avoid repeated freezing and thawing of samples.

PROCEDURE

Test Method

A. Screening

1. Dilute each patient serum 1:10 (with the Buffered Diluent provided 20 µl serum + 180 µl diluent) or 1:40 (10 µl serum + 390 µl diluent). Do not dilute Positive or Negative Controls. Save the undiluted sera to determine antibody titers if screening tests are found to be positive.

EN

2. Allow pouches containing substrate slides to equilibrate to room temperature for 10-15 minutes. Carefully remove the slides without touching the substrate.
3. Label the slides and place them in an incubation chamber lined with paper towels moistened with water to prevent drying.
4. Invert dropper vial and gently squeeze to apply 1 drop (approximately 50 μ l) of the Negative Control to well 1. Similarly apply 1 drop of ANA Positive Control to well 2. If applicable, 1 drop of AMA Positive Control to well 3 (1125, 1134). Avoid overfilling the wells.
5. Using a micropipette or Pasteur Pipette, apply 1 drop of patient's diluted serum (approximately 50 μ l) to the other wells. Avoid overfilling the wells.
6. Place the lid on the incubation chamber and incubate slides 30 minutes at room temperature.
7. Remove a slide from the incubation chamber. Hold slide at tab end and rinse gently with approximately 10 ml PBS using a pipette, or rinse slide in beaker filled with PBS. Do not use wash bottle. Transfer slide immediately into Coplin jar and wash 10 minutes. Repeat process with all remaining slides.
8. Remove slide(s) from Coplin jar. Blot the edge of the slide on a paper towel to remove excess PBS. Place the slide in the incubation chamber. Immediately invert the Conjugate dropper vial and gently squeeze to apply 1 drop (approximately 50 μ l) to each well.
9. Repeat steps **7 and 8** for each slide.
10. Replace the lid on the incubation chamber. Incubate 30 minutes at room temperature.
11. Remove a slide from incubator. Hold the slide at the tab end and dip the slide in a beaker containing PBS to remove excess conjugate. Place slide(s) in a staining dish filled with PBS for 10 minutes. If optional conjugate without counterstain is used (see optional components in Materials Provided Section), 2-3 drops of Evan's Blue counterstain may be added to the final wash. Repeat for the remaining slides. NOTE: Improper washing may lead to increased background fluorescence.
12. Remove a slide from the staining dish. Blot the edge of the slide on a paper towel to remove excess PBS. **To prevent slide from drying, proceed immediately with next step while slide is still wet.**
13. Mount the coverslip by applying **3 drops** of Mounting Medium evenly on the coverslip and place coverslip over slide. Avoid applying undue pressure and prevent lateral movement of the coverslip.
14. Repeat steps 12 and 13 for each slide.
15. Examine for specific fluorescence under a fluorescence microscope at a magnification of 200x or greater.

Slides may be read as soon as prepared. However, because of the presence of antifading agent in the mounting medium, no significant loss of staining intensity occurs if reading is delayed for up to 48 hours. Slides should be stored in the dark at 2-8°C.

B. Endpoint Determination (titration)

A serum positive in the screening test may be further tested following steps 5 through 13 to determine the titer. Each test run should include the Positive and Negative Controls. Make serial two-fold dilutions starting at 1:10. The reciprocal of the highest dilution producing a positive reaction is the titer.

Preparation of Serial Dilutions Starting at 1:10

Number six tubes 1 through 6. Add 0.9 ml of Sample Diluent to tube 1 and 0.2 ml to tubes 2 through 6. Pipette 0.1 ml of undiluted serum to tube 1 and mix thoroughly. Transfer 0.2 ml from tube 1 to tube 2 and mix thoroughly. Continue transferring 0.2 ml from one tube to the next after mixing to yield the dilutions depicted in the following table:

Tubes	1	2	3	4	5	6
Serum	0.1 ml					
	+					
Buffered Diluent	0.9 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml
		↗	↗	↗	↗	↗
Transfer		0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml
Final dilution	1:10	1:20	1:40	1:80	1:160	1:320 etc.

QUALITY CONTROL

Both a Positive and Negative Control should be included with each test run. The Negative Control should show no specific fluorescence of the nuclei, smooth muscle, tubules of the kidney or gastric parietal cells. The AMA Positive Control should have 2+ or greater staining intensity of the tubules of the kidney. The ANA Positive Control should have 2+ or greater staining intensity of the nuclei of the kidney with a predominantly homogeneous pattern.

If expected results are not obtained, the run should be repeated. If inadequate results continue to occur with the controls, these may be due to:

- Turbidity. Discard and use another control
- Problems with the optical system of the fluorescence microscope. These may include : improper alignment, bulb beyond useful life expectancy, etc.
- Allowing the slide to dry during the procedure.

INTERPRETATION OF RESULTS

The results of the tests for ANA, AMA, ASMA, and AGPA antibodies should be reported as negative (<10) on kidney stomach sections (1125, 1134) negative (<40) on HEp-2 cells (1102, 1102-60, 1102-120, 1130, 1103-240, 1103-480, 1103-512, 1125, 1134), positive (with titer greater than or equal to 320) or, alternatively, positive with specific endpoint titer.

Read only fields which contain specific staining of the nuclei of the kidney and HEp-2 cells and the pattern observed for ANA, the kidney tubules for AMA, the kidney blood vessel walls for ASMA and gastric parietal cells only for AGPA. All other reactions should be reported as negative for ANA, AMA, ASMA and/or AGPA.

ANA can be detected on all substrates but should be quantified on the kidney or HEp-2 cells. The nuclear staining patterns observable with the kidney substrate or HEp-2 cells provided include homogeneous, peripheral (rim), speckled and nucleolar. The centromere staining pattern (including mitotic figures) is seen most easily on HEp-2 cells. These nuclear staining patterns are described below. They may be one or a combination of several staining patterns. The latter are due to reactions to several different nuclear antigens.

Homogeneous: The entire nucleus fluoresces evenly with a diffuse staining pattern.

Nuclear Membranous: The nuclear membrane stains most intensely as fine linear pattern with decreasing staining intensity of the nucleoplasm towards the center of the nucleus.

Speckled: Discrete coarse to fine round speckles fluoresce throughout the nucleus.

Nucleolar: The nucleoli stain as multiple solid bodies within the nucleus.

Centromere: Large speckles of finite number. Reactive antigen segregates with condensed chromosomes in cells undergoing mitosis.⁵

The specificity of some of the antibodies giving the above staining patterns may be further identified by tests for antibodies to nDNA and to various extractable nuclear antigens. These may be of diagnostic significance as listed in Table 1 at the end of this document.

AMA may be observed on both the distal and proximal tubules of the kidney with the distal tubules staining more brightly. Even though the cytoplasm of the gastric parietal cells also stains, AMA should be quantitated on the kidney.

Staining of the stomach muscularis and kidney glomeruli may also be observed with ASMA, but only ASMA seen on the blood vessel walls of the kidney should be reported.

On HEp-2 cells, detectable cytoplasmic antibodies include anti-mitochondrial antibodies (AMA) and anti-smooth muscle antibodies (ASMA). In an AMA pattern, the cytoplasm appears granular, whereas the ASMA pattern is a

fibrillar network of staining throughout the cytoplasm. Both patterns should be reported as negative for ANA. AGPA only reacts on the parietal cells of the stomach and provides cytoplasmic reactions. Negative reactions on kidney with positive reactions on stomach are indicative of AGPA.

LIMITATION OF THE PROCEDURE

In some cases, sera positive for ANA may either be very weak or negative at the initial screening dilution (prozone phenomenon). In such doubtful cases the sera should be screened at higher dilutions and, if positive, antibody titers determined.

In some cases the presence of two or more antibodies in a serum which are reactive with the same substrate may cause an interference in their detection by immunofluorescence. This interference may cause either failure to detect ANA or suppression of its titer if the interfering antibody has a higher titer than ANA. All ANA reactions should be reported.

The goat anti-human IgG F (TC) Conjugate supplied in this kit is primarily heavy chain specific but has some light chain activity. It reacts primarily with IgG class autoantibodies, but may, to a lesser degree, react with light chains of other classes such as IgM.

A positive ANA should not be considered diagnostic of SLE by itself. They also occur in patients with other connective tissue diseases and certain drugs such as procainamide and hydralazine may induce a positive ANA¹. Moreover, sera of patients with malignancies and infectious diseases may also have positive ANA²⁰.

The clinician should consider the results of all positive indirect immunofluorescence tests along with the results of other laboratory tests and the clinical condition of the patient when making a diagnosis.

EXPECTED VALUES

As seen in Tables 2, 3, 4 and 6 at the end of this document, tests for nuclear antibodies are used to screen for SLE and certain other immunologic disturbances. AMA occur in over 90% of cases of primary biliary cirrhosis and 3-11% of cases of chronic hepatitis. ASMA occur in the majority of cases of chronic active hepatitis and AGPA are commonly associated with pernicious anemia and chronic atrophic gastritis.

PERFORMANCE CHARACTERISTICS

The ImmuGlo™ Autoantibody Test System was compared with another commercially available fluorescent antibody test using HEp-2 cells as a substrate. The comparison included 15 serum samples from normal subjects as well as sera from patients with the diagnosis of SLE, subacute cutaneous lupus erythematosus, scleroderma or rheumatoid arthritis. Sera were tested according to the procedure and screening dilution recommended by the manufacturer. These yielded comparable results as summarized below:

Comparison of Kits Using HEp-2 Cell Substrate for the Detection of Antinuclear Antibodies

Clinical Condition	No. of Sera	% Positive	
		Immco™	Other
SLE	12	100	100
Subacute Cutaneous LE (SCLE)	7	85	85
Scleroderma	6	100	100
Rheumatoid Arthritis	10	50	30
Normal Controls	15	0	0

The ImmuGlo™ Autoantibody Test System (Mouse Kidney/Stomach Sections) was compared with another commercially available fluorescent antibody test using mouse kidney/stomach as a substrate. The comparison included: 20 samples of ANA positive sera, 19 samples of AMA positive sera, 19 samples of ASMA positive sera, 20 samples of AGPA positive sera and 38 serum samples from normal subjects. Sera were tested starting at a 1:10 dilution with the procedure recommended by the manufacturer. These yielded comparable results as summarized in Tables 7 and 8 at the end of this document.

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REFERENCES • ΒΙΒΛΙΟΓΡΑΦΙΑ • LITERATUR • BIBLIOGRAPHIE • BIBLIOGRAFIA

1. Tan EM. Autoantibodies to nuclear antigens (ANA): Their immunobiology and medicine. *Adv Immunol* 33: 167-240, 1982.
2. Kumar V, Beutner EH and Chorzelski TP. Autoimmunity of the skin. In "Concepts in Immunopathology", Vol 1, Cruse JM and Lewis RE Jr, Eds, Karger, Basel, 318-353, 1985.
3. Reimer G, Cornell RC and Tan EM. The biochemical nature of nuclear antigens reactive with antinuclear antibodies. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 519-531, 1987.
4. Beutner EH, Kumar V, Krusny SA and Chorzelski TP. Standardization of antinuclear antibody and other immunofluorescent tests used in immunopathologic studies of the skin. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 41-64, 1987.
5. Tan EM, Chen EKL, Sullivan KF and Rubin RL. Antinuclear antibodies (ANAs): Diagnostically specific immune markers and clues toward the understanding of systemic autoimmunity. *Clin Immunol Immunopathol* 47: 121-141, 1988.
6. Manns M, Gerken G, Kyriakoulis A and Meyer zum Büschenfelde KH. Significant autoimmune markers of autoimmune liver disorders: Current status. *J Clin Lab Anal* 1: 362-370, 1987.
7. Mackay R. Autoimmunity and the liver. *Clin Aspects Immunity* 2: 8-17, 1988.
8. McMillan SA, Alderdice JM, McKee CM et al. Diversity of autoantibodies in patients with anti-mitochondrial antibody and their diagnostic value. *J Clin Pathol* 4: 232-236, 1987.
9. Gershwin ME, Coppel RL and Mackay R. Primary biliary cirrhosis and mitochondrial autoantigens - insights from molecular biology. *Hepatology* 8: 147-151, 1988.
10. Berg PA and Klein R. Mitochondrial antigens and autoantibodies from anti-M1 to anti-M9. *Klin Wochenschr* 64: 897-909, 1986.
11. Popper H and Ponzetto F. Clinic, histologic and immunopathologic features of primary biliary cirrhosis. *Springer Semin Immunopathol* 3: 339-354, 1980.
12. Berg PA and Bacon H. Serology of primary biliary cirrhosis. *Springer Semin Immunopathol* 3: 355-373, 1980.
13. Anderson P, Smith JV and Sobieszek A. Studies on the specificity of smooth muscle antibodies. *Clin Exp Immunol* 22: 22-29, 1975.
14. Kurki P, Miettinen A, Linder E, Pikkariinen P, Vuorio M and Silspuro MP. Different types of smooth muscle antibodies in chronic active hepatitis and primary biliary cirrhosis: Their diagnostic and prognostic significance. *Gut* 21: 878-884, 1980.
15. Fisher JB and Taylor KB. The significance of gastric antibodies. *Brit J Hematol* 20: 1-7, 1971.
16. Chisholm M. Immunology of gastritis. *Clin Gastroenterol* 5: 419-428, 1976.
17. Bigazzi PE, Burek CL and Rose NR. Antibodies to tissue-specific endocrine, gastrointestinal and neurologic antigens. In "Manual of Clinical Laboratory Immunology". Rose NR, Friedman H and Frey JL, Eds, American Society for Microbiology, Washington DC, 762-770, 1986.
18. Beutner EH, Kumar V, Krusny SA and Chorzelski TP. Defined immunofluorescence in immunodermatology. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 3-40, 1987.
19. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health, 1999 [HHS Pub. No. (CDC) 93-8395].
20. Nisengard RJ. Antinuclear antibodies: Significance of titers. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Ben S, Eds, John Wiley and Sons, New York, 2nd Ed, 387-398, 1979.

21. Meyer zum Büschenfelde KH, Manns M and Trutmann F. Autoimmunity in chronic liver diseases - relationship to SLE? In "Recent Advances in Systemic Lupus Erythematosus". Lambert PH, Perrin L, and Zui S, Academic Press, New York, 259-269, 1984.
22. Walker JG, Doniach D, Roitt M and Sherlock S. Serologic tests in the diagnosis of primary biliary cirrhosis. *Lancet* 1: 827, 1965.
23. Pironetto F and Popper H. Hetero-isotopic autoimmune phenomenon in the liver. In "Textbook of Immunopathology", Miescher PA and Müller-Eberhard HJ, Eds, Grune and Stratton, New York, 2nd Ed, 789-817, 1976.
24. Leung PSC, Manns MP, Coppel RL, Gershwin ME. Detection of anti-mitochondrial antibodies in primary biliary cirrhosis and liver-kidney microsomal antibodies in autoimmune hepatitis. In "Manual of Clinical Laboratory Immunology", Rose NR, Hamilton RG, Detrick B, Eds, ASM Press, Washington DC, 6th Ed, 1023-1031, 2002.
25. Muratori P et al. Smooth muscle antibodies and type 1 autoimmune hepatitis. *Autoimmunity*. 35 (8): pp. 497-500. 2002.
26. Gtselis NK et al. Autoantibodies in HCV-treated patients. *World J Gastroenterol*. 11(4):482-487.2005.
27. Miller MH et al. Clinical comparison of cultured human epithelial cells and rat liver as substrates for the fluorescent antinuclear antibody test. *J Rheumatol*. 12 (2): 265-9. 1985.

Table 1. Diagnostic Significance of Antinuclear Antibodies

IF Staining Pattern	Nature of Antigen	Associated Disease
Homogeneous	dsDNA/Histones	SLE
Nuclear membranous	Laminins	SLE, vasculitis or chronic hepatitis
Speckled	RNP	SLE or MCTD*
	Sm	SLE
	SS-A/SS-B	SLE or Sjögren's Syndrome
	ScI-70	Scleroderm
Nucleolar probably U3 RNA	RNAP-I Pm-ScI RNA	Scleroderm
Centromere/Kinetochore	inner and outer plates of kinetochore	CREST syndrome

*Mixed Connective Tissue Disease

Table 2: Incidence of Antinuclear Antibodies (ANA) Detected by Indirect Immunofluorescence on HEp-2 Cells

Clinical Condition	No. of Sera	% Positive
SLE	12	100
Subacute Cutaneous LE (SCLE)	7	86
Scleroderm	6	100
Rheumatoid Arthritis	10	50
Normal Controls	15	0

Table 3: Incidence of Antinuclear Antibodies (ANA) Detected by Indirect Immunofluorescence on Mouse Kidney Substrate

Clinical Condition	No. of Sera	% Positive
SLE	21	95
Scleroderm	17	82
Rheumatoid Arthritis	20	5
Normal Controls	96	0

Table 4: Incidence of Anti-Mitochondrial Antibodies (AMA) Detected by IFA on Mouse Kidney Substrate

Clinical Condition	% Positive
Primary Biliary Cirrhosis	100
Autoimmune Chronic Active Hepatitis	8
HBsAg and Chronic Active Hepatitis	0
Extrahepatic Jaundice and Other Liver Diseases	0
Systemic Lupus Erythematosus (SLE)	3
Rheumatoid Arthritis	0
Normal Controls	0

Adapted from Meyer zum Büschenfelde KH, et al.²¹; Walker JG, et al.²² and Pironetto F and Popper H²³.

Table 5: Incidence of Anti-Smooth Muscle Antibodies (ASMA) as Detected by IFA on Mouse Kidney Substrate

Clinical Condition	% Positive
Chronic Active Hepatitis (Type A)	50-87
Primary Biliary Cirrhosis	25
Acute Viral Hepatitis	87
Infectious Mononucleosis	87
Burkitt's Lymphoma	73
Nasopharyngeal Carcinoma	23
Hodgkin's Disease	23
Myeloproliferative Disorder	5
Warts	4
Normal Controls	3-18

Adapted from Anderson P, et al.¹³

Table 6. Incidence of Anti-Gastric Parietal Cell Antibodies (AGPA) as Detected by Indirect Immunofluorescence on Mouse Stomach Substrate

Clinical Condition	% Positive
Pernicious Anemia (PA)	85-95
Chronic Atrophic Gastritis without PA	30-60
Gastric Ulcer	25-30
Autoimmune Endocrinopathies	25-33
Sjögren's Syndrome	30
First Degree Relatives of PA Patients	30
Normal Controls	
< 20 years old	2
20-60 years old	6-8
> 60 years old	16

Table 7. Findings in Positive Sera

	n	Negatives			Positives Titer		
		< 10	10-20	40-80	160-320	640-2560	
ANA Positive Sera							
Immco™	20	0	0	7	6	7	
Other	20	0	0	7	8	5	
AMA Positive Sera							
Immco™	19	4	1	3	1	10	
Other	19	4	1	3	4	7	
ASMA Positive Sera							
Immco™	19	3	4	8	4	0	
Other	19	2	5	6	5	1	
AGPA Positive Sera							
Immco™	20	0	1	6	4	9	
Other	20	0	2	6	7	5	

Table 8. Findings in Normal Controls

	n	Negatives		Positives Titer	
		< 10	10-20	40-80	>160
ANA Positive Sera					
Immco™	38	30	3	5	0
Other	38	36	0	2	0
AMA Positive Sera					
Immco™	38	38	0	0	0
Other	38	38	0	0	0
ASMA Positive Sera					
Immco™	38	35	2	1	0
Other	38	29	8	1	5
AGPA Positive Sera					
Immco™	38	37	0	1	0
Other	38	37	0	1	0

For technical assistance please contact:

IMMCO Diagnostics, Inc.

60 Pineview Drive

Buffalo, NY 14228-2120

Telephone: (716) 691-0091

Fax: (716) 691-0466

Toll Free USA/Canada: 1-800-537-TEST

E-Mail: info@immco.com

or your local product distributor

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante
Autorizzato/Representante Autorizado/Représentant Autorisé

EMERGO Group, nc.

Molenstraat 15, 2513 BH, The Hague,
The Netherlands

Tel (+31) 345 8570, Fax (+31) 346 7299
www.emergogroup.com