ImmuGlo™

Anti-Neuronal Antibody IFA

IVD

PRODUCT INSERT

REF

Code 1111

48 Determinations

INTENDED USE

An indirect immunofluorescence (IFA) antibody test for the detection and semi-quantitation of anti-neuronal (*paraneoplastic*) autoantibodies; anti-Yo, anti-Hu and anti-Ri in human serum.

SUMMARY AND EXPLANATION

Autoimmune responses of the central nervous system (CNS), recognized as *paraneoplastic neurologic disorders* are manifestations of an antitumor immune response. These consist of a variety of neurological disorders like *Paraneoplastic encephalomyelitis (PE)*, *Paraneoplastic sensory neuropathy (PSN)*, *Paraneoplastic cerebellar degeneration (PCD)*, *Paraneoplastic opsoclonus myoclonus ataxia (POMA)* and *Stiffmann syndrome*¹⁻³. Reliable diagnosis of such conditions and detection of the underlying tumor is difficult. In a significant number of cases in fact, the underlying tumor is not discovered until the patient presents with neurological symptoms^{4,5}. *Paraneoplastic disorders* are characterized by the presence of neuronal autoantibodies in patient serum. The detection of these autoantibodies is useful for the clinician as it represents the presence of an underlying tumor. Tumors that have been known to initiate paraneoplastic disorders are *small-cell lung cancer, neuroblastoma, breast, ovarian and testicular cancers*. The following autoantibodies are found in *paraneoplastic syndromes*⁶⁻⁸.

- a) Anti-Hu, Type I anti-neuronal nuclear antibody (ANNA-1) is associated with small cell lung cancer (SCLC) resulting in PE.
- b) Anti-Yo, anti-purkinje cell antibodies (PCA-1) is associated with ovarian and breast carcinomas resulting in PCD.
- c) Anti-Ri, Type II anti-neuronal nuclear antibody (ANNA-2) is associated with neuroblastoma in children and fallopian or breast cancer resulting in POMA.

The presence of one of these antibodies supports the clinical diagnosis of *paraneoplastic syndrome* and leads to a focused search for underlying neoplasm. These markers also help in discriminating between true paraneoplastic disorders and other inflammatory disorders of the nervous system that mimic a paraneoplastic syndrome.

Indirect immunofluorescence provides a sensitive method of detecting these autoantibodies. Anti-Hu/anti-Ri autoantibodies, which characteristically stain the granular cell nucleus, are easily distinguished from the Purkinje cell cytoplasm staining anti-Yo antibodies. However, to facilitate distinction between Hu and Ri autoantibodies we have designed a triple substrate. This triple substrate presentation also identifies and eliminates anti-nuclear antibodies from the specific screening of the paraneoplastic autoantibodies. If the specimen yields no immunoreactivity by immunofluorescence the result is reported as negative. The specificity of those determined positive by IFA could then be confirmed by Western immunoblotting method.

PRINCIPLES OF PROCEDURE

Antibodies to Hu (ANNA 1), Yo1 (PCA-1), and Ri (ANNA2) are detected by indirect IF using a substrate composite of sections of monkey cerebellum, rodent gut and liver. In this IFA

method, patient sera are incubated on the sections to allow binding of antibodies to the specific antigens of the substrate. Any immunoglobulins and other serum proteins not bound to the tissue sections are removed by rinsing. Bound antibodies of the IgG class are detected by incubation of the substrate with fluorescein-labeled antihuman IgG conjugate. Reactions are observed under a fluorescence microscope equipped with appropriate filters. The presence of antibodies of a specificity is determined by their specific reactions to the cerebellar neurons (Purkinje, molecular or granular cells) and myenteric plexus neurons in the gut.

PRODUCT INFORMATION

Storage and Preparation

Store all reagents at 2-8°C. Reagents are ready for use after equilibration to room temperature.

Materials Provided

ImmuGlo™ anti-Neuronal IFA **REF** 1111

Kits contain sufficient reagents to perform 48 determinations each.

| | | |
|------------|--|---|
| 8 x | SORB SLD 6 | 6 well cerebellum-rodent intestine/ liver substrate slides |
| 1 x 0.5 ml | CONTROL + HU * | Hu (ANNA-1) Positive Control |
| 1 x 0.5 ml | CONTROL - * | Negative Control. Contains human serum. |
| 1 x 5 ml | IgG-CONJ FITC EVANS * | Goat antihuman IgG FITC Conjugate. Ready to use. Protect from light. |
| 1 x 60 ml | BUF * | Sample Diluent. |
| 2 vials | BUF WASH | Phosphate Buffered Saline (PBS). Dissolve each vial to 1 liter. |
| 1 x 5.0 ml | MOUNTING MEDIUM * | Mounting Medium. Do not freeze. |
| 1 x 12 | COVER SLD | Coverslips. Reconstitute to one liter each. |

* Contains < 0.1% NaN₃

Optional Components

| | | |
|------------|-------------------------------|---|
| 1 x 5 ml | IgG-CONJ FITC * | Goat antihuman IgG FITC Conjugate. Ready to use. Protect from light. (2100X) |
| 1 x 1.0 ml | EVANS | Evan's Blue Counterstain. (2510) |

Materials Required But Not Provided

Fluorescence microscope
 Micropipette or Pasteur pipette
 Serological pipettes
 Staining dish (e.g. Coplin jar)
 Small test tubes (e.g. 13 x 75 mm) and test tube rack
 Distilled or deionized water
 1 liter container
 Wash bottle
 Paper towels
 Incubation chamber

WARNINGS AND PRECAUTIONS

For *in vitro* Diagnostic Use. All human derived components used have been tested for HbsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. All human serum specimens and human derived products should be treated as potentially hazardous, regardless of their origin. Follow good laboratory practices in storing, dispensing and disposing of these materials⁹.

WARNING - Sodium azide (NaN₃) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

Instructions should be followed exactly as they appear in this insert to ensure valid results. Do not interchange kit components with those from other sources other than the same catalog number from IMMCO. Do not use beyond expiration date.

SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used for this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of this test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum should be frozen at -20°C. Avoid repeated freezing and thawing of samples.

PROCEDURE

Test Method

A. Screening

1. Dilute each patient serum **1:10** with the Buffered Diluent provided (0.1 ml serum + 0.9 ml Diluent). **Do not** dilute Positive or Negative Controls. Save the undiluted sera to determine antibody titers if screening tests are found to be positive.
2. Allow pouches containing substrate slides to equilibrate to room temperature for **10-15 minutes**. Carefully remove the slides without touching the substrate.
3. Label the slides and place them in an incubation chamber lined with paper towels moistened with water to prevent drying.
4. Invert dropper vial and gently squeeze to apply **1 drop** (approximately 50 µl) of the **Negative Control** to well #1. Similarly apply **1 drop** of **Positive Control** to well #2. Avoid overfilling the wells.
5. Using a micropipette or Pasteur pipette, apply **1 drop** of patient's diluted serum (approximately 50 µl) to the other wells. Avoid overfilling the wells.
6. Place the lid on the incubation chamber and incubate slides **30 minutes** at room temperature.
7. Remove a slide from the incubation chamber. Hold slide at tab end and rinse gently with approximately **10 ml** of PBS using a pipette, or rinse slide in a beaker filled with PBS. Do not use a wash bottle. Transfer slide immediately into a Coplin jar filled with PBS and wash **10 minutes**. Repeat process with all remaining slides.
8. Blot the edge of the slide on a paper towel to remove excess PBS. Place the slide in the incubation chamber. Immediately invert the **Conjugate** dropper vial and gently squeeze to apply **1 drop** (approximately 50 µl) to each well. Repeat process with all remaining slides.
9. Replace the lid on the incubation chamber. Incubate **30 minutes** at room temperature.
10. Remove a slide from the incubator. Hold the slide at the tab end and dip the slide in a beaker containing PBS to remove excess conjugate. Place slide(s) in a staining dish filled with PBS for **10 minutes**. Repeat process with all remaining slides. If desired, **2-3 drops** of Evans blue counterstain may be added to the final wash.

NOTE: Improper washing may lead to increased background fluorescence.

11. Remove a slide from the staining dish. Blot the edge of the slide on a paper towel to remove excess PBS. **While slide is still wet mount the coverslip.** Place **3 drops** of the Mounting Medium evenly spaced on a coverslip placed on a paper towel and invert the slide onto the coverslip. To remove any air bubbles gently apply pressure along the edge of the coverslip. Avoid any movement of the coverslip. Repeat process with all remaining slides
12. Examine for specific fluorescence under a fluorescence microscope at a magnification of **200x** or greater.

Slides may be read as soon as prepared. However, because of the presence of an antifading agent in the mounting medium, no significant loss of staining intensity occurs if reading is delayed. Slides should be stored in the dark at 2- 8°C.

B: End Point Determination (Titration)

A serum positive in the screening test may be further tested following **Steps 4 through 15** to determine the titer. Each test run should include the Positive and Negative Controls. Make serial twofold dilutions starting at 1:10 (see below). Using one slide, a serum may be tested at dilutions ranging from 1:10 to 1:320. If positive at a 1:320 dilution, the titer is reported as greater or equal to 320. Additional slides may be used to obtain endpoints for those sera still positive at a 1:320 dilution. The reciprocal of the highest dilution producing a positive reaction is the titer.

Preparation of Serial Dilutions

Number six tubes 1 through 6. Add 0.9 ml of Buffered Diluent to tube 1 and 0.2 ml to tubes 2 through 6. Pipette 0.1 ml of undiluted serum to tube 1 and mix thoroughly. Transfer 0.2 ml from tube 1 to tube 2 and mix thoroughly. Continue transferring 0.2 ml from one tube to the next after mixing to yield the dilutions depicted in the following table:

| Tubes | 1 | 2 | 3 | 4 | 5 | 6 |
|------------------|--------|--------|--------|--------|--------|------------|
| Serum | 0.1 ml | | | | | |
| | + | | | | | |
| Buffered Diluent | 0.9 ml | 0.2 ml | 0.2 ml | 0.2 ml | 0.2 ml | 0.2 ml |
| | | | | | | |
| Transfer | 0.2 ml | 0.2 ml | 0.2 ml | 0.2 ml | 0.2 ml | 0.2 ml |
| Final dilution | 1:10 | 1:20 | 1:40 | 1:80 | 1:160 | 1:320 etc. |

Quality Control

Both Positive and Negative Controls should be included with each test run. The negative control should show no specific fluorescence of the neurons, whereas the positive control should have 1+ or greater staining intensity.

If expected results are not obtained, the run should be repeated. If inadequate results continue to occur with the controls, these may be due to:

- Turbidity. Discard and use another control.
- Problems with the optical system of the fluorescence microscope. These may include: improper alignment, use of the bulb beyond useful life expectancy, etc.
- Allowing the slide to dry during the procedure.

RESULTS

Results of the test for anti-neuronal antibodies should be reported as negative (<10) or positive with titer. The specificity of the antibody reaction is determined as presented in

Table 1 and Figures 1 and 2 at the end of this document

LIMITATIONS OF THE PROCEDURE

In some cases, sera positive for anti-neuronal antibodies may either be very weak or negative at the initial screening dilution (prozone phenomenon). In such doubtful cases the sera should be retested at higher dilutions and, if positive, antibody titers determined.

The presence of two or more antibodies in a serum reacting with the substrate may cause an interference in their detection by immunofluorescence and may cause either a failure to detect anti-neuronal antibodies or a suppression of the titer. In such cases, tests could be performed by Western Immunoblotting. It is recommended that all positive samples be confirmed for antibody specificity by western immunoblotting techniques.

EXPECTED VALUES

There is a strong association among *paraneoplastic syndromes*, anti-neuronal antibody specificities and the associated tumor type. However, in a minority of patients with paraneoplastic syndrome with anti-neuronal antibodies, no tumor may be found. Certain patients with *paraneoplastic syndromes* may not have detectable levels (by IFA) of anti-neuronal antibodies. Most cases who develop PCD in association with Hodgkin's disease, non-SCLC, gastrointestinal carcinomas, do not have demonstrable anti-Yo antibodies. There also exist anti-Yo negative but atypical anti-Purkinjee cell cytoplasm autoantibodies (atypical APCA's). Such antibodies can occur in association with Hodgkin's disease, *adenocarcinoma* of the lung, colon, or prostate¹⁰. 5-6% patients with ovarian cancer have circulating anti-Yo or anti-Ri antibodies in the absence of any paraneoplastic neurological syndromes¹¹. Table 2 at the end of this document indicates the utility of antibody tests for anti-neuronal antibodies.

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Matériel fourni

ImmUGlo™ Anti-Neuronal Antibody IFA **REF** 1111

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REFERENCES • REFERENCIAS • LITERATUR • RIFERIMENTI


1. Posner, JB. Paraneoplastic syndromes. *Neuro Clinics*; 9:919-936, 1991.
2. Dalmau JO, Posner JB. Paraneoplastic syndromes. *Arch Neurol*; 56: 405-408, 1999.
3. Dropcho EJ. Principles of Paraneoplastic syndromes. *Ann NY Acad Sci*; 841:246-261, 1998.
4. Moll JWB, Henzen-Logmans SC, Splinter TAW et al. Diagnostic value of anti-neuronal antibodies for paraneoplastic disorders of the nervous system. *J Neurol, Neurosurg Psych*; 53:940-943, 1990.
5. Moll JWB, Vecht ChJ. Immune diagnosis of paraneoplastic neurological disease. *Clin Neurol Neurosurg*; 97:71-81, 1995.
6. Lennon VA. Paraneoplastic autoantibodies: The case for a descriptive generic nomenclature. *Neurology*; 44:2236-2240, 1994.
7. Lennon VA. The case for a descriptive generic nomenclature: Clarification of immunostaining criteria for PCA-1, ANNA-1 and ANNA-2 autoantibodies. *Neurology*; 44:2412-2415, 1994.
8. Dalmau J, Posner JB. Neurologic paraneoplastic antibodies (anti-Yo; anti-Hu; anti-Ri): The case for a nomenclature based on antibody and antigen specificity. *Neurology*; 44:2241-2246, 1994.
9. Biosafety in Microbiological and Biomedical Laboratories. Center for Disease Control, National Institute for Health. (HHS Pub. No {CDC} 93-8395) 1993.
10. Dropcho EJ. Autoimmune central nervous system paraneoplastic disorders: Mechanisms, diagnosis, and therapeutic options. *Ann. Neurol*; 37:S102-S113, 1995.
11. Drlicek M, Bianchi G, Boglium G, et. al. Antibodies of the anti-Yo and anti-Ri type in the absence of paraneoplastic neurological syndromes: a long term survey of ovarian cancer patients. *J. Neurol*; 244: 85-89, 1997.
12. Dalmau J, Graus F, Rosenblum MK and Posner JB. Anti-Hu-associated paraneoplastic encephalomyelitis/sensory neuronopathy. A clinical study of 71 patients. *Medicine*; 71:59-72, 1992.
13. Greenlee JE, Brashear HR. Antibodies to cerebellar purkinje cells in patients with paraneoplastic cerebellar degeneration and ovarian carcinoma. *Ann Neurol*; 14:609-613, 1983.
14. Altermatt HJ, Rodriguez M, Scheithauer BW and Lennon VA. Paraneoplastic anti-purkinje and type I anti-neuronal nuclear autoantibodies bind selectively to central, peripheral, and autonomic nervous system cells. *Lab Invest*; 65:412-420, 1991.
15. Luque FA, Furmeaux HM, Ferziger R, et al. Anti-Ri: An antibody associated with paraneoplastic opsoclonus and breast cancer. *Ann Neurol*; 29:241-251, 1991.

Table 1: Immunofluorescent Reaction Pattern Identification


| Antibody Specificity | Cerebellar Purkinje | Neurons Granular | Gut Myenteric Plexus | Liver |
|----------------------|---|---|---|-------|
| Hu (ANNA-1) | Nucleus + Nucleolus - Cytoplasm + | Nucleus + Nucleolus - Cytoplasm + | Nucleus + Nucleolus - Cytoplasm + | - |
| Ri (ANNA-2) | Nucleus - Nucleolus - Cytoplasm - | Nucleus + Nucleolus - Cytoplasm + | Nucleus - Nucleolus - Cytoplasm - | - |
| Yo (PCA-1) | Nucleus - Cytoplasm + | Nucleus - Cytoplasm - | Nucleus - Cytoplasm - | - |

- = Negative; + = Positive

Fig. 1: Hu (ANNA-1) Antibody


Indirect IF on primate cerebellum - staining of nuclei and cytoplasm of Purkinje and granular layer cells.

Indirect IF on rodent gut myenteric plexus

Fig. 2: Yo (PCA-1) Antibody

Indirect IF on primate cerebellum - staining of cytoplasm only.

*Table 2: Clinical Significance of Paraneoplastic anti-Neuronal Antibodies*¹²⁻¹⁵

| Antibody Specificity | Paraneoplastic Neurological Syndrome | Most frequently associated tumors | Western Immunoblot criteria for detection |
|----------------------|---|---|--|
| Anti-Hu (ANNA-1) | Encephalomyelitis, sensory neuropathology, autonomic neuropathy | Small cell lung cancer, neuroblastoma; rarely non small cell lung cancer, prostate cancer, seminoma | 35-40 kD reactive bands on extracts of isolated CNS neurons or with recombinant proteins |
| Anti-Ri (ANNA-2) | Opsoclonus, ataxia, nystagmus, dizziness, dysarthria | Breast cancer, small cell lung cancer | 55 and 80 kD bands on neuronal proteins or with recombinant proteins |
| Anti-Yo (PCA-1) | Subacute cerebellar syndrome, dysarthria and nystagmus | Ovarian cancer, breast cancer, other gynecological cancers | 62 kD band and mostly 24 kD band on purified Purkinje cell extracts or with recombinant proteins |

For technical assistance please contact:

IMMCO Diagnostics, Inc.
60 Pineview Drive
Buffalo, NY 14228-2120
Telephone: (716) 691-0091
Fax: (716) 691-0466
Toll Free USA/Canada: 1-800-537-TEST
E-Mail: info@immcodiagnostics.com

or your local product distributor

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante
Autorizzato/Representante Autorizado/Représentant Autorisé
EMERGO Group, Inc.
Molenstraat 15, 2513 BH, The Hague,
The Netherlands
Tel (+31) 345 8570, Fax (+31) 346 7299
www.emergogroup.com

REV.APR2009

Document No. PI4111 CE