ANTI-NEUTROPHIL CYTOPLASMIC ANTIBODY (ANCA) TEST SYSTEM

IVD

PRODUCT INSERT

REF	1116 ANCA Kit (ethanol) 24 Determinations
REF	1140 ANCA Kit (ethanol) 48 Determinations
REF	1140-2 ANCA Kit (ethanol) 96 Determinations
REF	1140-250 ANCA Kit (ethanol) 250 determinations
REF	1141 ANCA Kit (formalin) 48 determinations
REF	1141-250 ANCA Kit (formalin) 250 determinations
REF	1142 ANCA Kit (ethanol+formalin) 48 determinations

An indirect immunofluorescence antibody test for the detection and semi-quantitation of anti-neutrophil cytoplasmic antibodies (ANCA) in human serum. ANCA are found in the sera of patients with necrotizing vasculitides and hence, serve as an aid to the clinical and other laboratory findings in the diagnosis of these disorders.

SUMMARY AND EXPLANATION

ANCA occur in patients with Wegener's granulomatosis, microscopic polyarteritis, necrotizing or crescentic glomerulonephritis, other vasculitides and inflammatory bowel disorders (primarily ulcerative colitis)¹⁻¹¹. Indirect immunofluorescence staining of ethanol fixed neutrophils may exhibit different types of fluorescent staining patterns. These include:

1. Autoantibodies against cytoplasmic antigens of neutrophils giving a diffuse cytoplasmic staining (cANCA).
2. Autoantibodies against neutrophil antigens giving a perinuclear reaction pattern (pANCA).
3. Common antibodies against nuclear antigens (ANA).

The significance of these two types of ANCA differs¹². cANCA are found primarily in patients with Wegener's granulomatosis and microscopic polyarteritis, whereas pANCA occurs in various vasculitic disorders, ulcerative colitis (UC) and primary sclerosing cholangitis (PSC). Eighty percent of patients in whom pANCA are found have histologic evidence of vasculitis, UC or PSC. ANCA occur in more than 90% of patients with active generalized Wegener's granulomatosis and in 67% of patients with active limited disease. The incidence of ANCA varies in patients with clinical remission. Patients with active Wegener's granulomatosis may occasionally be ANCA negative.

However, in such cases, repeat testing should yield positive ANCA reactions¹. The indirect immunofluorescence microscope method is considered to be the gold standard for detecting ANCA. Antigen fixation of PMN is usually performed in ethanol. When fixed in ethanol, two types of ANCA staining reaction patterns have been identified: cytoplasmic (cANCA) and perinuclear (pANCA). The pANCA reactions can be confirmed by repeat testing on formalin fixed slides where pANCA reactions convert to cANCA whereas ANA reactions either remain nuclear or become negative on formalin fixation.

PRINCIPLES OF PROCEDURE

In the indirect immunofluorescence method used in this kit, patients' sera are incubated on optimized preparations of human neutrophils to allow binding of antibodies to the substrate. Any antibodies not bound are removed by rinsing the slide. Bound antibodies of the IgG class are detected by incubation of the substrate with fluorescein-labeled, anti-human IgG conjugate. Reactions are observed under a fluorescence microscope equipped with appropriate filters. The presence of ANCA is demonstrated by an apple green fluorescence either of the cytoplasm with a diffuse granular cytoplasmic staining (cANCA) or a perinuclear staining (pANCA). The titer (the reciprocal of the highest dilution giving a positive reaction) is then determined by testing serial dilutions¹³.

PRODUCT INFORMATION

Storage and preparation

Store all reagents at 2-8°C. Reagents are ready for use after equilibration to room temperature.

EN

Materials provided

REF 1116	ANCA Kit (ethanol)	24 Determinations
REF 1140	ANCA Kit (ethanol)	48 Determinations
REF 1140-2	ANCA Kit (ethanol)	96 Determinations
REF 1140-250	ANCA Kit (ethanol)	250 Determinations
REF 1141	ANCA Kit (formalin)	48 Determinations
REF 1141-250	ANCA Kit (formalin)	250 Determinations
REF 1142	ANCA Kit (ethanol+formalin)	48 Determinations
4 x	SORB SLD 6	6-well Human Neutrophil Substrate Slides [ethanol fixed] (REF 1116), 8 slides (1140) 16 slides (1140-2)
25 x	SORB SLD 10	10-well Human Neutrophil Substrate Slides [ethanol fixed] (1140-250)
8 x	SORB SLD 6	6-well Human Neutrophil Substrate Slides [formalin fixed], (REF 1141)
25 x	SORB SLD 10	10-well Human Neutrophil Substrate Slides [formalin fixed] (1141-250)
8 x	SORB SLD 6+6	12-well Human Neutrophil Substrate Slides [6 ethanol/6 formalin] (REF 1142)
1 x 0.5 ml	CONTROL + ANCA *	cANCA Positive Control, contains human serum with BSA (REF 1140, 1142, 1116, 1140-2, 1140-250, 1141)
1 x 0.5 ml	CONTROL + pANCA *	pANCA Positive Control, contains human serum with BSA (REF 1141, 1142, 1141-250)
1 x 0.5 ml	CONTROL - *	Negative Control, contains human serum with BSA (REF 1140, 1141, 1142, 1116), 2 vials (1140-2) 3 vials (1140-250, 1141-250)
1 x 5.0 ml	IgG-CONJ FITC *	Anti-human IgG FITC Conjugate with Evan's Blue, contains BSA. Protect from light (REF 1116, 1140, 1141, 1142) 2 vials (1140-2), 3 vials (1140-250, 1141-250)
1 x 60 ml	BUF *	Sample Diluent, contains BSA (1116, 1140-2, 1140, 1141, 1142) 2 vials (1140-250, 1141-250)
2 vials	BUF WASH	Phosphate Buffered Saline (PBS), Dissolve each vial to 1liter. (REF 1116, 1140-2, 1140, 1141, 1142) 3 vials (1140-250, 1141-250)
1 x 5.0 ml	MOUNTING MEDIUM *	Mounting Medium. Do not freeze (REF 1116, 1140, 1140-2, 1141, 1142) 3 vials (1140-250, 1141-250)
1 x 12 ml	COVER SLD	Coverslips (REF 1116, 1140, 1141, 1142) 2x12 (1140-2) 3x12 (1140-250, 1141-250)

* Contains <0.1% NaN₃

Optional materials available from Immco

	REF
6 well Human Neutrophil Substrate Slide (ethanol fixed)	2162
16 well Human Neutrophil Substrate Slide (ethanol fixed)	2162-16
6 well Human Neutrophil Substrate Slide (formalin fixed)	2186
6 + 6 well Ethanol+Formalin Human Neutrophil Substrate Slide (ethanol & formalin fixed)	2189
10 well Human Neutrophil Substrate Slide (ethanol fixed)	2162-10
10 well Human Neutrophil Substrate Slide (formalin fixed)	2186-10
1 x 0.5 ml cANCA Positive Control*, contains human serum with BSA	2252
1 x 0.5 ml pANCA Positive Control*, contains human serum with BSA	2240
1 x 1.0 ml Evan's Blue Counterstain	2510

EN

Symbols used on labels:

-  Lot number
-  Catalog number
-  Use by
-  Storage temperature
-  Read instructions for use
-  In vitro diagnostic use
-  Manufacturer
-  Number of Tests

Material required but not provided

- Fluorescence microscope
- Micropipette or Pasteur pipette
- Serological pipettes
- Staining dish (e.g. Coplin jar)
- Small test tubes (e.g. 13 x 75 mm) and test tube rack
- Distilled or deionized water
- 1 liter container
- Wash bottle
- Paper towels
- Incubation chamber

WARNINGS AND PRECAUTIONS

For *in vitro* Diagnostic Use. All human derived components used have been tested for HbsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. All human serum specimens and human derived products should be treated as potentially hazardous, regardless of their origin. Follow good laboratory practices in storing, dispensing and disposing of these materials¹⁴.

WARNING - Sodium azide (NaN₃) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

Instructions should be followed exactly as they appear in this insert to ensure valid results. Do not interchange kit components with those from other sources other than the same catalog number from Immco Diagnostics Inc. Do not use beyond expiration date.

SPECIMEN COLLECTION AND PREPARATION

Only serum specimens should be used for this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of this test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum should be frozen at -20°C. Avoid repeated freezing and thawing of samples.

PROCEDURE

Test Method

A. Screening

1. Dilute each patient serum 1:20 with the Sample Diluent provided (50 μ l serum + 1.0 ml Diluent). Do not dilute Positive or Negative Controls. Save the undiluted sera to determine antibody titers if screening tests are found to be positive.

EN

2. Allow pouches containing substrate slides to equilibrate to room temperature for 10-15 minutes. Carefully remove the slides without touching the substrate.
3. Label the slides and place them in an incubation chamber lined with paper towels moistened with water to prevent drying.
4. Invert dropper vial and gently squeeze to apply 1 drop (approximately 50 μ l) of the Negative Control to well #1. Similarly apply 1 drop of Positive Control to well #2. Avoid overfilling the wells.
5. Using a micropipette or Pasteur Pipette, apply 1 drop of patient's diluted serum (approximately 50 μ l) to the other wells. Avoid overfilling the wells.
6. Place the lid on the incubation chamber and incubate slides 30 minutes at room temperature.
7. Remove a slide from the incubation chamber. Hold slide at tab end and rinse gently with approximately 10 ml PBS using a pipette, or rinse slide in beaker filled with PBS. Do not use wash bottle. Transfer slide immediately into Coplin jar and wash 10 minutes. Repeat process with all remaining slides.
8. Remove slide(s) from Coplin jar. Blot the edge of the slide on a paper towel to remove excess PBS. Place the slide in the incubation chamber. Immediately invert the Conjugate dropper vial and gently squeeze to apply 1 drop (approximately 50 μ l) to each well.
9. Repeat steps 7 and 8 for each slide.
10. Replace the lid on the incubation chamber. Incubate 30 minutes at room temperature.
11. Remove a slide from incubator. Hold the slide at the tab end and dip the slide in a beaker containing PBS to remove excess conjugate. Place slide(s) in a staining dish filled with PBS for 10 minutes. If desired, 2-3 drops of Evan's Blue counterstain may be added to the final wash. Repeat for the remaining slides. NOTE: Improper washing may lead to increased background fluorescence.
12. Remove a slide from the staining dish. Blot the edge of the slide on a paper towel to remove excess PBS. **To prevent slide from drying, proceed immediately with next step while slide is still wet.**
13. Mount the coverslip by applying **1 drop** of Mounting Medium gently onto each well and place coverslip over slide. Avoid applying undue pressure and prevent lateral movement of the coverslip.
14. Repeat steps 12 and 13 for each slide.
15. Examine for specific fluorescence under a fluorescence microscope at a magnification of 200x or greater.

Slides may be read as soon as prepared. However, because of the presence of antifading agent in the mounting medium, no significant loss of staining intensity occurs if reading is delayed for up to 48 hours. Slides should be stored in the dark at 2-8°C.

B. End point Determination (titration)

A serum positive in the screening test may be further tested following steps 5 through 13 to determine the titer. Each test run should include the Positive and Negative Controls. Make serial two-fold dilutions starting at 1:20. The reciprocal of the highest dilution producing a positive reaction is the titer.

Preparation of Serial Dilutions

Number four tubes 1 through 4. Add 1.0 ml of Sample Diluent to tube 1 and 0.2 ml to tubes 2 through 4. Pipette 0.05 ml of undiluted serum to tube 1 and mix thoroughly. Transfer 0.2 ml from tube 1 to tube 2 and mix thoroughly. Continue transferring 0.2 ml from one tube to the next after mixing to yield the dilutions depicted in the following table:

Tubes	1	2	3	4
Serum	0.05 ml			
	+			
Buffered Diluent	1.0 ml	0.2 ml	0.2 ml	0.2 ml
Transfer		↻ 0.2 ml	↻ 0.2 ml	↻ 0.2 ml
Final Dilution	1:20	1:40	1:80	1:160 etc.

EN

QUALITY CONTROL

Both a Positive and Negative Control should be included with each test run. The Negative Control should show no specific fluorescence of the neutrophil, whereas the Positive Control should have 2+ or greater staining intensity of the cytoplasm of the neutrophil with cANCA positive control and perinuclear with pANCA positive control on ethanol fixed slides. On formalin fixed slides, the cANCA Positive Control remains cytoplasmic whereas pANCA reaction will become cytoplasmic.

If expected results are not obtained, the run should be repeated. If inadequate results continue to occur with the controls, these may be due to:

- Turbidity. Discard and use another control
- Problems with the optical system of the fluorescence microscope. These may include : improper alignment, bulb beyond useful life expectancy, etc.
- Allowing the slide to dry during the procedure.

INTERPRETATION OF RESULTS

The results of the tests for ANCA should be reported as negative (<20), positive greater or equal to 160, or preferably, positive with titer and pattern. Read for specific diffuse, granular cytoplasmic staining (cANCA) or perinuclear staining (pANCA). Other detectable antibodies include antinuclear antibodies (ANA), which sometimes may mimic pANCA reactions. ANCA patterns on ethanol and formalin fixed slides are presented at the end of the document.

LIMITATION OF THE PROCEDURE

In some cases, sera positive for ANCA may either be very weak or negative at the initial screening dilution (prozone phenomenon).

In such doubtful cases the sera should be screened at higher dilutions and, if positive, antibody titers determined.

In some cases the presence of two or more antibodies in a serum which are reactive with the same substrate may cause an interference in their detection by immunofluorescence. This interference may cause either failure to detect ANCA or suppression of its titer if the interfering antibody has a higher titer than ANCA. The most common cause of the interference phenomenon in ANCA tests is the coexistence of ANA.

In some patients with Wegener's granulomatosis, ANCA tests may be negative. In such cases repeat testing may yield positive results. Patients with Wegener's granulomatosis on treatment invariably are negative for ANCA. ANA reactions may sometimes be confused with or mimic pANCA staining. To confirm pANCA reactivity, sera providing pANCA reactions should be retested either on formalin fixed slides or for ANA on HEp-2 slides. pANCA specimens on formalin fixed ANCA substrate should give cANCA reactivity, whereas reactions to ANA should be either negative or remain nuclear. Anti-cytokeratin antibodies may result in false positive cANCA reaction¹⁵.

In such cases, indirect IF tests on HEp-2 cells may help in distinguishing pseudo-ANCA from true cANCA reaction.

The antibody titers do not necessarily associate with disease activity. The results of the ANCA test should be evaluated in light of clinical findings as the presence or absence of ANCA may not be directly associated with a vasculitis disorder.

The positive ANCA results obtained by immunofluorescence should be confirmed by ELISA. ANCA of certain antigen specificities are more indicative of a specific vasculitic disorder. Also, ANCA have been associated with immunological disturbances other than vasculitis disorders such as ulcerative colitis^{7,8}.

EXPECTED VALUES

Sixty-four (64) normal samples were tested for ANCA. All samples were negative for ANCA at 1:20 dilution.

A positive ANCA in the appropriate clinical setting is useful in the diagnosis of systemic vasculitides and inflammatory bowel disorders⁷⁻¹¹. cANCA staining occurs most commonly in Wegener's granulomatosis and

pANCA staining in microscopic polyangitis, pauci-immune crescentic glomerulonephritis and ulcerative colitis⁷⁻¹¹. In other vasculitides (polyarteritis nodosa, takayasu disease, giant cell arteritis, Behcet's disease) ANCA are rare or absent. The incidence of ANCA in Wegener's granulomatosis and other vasculitides as abstracted from literature is summarized in Table 1, at the end of the document. In the European ANCA collaborative assay standardization project, Hagan and his collaborators¹⁶ evaluated the usefulness of the indirect IF test towards the diagnosis of idiopathic system vasculitides. In this study, ANCA was found to be present in 85% of patients with Wegener's granulomatosis, of these, 64% were positive for cANCA and 21% for pANCA. The pANCA was more prevalent than cANCA in microscoping polyangitis (58% vs 23%) with sensitivity of 81% for microscoping polyangitis and 82% for idiopathic rapidly progressive glomerulonephritis.

Of the disease and healthy controls, ANCA was present in 19% and 6% respectively, thus providing a specificity of 76% in disease controls and 94% in healthy individuals. (Table 2, at the end of the document).

However, the positive predictive value of the ANCA tests is much better and more significant if evaluated in conjunction with clinical signs and symptoms. In an editorial, Jennette Wilmant and Falk¹⁷ reported a positive predictive value of 92% in a patient with serum creatinine >3 mg/dl. Similarly, in patients in which the initial positive predictive value is low, a positive ANCA result increases the likelihood of a disease to a level that may warrant further evaluation. The authors conclude that ANCA testing in patients with strong clinical evidence for pauci-immune crescentic glomerulonephritis is most useful for substantiating the diagnosis whereas ANCA testing in patients with weak clinical evidence is most useful for ruling out the diagnosis.

SPECIFIC PERFORMANCE CHARACTERISTICS

The Anti-Neutrophil Cytoplasmic Antibody (ANCA) Test was compared with another indirect immunofluorescence ANCA test. The comparison included a total of 129 serum samples obtained from a diagnostic reference laboratory specializing in the detection of autoimmune diseases. These samples were tested according to the procedures recommended by the manufacturers. The results are as follows:

		Immco™		Total
		Positive	Negative	
Other	Positive	49	0	49
	Negative	7	73	80
	Total	56	73	129

Relative Specificity: 91%
Relative Sensitivity: 100%
Relative Agreement: 95%

Seven discrepant samples observed by indirect immunofluorescence method above were tested by ELISA for antibodies to myeloperoxidase (MPO) or proteinase 3 (PR3) antigens, the two major antigens associated with ANCA. Of the seven positives on ELISA and negative on the other test kit, all but one tested positive by ELISA, suggesting thereby, the true positivity of these samples.

Cross Reactivity:

Antinuclear antibodies (ANA) may exhibit positive reactions on neutrophils. To determine reactivities on neutrophils of ANA, which may be confused for ANCA reactivity, we tested a total of 25 ANA positive samples of varying antibody specificities on ethanol, formalin and COMVI slides. All ANA positive sera with the exception of SS-A (Ro) and SS-B (La) antibody specificities, showed nuclear reactions. These ANA reactions either remain nuclear or become negative on formalin fixed slides.

Reproducibility

Studies were performed to demonstrate intra and inter-assay variability. Four ANCA positive (two each of c and p ANCA) and one ANCA negative sera were tested starting at a 1:20 dilution to endpoint. They were tested on 3 different lots each of ethanol, formalin and COMVI slides for four days to determine intra as well as inter-reproducibility. Negative samples remained negative and positive samples provided the expected titer.

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

REFERENCES • ΒΙΒΛΙΟΓΡΑΦΙΑ • LITERATUR • BIBLIOGRAPHIE • BIBLIOGRAFIA

1. Nölle B, Specks U, Lüderman J, Rohrbach M, DeRemee RA and Gross WL. Anticytoplasmic antibodies: Their immunodiagnostic value in Wegener's granulomatosis. *Ann Int Med* 111: 28-40, 1989.
2. Venning MC, Quinn A, Broomhead V and Bird AG. Antibodies directed against neutrophils (cANCA and pANCA) are of distinct diagnostic value in systemic vasculitis. *Quart J Med* 77: 1287-1 296, 1990.
3. Van der Woude FJ, Daha MR and Van Es LA. The current status of neutrophil cytoplasmic antibodies. *Clin Exp Immunol* 78: 143-1 48, 1989.
4. Specks U, Wheatley CL, McDonald TJ, Rohbrach MS and DeRemee RA. Anticytoplasmic autoantibodies in the diagnosis and follow up of Wegener's granulomatosis. *Mayo Clin Proc* 64: 28-36, 1989.
5. Tervert JWC, van der Woude FJ, Fauci AS and Ambrus JL. Association between active Wegener's granulomatosis and anticytoplasmic antibodies. *Arch Int Med* 149: 2461-2465, 1989.
6. Cross CE and Lillington GA. Serodiagnosis of Wegener's granulomatosis: Pathobiologic and clinical implications. *Mayo Clin Proc* 64: 119-1 22, 1989.
7. Seibold F, Slametschka D, Gregor X and Weber P. Neutrophil Autoantibodies: A genetic marker in primary sclerosing cholangitis and ulcerative colitis. *Gastroenterol* 107:532-536, 1994.
8. Claise C, Johanet C, Bouhnik Y et al. Antineutrophil cytoplasmic autoantibodies in autoimmune liver and inflammatory bowel diseases. *Liver* 16:28-34, 1996.
9. Gigase P, DeClerck LS, Van Cotthem KA et al. Anti-Neutrophil cytoplasmic antibodies in inflammatory bowel disease with special attention for IgA-class antibodies. *Dig Dis and Sci* 42:2171 -2174, 1997.
10. Shanahan F. Neutrophil autoantibodies in inflammatory bowel disease: Are they important? *Gastroenterol* 107:586-589, 1994.
11. Shanahan F and Bernstein CN. ANCAs aweigh in colitis. *Gastroenterol* 105:946-947, 1993.
12. Lüdemann J, Utecht B and Gross WL. Laboratory methods for detection of antineutrophil cytoplasm antibodies. *Clin Immunol Newsletter* 10:159-166, 1990.
13. Beutner EH, Kumar V, Krasny SA and Chorzelski TP. Defined immunofluorescence in immunodermatology. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 3-40, 1987.
14. Biosafety in Microbiological and Biomedical Laboratories. Center for Disease Control, National Institute for Health, HHS Pub. No CDC 93-8395) 1999.
15. Streicher J, Fabian B, Herkner K et al. Anti-cytokeratins are a potential source of false positive indirect immunofluorescence assays for cANCA. *J Clin Lab Analysis*. 12:54-59, 1998.
16. Hagen CF, Daha MR, Hermand J et al. Diagnostic value of standardized assays for anti-neutrophil cytoplasmic antibodies in idiopathic systemic vasculitis. *Kidney Int* 53:743-753, 1998.
17. Jennette JC, Wilkman AS, Falk RJ. Diagnostic predictive value of ANCA serology. *Kidney Int* 53:796-798, 1998.

ANCA Reactions on Ethanol Fixed Slides

ANCA Reactions on Formalin Fixed Slides



Table 1: Prevalence of ANCA¹

Clinical Condition	Positive
Wegener's granulomatosis	
Generalized	
Active	96
Partial remission	71
Full remission	41
Local recurrence	80
Localized	
Active	67
Partial remission	54
Full remission	32
Inflammatory Bowel Disorders	
Ulcerative colitis	70
Primary sclerosing cholangitis	82
Crohn's disease*	27
Disease Controls	
Blood Donors	0
Connective tissue	
Autoimmune disorders**	5
Misc. medical conditions	0
Granulomatosis disease	0
Primary renal disease	1

* ANCA titers are usually low ** ANCA Reactivity is pANCA

Table 2: Sensivity and Specificity of the IF Test in Patients with Systemic Vasculitides (Adapted from Reference 16)

	N	Sensitivity		
		cANCA	pANCA	c or pANCA
Patients				
Wegener's granulomatosis	97	64	21	85
Microscopic polyangitis	44	23	58	81
Idiopathic RPGN	12	36	45	81
Classical polyarteritis nodosa	10	10	30	40
Churg-Strauss syndrome	6	33	33	66
Controls				
Specificity				
Disease Controls	184	95	81	76
Healthy controls	740	98	96	94

For technical assistance please contact:

IMMCO Diagnostics, Inc.

60 Pineview Drive

Buffalo, NY 14228-2120

Telephone: (716) 691-0091

Fax: (716) 691-0466

Toll Free USA/Canada: 1-800-537-TEST

E-Mail: info@immco.com

or your local product distributor

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante
Autorizzato/Representante Autorizado/Représentant Autorisé

EMERGO Group, Inc.

Molenstraat 15, 2513 BH, The Hague,
The Netherlands

Tel (+31) 345 8570, Fax (+31) 346 7299
www.emergogroup.com