Anti-Gliadin Antibody (AGA) ELISA

IVD

CLIA Complexity: High
 CDC Analyte Identification Code: 0528
 CDC Test System Identification Code: 28513/28514

PRODUCT INSERT

REF 1117A IgA-AGA ELISA 96 Determinations

REF 1117G IgG-AGA ELISA 96 Determinations

An enzyme linked immunosorbent assay (ELISA) for the detection and semi-quantitation of IgG- or IgA-anti-gliadin antibodies in human serum to aid in the diagnosis of patients with *celiac disease* and *dermatitis herpetiformis*.

SUMMARY AND EXPLANATION

Gluten sensitive enteropathy (GSE), i.e. *celiac disease* (CD) and *dermatitis herpetiformis* (DH), is a common clinically heterogeneous gastrointestinal disorder with sensitivity to gluten, and can exhibit with non-classic or minimal symptoms¹. There is some genetic component associated with GSE which is illustrated by showing, that approximately 5-10% of first degree relatives have symptomatic or asymptomatic CD. Children with short stature and patients with insulin dependent diabetes and other autoimmune disorders also demonstrate a greater likelihood of developing GSE. Strict avoidance of gluten in the diet is recommended to control the disease activity and early diagnosis in such patients may improve their overall prognosis².

Published literature suggests the use of serological testing for patients suspected of GSE and to monitor dietary compliance³⁻⁶. The European Society of Pediatric Gastroenterology and Nutrition (ESPGAN) has recommended the inclusion of various serological tests to reduce the number of intestinal biopsies needed to make a diagnosis⁷. These include tests for anti-gliadin (AGA), anti-reticulin (ARA) and anti-endomysial (EMA) antibodies. Of these, AGA have been studied the most⁸⁻²⁰. Using ELISA, both IgA and IgG class AGA are detected in sera of patients with GSE. Of these, IgG class AGA seem more sensitive but less specific indicators of GSE than IgA class AGA. IgA class AGA on the other hand are less sensitive but more specific for GSE.

PRINCIPALS OF THE PROCEDURE

Gliadin antigen is bound to the wells of a polystyrene microwell plate followed by blocking the unreacted sites to reduce non-specific binding. Controls, calibrators and diluted patient sera are added to separate wells, allowing any gliadin antibodies present to bind to the immobilized antigen. Unbound sample is washed away and an enzyme labeled anti-human IgG or IgA conjugate is added to each well. These enzyme conjugated antibodies bind specifically to the human immunoglobulin of the appropriate class. After washing away any unbound conjugate, specific enzyme substrate (pNPP) is then added to the wells. After stopping the enzymatic reaction, the intensity of color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 405 nm. Results are expressed in ELISA units per milliliter (EU/ml).

REAGENTS

Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.** Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date. Reconstitute the wash buffer to 1 liter with distilled or deionized water. Coated microwell strips are for one time use only.

Precautions

For *in vitro* Diagnostic Use. All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials²¹. **WARNING** - Sodium azide (NaN₃) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup.

EN

Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

Instructions should be followed exactly as they appear in this kit insert to en-sure valid results. Do not interchange kit components with those from other sources other than the same catalog number from Immco Diagnostics Inc.

Use good laboratory techniques to minimize microbial and chemical contamination. Do not use after expiration date.

Materials provided

ImmuLisa™ IgA-AGA ELISA **REF** 1117A

ImmuLisa™ IgG-AGA ELISA **REF** 1117G

12 x 8	MICROPLATE AGA	Microplate with individual breakaway microwells coated with gliadin antigen.
1 x 1.5 ml	CALIBRATOR A AGA * †	Ready to use Calibrator A (<i>green cap</i>). Human serum containing antibodies to gliadin.
1 x 1.5 ml	CALIBRATOR B AGA * †	Ready to use Calibrator B (<i>violet cap</i>). Human serum containing antibodies to gliadin.
1 x 1.5 ml	CALIBRATOR C AGA * †	Ready to use Calibrator C (<i>blue cap</i>). Human serum containing antibodies to gliadin.
1 x 1.5 ml	CALIBRATOR D AGA * †	Ready to use Calibrator D (<i>yellow cap</i>). Human serum containing antibodies to gliadin.
1 x 1.5 ml	CONTROL + AGA * †	Ready to use Positive Control (<i>red cap</i>). Contains human serum positive for AGA.
1 x 1.5 ml	CONTROL - * †	Ready to use Negative Control (<i>white cap</i>). Contains human serum.
1 x 12 ml	IgA-CONJ ALKPHOS * †	Ready to use anti-human Alk. Phos. Conjugate . Color coded pink.
1 x 12 ml	IgG-CONJ ALKPHOS * †	Ready to use anti-human Alk. Phos. Conjugate . Color coded pink.
1 x 60 ml	DIL *	Ready to use Serum Diluent . Color coded blue.
1 x 12 ml	SUBSTRATE *	Ready to use Enzyme Substrate . Contains pNPP. Protect from light.
1 x 12 ml	STOP	Ready to use Stop Solution .
2 x	BUF WASH	Powder Wash Buffer . Reconstitute to one liter each.

* Contains < 0.1% NaN₃

† **REF** 1117A contains IgA-AGA calibrators, controls and IgA conjugate

REF 1117G contains IgG-AGA calibrators, controls and IgG conjugate

Symbols used on labels:

-  Lot number
-  Catalog number
-  Use by
-  Storage temperature
-  Read instructions for use
-  In vitro diagnostic use
-  Manufacturer
-  Number of Tests

Materials Required But Not Provided

- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Deionized or distilled water
- Microplate reader capable of reading absorbance values at 405 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm.
- Squeeze bottle to hold diluted wash buffer
- Timer
- Absorbent paper
- Automated microplate washer capable of dispensing 200 µl

SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2° - 8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples.

PROCEDURE**Procedural Notes**


- Before starting the assay read these instructions carefully.
- Bring all reagents and samples to room temperature (20-26°C) for 30 minutes. Return materials to refrigerator immediately after use.
- Prepare all dilutions of the patient specimens before starting the test.
- **Immediately return unused strips to the pouch containing desiccants and seal securely to minimize exposure to water vapor.**
- Wash step: Good technique is critical. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 wells simultaneously. This speeds the process and provides for a more uniform incubation time.
- Careful timing is important. Incubation periods begin after dispensing reagents.

Assay procedure

1. ALL REAGENTS MUST BE BROUGHT TO ROOM TEMPERATURE (20-26°C) PRIOR TO BEGINNING THE ASSAY.
2. Label protocol record to indicate specimen placement in the microplate. It is good laboratory practice to test specimens in duplicate.

EN

3. **Qualitative determination:** use only Calibrator D.
Semi-quantitative determination: use Calibrators A - D as shown in the example below.

4. **Prepare a 1:51** dilution of the patient specimen by mixing 10 μ l of the patient specimen with **0.5 ml** of Serum Diluent.
5. Add **100 μ l** of Calibrators, Positive and Negative controls and diluted patient specimens to the appropriate microwells indicated on the protocol record.
Note: Include one well with **100 μ l** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank. The absorbance of this well should not be greater than 0.3.
6. Incubate **30 minutes** (\pm 5 minutes) at room temperature on a level surface.
7. Wash step: Thoroughly aspirate the contents of each well. Add 200-300 μ l of the **reconstituted** wash buffer to all wells then aspirate. Repeat this sequence thrice more for a total of four washes. Invert the plate and tap it on absorbent material to remove any residual fluid after the last wash. Do not dry wells completely.
8. Add 100 μ l of the Conjugate to each well.
9. Incubate the wells for **30 minutes** (\pm 5 min) at room temperature.
10. Wash step: Repeat step 7.
11. Add 100 μ l of Enzyme Substrate to each well.
12. Incubate the wells for **30 minutes** (\pm 5 min) at room temperature.
13. Add 100 μ l of Stop Solution to each well. Maintain the same sequence and timing of Stop Solution addition as was used for the Enzyme Substrate. Read the absorbance (OD) of each well at 405nm within one hour of stopping the reaction.
14. Read the absorbance (OD) of each well at 405nm using a single or dual wave-length microplate reader against the reagent blank set at zero absorbance.

QUALITY CONTROL

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be <0.3 . The Calibrator A should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be <20 EU/ml. If the test is run in duplicate, take the mean of the two readings to determine the EU/ml. While performing Qualitative determinations, the optical density of the Calibrator D must be greater than that of the negative control and lesser than the absorbance of the positive control. For semi-quantitative determinations, the positive control must give values in the range stated on the vial.

EN

RESULTS

Calculations

The concentrations of the patient samples can be determined by either of two methods:

1. QUALITATIVE DETERMINATION

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator D}} \times \text{EU/ml of Calibrator D} = \text{EU/ml Test Sample}$$

2. SEMI-QUANTITATIVE DETERMINATION

Plot absorbance of Calibrator A through D against their respective concentration on a linear-linear graph paper. Plot the concentration in EU/ml on the X-axis against the absorbance on the Y-axis and draw the best fit curve. Determine the concentrations of the patient samples from the curve against its corresponding absorbance value.

Calibrator

The calibrators are included to provide semi-quantitation and must be used with each test. Patient specimens containing higher antibody levels may give absorbance values greater than that of the Calibrator A. For determining accurate semi-quantitative values such specimens should be further diluted so they fall within the range of the calibrator curve when retested. For determining EU/ml, multiply the units obtained by the dilution factor. See sample standard curves (Figure 1) at the end of this document.

Interpretation

The following information serves only as a guide in the interpretation of the laboratory results. Each laboratory must determine its own normal values.

AGA-IgG Value		Interpretation
Child	Adult	
<28 EU/ml	<20 EU/ml	Neg (-)
>28 EU/ml	>20 EU/ml	Pos (+)
AGA-IgA Value		Interpretation
Child	Adult	
<23 EU/ml	<20 EU/ml	Neg (-)
>23 EU/ml	>20 EU/ml	Pos (+)

LIMITATIONS OF THE PROCEDURE

This test should not be performed on grossly hemolyzed, microbially contaminated or lipemic samples. This method should be used for testing human serum samples only. The results obtained serve only as an aid in the diagnosis and should not be interpreted as diagnostic in themselves. In IgA deficient CD patients IgA-AGA may be negative.

Expected Values

The expected values in a normal population are negative (<20 EU/ml for adults and <23-28 EU/ml for children). However, it has been determined that some apparently healthy, asymptomatic individuals may test positive for IgA or IgG anti-gliadin antibodies. **See Figure 2.**

Clinical Studies

The sensitivity and specificity of the IgG- and IgA-AGA tests were determined by testing patients with active CD (children=38), patients with DH (adults=10), patients with thyroid disorders (adults=200), patients with *pemphigus*

EN

and *pemphigoid* (adults=69), patients with other gastrointestinal disorders and 202 normal blood donors (children=113 and adults=89). Results from this study appear in Table 1 at the end of this document.

Using the ImmuLisa™ IgG and IgA anti-gliadin antibody test kits, the following incidences of AGA in normals, in patients with CD, and in other disease states were observed. IgA class AGA are specific for CD with a sensitivity of 90%. IgG-AGA, on the other hand are sensitive indicators of CD. See Table 2.

The algorithm in Table 3 serves as a guide in the interpretation of serological testing for CD.

Some investigators have taken the approach in Figure 3 to establish the diagnosis of CD using serological methods.

PERFORMANCE CHARACTERISTICS

Results obtained with the ImmuLisa™ IgG and IgA-AGA test were compared with the ImmuGlo™ anti-endomysial antibody (EMA) kit for diagnosing CD. The comparison included a total of 186 subjects. Included in the study were 34 patients on a gluten containing diet and 32 patients on gluten-free diet with a biopsy-confirmed diagnosis of CD as per ESPGAN criteria. **Results appear in Tables 4 and 5 at the end of this document.**

Precision:

Samples with known concentrations of AGA were assayed in 10 replicates over a period of two weeks. Intra-and inter-assay coefficient of variation (CV) appear in **Tables 6 and 7.**

Recovery

Samples with known AGA concentrations were mixed with appropriate dilutions of another positive sample with known amounts of AGA. Anti-gliadin antibody levels of the mixed samples were determined and from the values obtained the percent recovery calculated. The results appear in **Table 8.**

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REFERENCES • ΒΙΒΛΙΟΓΡΑΦΙΑ • LITERATUR • BIBLIOGRAPHIE • BIBLIOGRAFIA

1. Pare P, Douville P, Coron D, et al. Adult celiac sprue: changes in the pattern of clinical recognition. *J Clin Gastroenterol*; 10:395-400,1988.
2. Hall MJ, Cooper BT, Rooney H et al. Coeliac disease and malignancy of the duodenum: diagnosis by endoscopy, successful treatment of malignancy and response to gluten free diet. *Gut* 32:90-92, 1991.
3. Rautonen J, Rautonen N and Sauifahti E. Antibodies to gliadin in children with celiac disease. *Acta Paediatr Scand*; 1991, 80:1200-1205.
4. Bürgin-Wolff A, Lentze MJ. Relation of anti-gliadin antibodies to gluten sensitive enteropathy B. Diagnostic significance of antibodies against gliadin in celiac disease in children. In serologic diagnosis of celiac disease eds. Chorzelski TP, Beutner EH, Kumar V and Zalewski TP. CRC Press, Boca Raton FL, 1990 pp. 59-75.
5. Lerner A, Kumar V, Iancu TC. Immunological diagnosis of childhood coeliac disease: comparison between anti-gliadin, anti-reticulin and anti-endomysial antibodies. *Clin Exp Immunol*; 1994, 95:78-82.
6. Patinen P, Björkstén F, Malmström M et al. Salivary and serum anti-gliadin antibodies in dermatitis herpetiformis. *Eur J Oral Sci*; 1995, 103.
7. Walker-Smith JA, Guandalini S, Schmitz J, et al. Revised criteria for diagnosis of celiac disease. Report of working group of European Society of Pediatric Gastroenterology and Nutrition. *Arch Dis Child* 1990; 65:909-911.
8. Savilahti E, Perkkio M, Kalimo K, et al. IgA anti-gliadin antibodies: a marker of mucosal damage in childhood celiac disease. *Lancet* 1973; i:320-322.
9. Kumar V, Jain N, Lerner A, Beutner EH, Chorzelski TP, Lebenthal E. Comparative studies of different gliadin preparations in detecting anti-gliadin antibodies. *J Pediatr Gastroenterol Nutr* 1986; 5:730-734.
10. Montgomery AMP, Goka AKJ, Kumar PJ, et al. Low gluten diet in the treatment of adult coeliac disease: effect on jejunal morphology and serum anti-gluten antibodies. *Gut* 1988; 29:1564-1568.
11. Weiss JB, Austin RK, Schanfield MS, et al. Gluten-sensitive enteropathy. IgG heavy-chain (Gm) allotypes and the immune response to wheat gliadin. *J Clin Invest* 1983; 72:96-101.
12. Mearin ML, Konincx CR, Biemond I, et al. Influence of genetic factors on the serum levels of anti-gliadin antibodies in celiac disease. *J Pediatr Gastroenterol Nutr* 1984; 3:373-377.
13. Levenson SD, Austin RK, Dietler MD, et al. Specificity of anti-gliadin antibody in celiac disease. *Gastroenterology* 1985; 89:1-5.
14. Lerner A, Lebenthal E. The controversy of anti-gluten antibody (AGA) as a diagnostic tool in celiac disease. *J Pediatr Gastroenterol Nutr* 1991; 12:407-409.
15. Lebenthal E, Heitlinger LA. Gliadin antibodies in celiac disease. *J Pediatr* 1983; 102:711-712.
16. Kumar V, Lerner A, Jain N, Beutner EH. Are anti-gliadin antibodies specific for celiac disease? *J Pediatr Gastroenterol Nutr* 1984; 3:815.
17. Tucker NT, Barghuthy FS, Prihoda TJ, Kumar V, Lerner A, Lebenthal E. Anti-gliadin antibodies detected by enzyme-linked immunosorbent assay as a marker of childhood celiac disease. *J Pediatr* 1988; 113:286-289.
18. Unsworth DJ, Kieffer M, Holborow EJ, et al. IgA anti-gliadin antibodies in coeliac disease. *Clin Exp Immunol* 1981; 46:286-293.
19. Bürgin-Wolff A, Bertele RM, Berger R et al. A reliable screening test for childhood celiac disease: fluorescent immunosorbent test for gliadin antibodies. A prospective multicenter study. *J Pediatr* 1983; 102:655-660.
20. Kelly J, O'Farrelly C, Rees JPR, et al. Humoral response to alpha gliadin as serological screening test for coeliac disease. *Arch Dis Child* 1987; 62:469-473.
21. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control National Institutes of Health; 1993 (HHS Pub No [CDC] 93-8395).
22. Unsworth DJ. Serological diagnosis of gluten sensitive enteropathy. *J Clin Pathol*, 1996, 49:704-711.

Figure 1. Standard Curve



Figure 2. Incidence of Anti-Gliadin Antibodies

		Clinical Specificity	Clinical Sensitivity
Children	IgG-AGA	79%	94%
	IgA-AGA	97%	61%
Adults	IgG-AGA	88%	70%
	IgA-AGA	98%	60%

Adult patients with *gluten sensitive enteropathy* in our study were patients with *dermatitis hepetiformis*.

Table 2. Incidence of Anti-Gliadin Antibodies

Disease State		Antibody Specificity	
		IgG	IgA
Confirmed CD	On Gluten Diet	100%	90%
	On Gluten-Free Diet	0%	0%
Disease Controls		0%	0%
Normals		5%	0%

Table 3. CD Diagnosis Algorithm

EMA	ARA	AGA		anti-tTG	Interpretation
		IgG	IgA		
+	+	+	+	+	Positive
+	+	-	+	+	Positive
+	+	-	-	+	Positive
+	-	+	+	+	Positive
+	-	-	-	+	CD prob-ble
-	+	+	+	+	CD prob-ble
-	-	+	+	+	CD possible
-	-	+	-	-	CD unlikely
-	-	-	+	-	CD unlikely
-	-	-	-	-	Negative

Figure 3. Laboratory Testing Algorithm


Table 4. Immco™ AGA-IgA vs. EMA

		Immco™ AGA-IgA		
		POS (+)	NEG (-)	TOT (=)
Immco™	POS (+)	58	21	79
EMA	NEG (-)	5	102	107
	TOT (=)	63	123	186

relative specificity: 86%
 relative sensitivity: 73%
 relative agreement: 95%

Table 5. Immco™ AGA-IgG vs. EMA

		Immco™ AGA-IgG			
		POS (+)	NEG (-)	TOT (=)	
Immco™	POS (+)	77	3	80	
EMA	NEG (-)	26	81	107	
		TOT (=)	103	84	187

relative specificity: 84%
 relative sensitivity: 96%
 relative agreement: 75%

Table 6. AGA-IgG Precision

Anti-AGA-IgG (EU/ml)	Intra-plate CV	Inter-plate CV
1. 68.4	11.9%	8.8%
2. 130.1	7.2%	12.0%

Table 7. AGA-IgA Precision

Anti-AGA-IgA (EU/ml)	Intra-plate CV	Inter-plate CV
1. 83.3	10.7%	12.0%
2. 197.4	4.5%	8.4%

Table 8. AGA-IgG Recovery

IgG	EU/ml added	EU/ml obtained	% Recovery
1.	88	90	102%
2.	79	81	103%
3.	73	76	104%

For technical assistance please contact:

IMMCO Diagnostics, Inc.

60 Pineview Drive

Buffalo, NY 14228-2120

Telephone: (716) 691-0091

Fax: (716) 691-0466

Toll Free USA/Canada: 1-800-537-TEST

E-Mail: info@immco.com

or your local product distributor

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante
Autorizzato/Representante Autorizado/Représentant Autorisé

EMERGO Group, Inc.

Molenstraat 15, 2513 BH, The Hague,

The Netherlands

Tel (+31) 345 8570, Fax (+31) 346 7299

www.emergogroup.com