## Anti-Liver/Kidney Microsomal-1 (LKM1) Antibody ELISA

IVD

### PRODUCT INSERT

REF 1168 LKM1 ELISA 96 Determinations

#### INTENDED USE

An enzyme-linked immuno-sorbent assay (ELISA) for the detection and semi-quantitation of anti-liver/kidney microsomal-1 (LKM1) antibodies in human serum. The presence of LKM1 antibody can be used as an adjunct to clinical and other laboratory findings in the diagnosis of *autoimmune hepatitis*.

#### SUMMARY AND EXPLANATION

*Autoimmune hepatitis* (AIH) is a distinct chronic inflammatory liver disease, characterized by the attack of the immune system directed against "self" antigens, especially those expressed in the liver<sup>1,2</sup>. It occurs in both sexes and all age groups, however, women are more likely victims of AIH than men. In women, 70 % of diagnosed cases of AIH occur between the ages of 15 and 40.<sup>1-4</sup>

*Hepatomegaly* and *splenomegaly* are the most common pathological findings associated with AIH. Abnormalities of the immune system that mark AIH include autoantibodies to liver antigens, *hyper-gammaglobulinemia*, and an increased CD4/CD8 ratio in peripheral blood and liver. Liver-Kidney Microsomal (LKM1) antibodies can be induced not only by autoimmune mechanisms, but also by drugs such as *tienic acid*, *dihydralazine*, *halothane*, *phenytoin*, *phenobarbital*, *carbamazepine* and by Hepatitis C and D infections.<sup>4-6</sup>

The International Autoimmune Hepatitis Group categorizes AIH into two separate disease groups: Type 1 and Type 2. This distinction is based on the presence of marker autoantibodies in serum of affected patients.<sup>5</sup> AIH Type 1 is characterized by antinuclear autoantibodies (ANA's) and smooth-muscle antibodies (SMA's). Type 1 is the more common type of AIH, accounting for 60-70% of patients with AIH. Type 2 is a somewhat rarer disease (prevalance of AIH Type 2 is about 10 cases per million) characterized by the presence of autoantibodies against microsomal antigens of liver and kidney (LKM) and the absence of ANA and SMA<sup>3</sup>. 70% of these patients are young, from 2 to 14 years old at onset. The consequences of Type 2 AIH are severe: typically, acute onset progresses rapidly to cirrhosis and liver failure.

The autoantigen associated with LKM1 antibodies is cytochrome P4502D6. Antibodies to LKM1 occur predominantly in patients with AIH but can also be detected in some patients with HCV infection<sup>7-10</sup>. The Immco™ anti-LKM1 assay uses specific antigen to assure a high degree of specificity for AIH.<sup>7-10</sup>

\* manufactured under patent # 5,830,067

#### PRINCIPALS OF THE PROCEDURE

LKM1 antigen is bound to the wells of a polystyrene microwell plate followed by blocking the unreacted sites to reduce non-specific binding. Controls, calibrators and diluted patient sera are added to separate wells, allowing any LKM1 antibodies present to bind to the immobilized antigen. Unbound sample is washed away and an enzyme labeled anti-human IgG conjugate is added to each well. These enzyme conjugated antibodies bind specifically to the human immunoglobulin of the IgG class. After washing away any unbound conjugate, specific enzyme substrate (pNPP) is then added to the wells. After stopping the enzymatic reaction, the intensity of color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 405 nm. Results are expressed in Enzyme Units per milliliter (EU/ml).

#### REAGENTS

##### Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.** Do not use if reagent is not clear or if a precipitate is present. All

EN

reagents must be brought to room temperature (20-25°C) prior to use. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date. Reconstitute the wash buffer to 1 liter with distilled or deionized water. Coated microwell strips are for one time use only.

### Precautions

For *in vitro* Diagnostic Use. All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials<sup>11</sup>. **WARNING** - Sodium azide (NaN<sub>3</sub>) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

**Instructions should be followed exactly as they appear in this kit insert to ensure valid results.** Do not interchange kit components with those from other sources other than the same catalog number from Immco Diagnostics Inc. Use good laboratory techniques to minimize microbial and chemical contamination. Do not use after expiration date.

### Materials provided

Immulin™ LKM1 ELISA **REF** 1168


Kits contain sufficient reagents to perform 96 determinations each.

12 x 8	<b>MICROPLATE</b> <b>LKM1</b>	<b>Microplate</b> with individual breakaway microwells coated with LKM1 antigen.
1 x 1.5 ml	<b>CALIBRATOR A</b> <b>LKM1</b> *	Ready to use <b>Calibrator A</b> ( <i>green cap</i> ). Human serum containing antibodies to LKM1.
1 x 1.5 ml	<b>CALIBRATOR B</b> <b>LKM1</b> *	Ready to use <b>Calibrator B</b> ( <i>violet cap</i> ). Human serum containing antibodies to LKM1.
1 x 1.5 ml	<b>CALIBRATOR C</b> <b>LKM1</b> *	Ready to use <b>Calibrator C</b> ( <i>blue cap</i> ). Human serum containing antibodies to LKM1.
1 x 1.5 ml	<b>CALIBRATOR D</b> <b>LKM1</b> *	Ready to use <b>Calibrator D</b> ( <i>yellow cap</i> ). Human serum containing antibodies to LKM1.
1 x 1.5 ml	<b>CONTROL +</b> <b>LKM1</b> *	Ready to use <b>Positive Control</b> ( <i>red cap</i> ). Contains human serum positive for LKM1.
1 x 1.5 ml	<b>CONTROL -</b> *	Ready to use <b>Negative Control</b> ( <i>white cap</i> ). Contains human serum.
1 x 12 ml	<b>IgG-CONJ</b> <b>ALKPHOS</b> *	Ready to use <b>anti-human Alk. Phos. Conjugate</b> . Color coded pink.
1 x 60 ml	<b>DIL</b> *	Ready to use <b>Serum Diluent</b> . Color coded blue.
1 x 12 ml	<b>SUBSTRATE</b> *	Ready to use <b>Enzyme Substrate</b> . Contains pNPP. <b>Protect from light.</b>
1 x 12 ml	<b>STOP</b>	Ready to use <b>Stop Solution</b> .
2 x	<b>BUF</b> <b>WASH</b>	Powder <b>Wash Buffer</b> . Reconstitute to one liter each.

\* Contains < 0.1% NaN<sub>3</sub>

EN

**Symbols used on labels:**

-  Lot number
-  Catalog number
-  Use by
-  Storage temperature
-  Read instructions for use
-  In vitro diagnostic use
-  Manufacturer
-  Number of Tests

**Materials Required But Not Provided**

- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Deionized or distilled water
- Microplate reader capable of reading absorbance values at 405 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm.
- Squeeze bottle to hold diluted wash buffer
- Timer
- Absorbent paper
- Automatic microplate washer capable of dispensing 200 µl

**SPECIMEN COLLECTION AND HANDLING**

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples.


**PROCEDURE**

**Procedural Notes**

- Before starting the assay read these instructions carefully.
- Bring all reagents and samples to room temperature (20-26°C) for 30 minutes. Return materials to refrigerator immediately after use.
- Prepare all dilutions of the patient specimens before starting the test.
- **Immediately return unused strips to the pouch containing desiccants and seal securely to minimize exposure to water vapor.**
- Wash step: Good technique is critical. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 wells simultaneously. This speeds the process and provides for a more uniform incubation time.
- Careful timing is important. Incubation periods begin after dispensing reagents.

### Assay procedure

- ALL REAGENTS MUST BE BROUGHT TO ROOM TEMPERATURE (20-26°C) PRIOR TO BEGINNING THE ASSAY.**
- Label protocol record to indicate specimen placement in the microplate. It is good laboratory practice to test specimens in duplicate.
- Qualitative determination:** use only Calibrator D. **Semi-quantitative determination:** use Calibrators A - D as shown in the example below.

- Prepare a **1:101** dilution of the patient specimen by mixing **5 µl** of the patient specimen with **0.5 ml** of Serum Diluent.
- Add **100 µl** of Calibrators, Positive and Negative controls and diluted patient specimens to the appropriate microwells indicated on the protocol record.  
**Note:** Include one well with **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank.
- Incubate **30 minutes** ( $\pm$  5 minutes) at room temperature on a level surface.
- Wash step: Thoroughly aspirate the contents of each well. Add 200-300µL of the **reconstituted** wash buffer to all wells then aspirate. Repeat this sequence thrice more for a total of four washes. Invert the plate and tap it on absorbent material to remove any residual fluid after the last wash. Do not dry wells completely.
- Add 100µL of the Conjugate to each well.
- Incubate the wells for 30 minutes.
- Wash step: Repeat step 7.
- Add 100µL of Enzyme Substrate to each well.
- Incubate for 30 minutes at room temperature.
- Add 100µL of Stop Solution to each well. Maintain the same sequence and timing of Stop Solution addition as was used for the Enzyme Substrate. Read the absorbance (OD) of each well at 405nm within one hour of stopping the reaction.
- Read the absorbance (OD) of each well at 405nm using a single or dual wavelength microplate reader against the reagent blank set at zero absorbance.

### Quality Control

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be  $<0.3$ . The Calibrator A should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be  $<20$  EU/ml. If the test is run in duplicate, take the mean of the two readings to determine the EU/ml. While performing Qualitative determinations, the optical density of the Calibrator D must be greater than that of the negative control and lesser than the absorbance of the positive control. For semi-quantitative determinations, the positive control must give values in the range stated on the vial.

EN

## RESULTS

### Calculations

The concentrations of the patient samples can be determined by either of two methods:

#### 1. QUALITATIVE DETERMINATION

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator D}} \times \text{EU/ml of Calibrator D} = \text{EU/ml Test Sample}$$

#### 2. SEMI-QUANTITATIVE DETERMINATION

Plot absorbance of Calibrator A through D against their respective concentration on a linear-linear graph paper. Plot the concentration in EU/ml on the X-axis against the absorbance on the Y-axis and draw the best fit curve. Determine the concentrations of the patient samples from the curve against its corresponding absorbance value.

### Calibrator

The calibrators are included to provide semi-quantitation and must be used with each test. Patient specimens containing higher antibody levels may give absorbance values greater than that of the Calibrator A. For determining accurate semi-quantitative values such specimens should be further diluted so they fall within the range of the calibrator curve when retested. For determining EU/ml, multiply the units obtained by the dilution factor.

### Interpretation

The following information serves only as a guide in the interpretation of the laboratory results. Each laboratory must determine its own normal values.

#### **anti-LKM1 value Interpretation**

<20 EU/ml	Negative
20-25 EU/ml	Indeterminate (Borderline)
>25 EU/ml	Positive

### LIMITATIONS OF THE PROCEDURE

The results obtained serve only as an aid in the diagnosis of AIH type 2 and should not be interpreted as diagnostic in themselves. Serum of some AIH type 2 patients may be negative for anti-LKM1 antibodies

### EXPECTED VALUES

There is no single sensitive and specific marker for AIH. Diagnosis is made from a combination of clinical, biochemical, serological and histological features. The presence of a variety of autoantibodies helps in the identification and diagnosis of AIH. Based on the type of antibodies present, AIH can be divided into subtypes. Autoantibodies to LKM1 are specific indicators of type 2 AIH (see tables 1 and 2 at the end of this document).

### Precision

Two LKM1 positive sera were tested with the Immco™ LKM1 to determine inter- and intra-assay variability. See table 3 at the end of this document.

### Recovery:

Samples with known LKM1 concentrations were mixed with appropriate dilutions of another positive sample with known amounts of LKM1. LKM1 levels of the mixed samples were determined and from the values obtained the percent recovery calculated. See table 4 at the end of this document.

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
## REFERENCES • ΒΙΒΛΙΟΓΡΑΦΙΑ • LITERATUR • BIBLIOGRAPHIE • BIBLIOGRAFIA

1. Meyer zum Büschenfelde KH, Dienes, HP. Autoimmune hepatitis, definition – classification – histopathology – immunopathogenesis. *Virchows Arch* 429:1-12, 1996.
2. Czaja AJ. Frequency and nature of the variant syndromes of autoimmune liver disease. *Hepatology* 28:360-365, 1998
3. Chazouilleres O, Wendum D, Serfaty L et al. Primary biliary cirrhosis – autoimmune hepatitis overlap syndrome: clinical features and response to therapy. *Hepatology* 28:297-301, 1998.
4. van den Berg AP. Autoimmune hepatitis: pathogenesis, diagnosis and treatment. *Scand J Gastroenterol* 225:66-9, 1998.
5. Czaja AJ. Diagnosis and therapy of autoimmune liver disease. *Med Clin North Am* 80:973994, 1996.
6. Yamamoto AM, Cresteil D, Homberg JC and Alvarez F. Characterization of anti-liverkidney microsome antibody (Anti-LKM1) from hepatitis C virus – positive and – negative sera. *Gastroenterol* 104:1762-1767, 1993.
7. Choudhuri K, Gregorio GV, Mieli-Vergani G, and Vergani D. Immunological cross-reactivity to multiple autoantigens in patients with liver kidney microsomal type 1 autoimmune hepatitis. *Hepatology* 28:1177-1181, 1998.
8. Ma Y, Peakman M, Lobo-Yeo A et al. Differences in immune recognition of cytochrome P4502D6 by liver kidney microsomal (LKM) antibody in autoimmune hepatitis and chronic hepatitis C virus infection. *Clin Exp Immunol* 97:94-99, 1994.
9. Manns MP, Griffin KJ, Sullivan KF and Johnson EF. LKM-1 autoantibodies recognize a short linear sequence in P450IID6, a cytochrome P-450 monooxygenase. *J. Clin Invest* 88:1370-1378, 1991.
10. Gueguen M, Boniface O, Bernard O et al. Identification of the main epitope on human cytochrome P450IID6 recognized by anti-liver kidney microsome antibody. *J Autoimm* 4:607-615, 1991.
11. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health, 1993 [HHS Pub. No. (CDC) 93-8395].

**Figure 1. Sample Standard Curve**

**Table 1. Antibodies in Autoimmune Hepatitis**

Type of Antibody	<i>Antibodies in Autoimmune Hepatitis</i>	
	Incidence (%)	
	Type 1	Type 2
Anti-nuclear Antibody (ANA)	60	0
Smooth Muscle Antibody	70	0
Liver Kidney Microsomal Antibody (LKM1)	0	100

From van Den Berg, Autoimmune Hepatitis, pathogenesis, diagnosis and treatment, Scand J Gastroenterol, 1998; 33 suppl 225; 66-69.

**Table 2. Autoantibodies to LKM 1 antigen by ELISA can be detected in almost all cases of AIH and in none with HCV infection induced hepatitis.**

	Indirect IF	ELISA
AIH	100%	>95%
HCV Infection	1.5%	0%

**Table 3: Precision**

	inter-assay	intra-assay
	%CV	%CV
Sample 1	7.2	5.0
Sample 2	8.5	7.0

**Table 4: Recovery**

	LKM1 conc. added (EU/ml)	LKM1 conc. obtained (EU/ml)	% Recovery
Sample 1	121.2	128.7	106.2
Sample 2	102.4	112.1	109.5
Sample 3	63.3	61.4	97.0

*For technical assistance please contact:*

**IMMCO Diagnostics, Inc.**

**60 Pineview Drive**

**Buffalo, NY 14228-2120**

**Telephone: (716) 691-0091**

**Fax: (716) 691-0466**

**Toll Free USA/Canada: 1-800-537-TEST**

**E-Mail: [info@immco.com](mailto:info@immco.com)**

*or your local product distributor*

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante  
Autorizzato/Representante Autorizado/Représentant Autorisé

EMERGO Group, Inc.

Molenstraat 15, 2513 BH, The Hague,  
The Netherlands

Tel (+31) 345 8570, Fax (+31) 346 7299

[www.emergogroup.com](http://www.emergogroup.com)