EN

## Anti-Thyroid Peroxidase (TPO) Antibody ELISA

IVD

CLIA Complexity: High  
CDC Analyte Identification Code: 6135  
CDC Test System Identification Code: 28442

### PRODUCT INSERT

REF 1132 TPO Antibody ELISA 96 Determinations

#### INTENDED USE

An enzyme linked immunosorbent assay (ELISA) for the qualitative detection and quantitation of autoantibodies against TPO in human serum to aid in the diagnosis of certain thyroid disorders - *Hashimoto's*, *Grave's disease* and non-toxic goiter.

#### SUMMARY AND EXPLANATION

The clinical spectrum of autoimmune thyroid disorders is broad and the affected patient may be either hyper, hypo or even euthyroid. There are two major forms of autoimmune thyroid disorders, *Grave's disease* and *Hashimoto's thyroiditis*. Thyroid autoimmune reactions can also occur in other thyroid abnormalities such as sporadic and endemic goiter, Plummer's disease and endocrine ophthalmopathy<sup>1-5</sup>.

Patients with these disorders are often associated with the presence of autoantibodies to thyroglobulin and thyroid peroxidase (TPO) antigens. Thyroglobulin (Tg) is a 660 kD homodimeric glycoprotein which functions as thyroid prohormone. TPO is a membrane bound enzyme of 105 kD that catalyses thyroid hormone biosynthesis. Thyroxine and tri-iodo thyronine are generated by the TPO catalysed iodination and coupling at specific homogenic tyrosines. Measurement of autoantibodies to thyroglobulin and TPO are important in the diagnosis of autoimmune thyroid diseases. These antibodies can be measured by various methods such as indirect immunofluorescence, passive hemagglutination and ELISA.

#### PRINCIPLES OF PROCEDURE

The ELISA is performed as a solid phase immunoassay. Microwells are coated with native TPO antigen followed by blocking the unreacted sites to reduce non-specific binding. Controls, calibrators and patient serum samples are incubated in the antigen coated wells which allows specific antibodies present in the serum to bind to the TPO antigen. Unbound antibody and other serum proteins are removed by washing the microwells.

Bound antibodies are detected by adding an enzyme labeled anti-human IgG conjugate to the microwells. Unbound conjugate is removed by washing.

Specific enzyme substrate (pNPP) is then added to the wells and the presence of antibodies is detected by a color change produced by the conversion of pNPP substrate to a colored reaction product. The reaction is stopped and the intensity of the color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 405 nm. Results are expressed in International units per milliliter (IU/ml).

#### REAGENTS

##### Storage and Preparation

Store all reagents at 2° - 8°C. **Do not freeze.**

Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20°-25°C) prior to use.

When stored at 2° - 8°C, the reconstituted wash buffer is stable until the kit expiration date. Reconstitute the wash buffer to 1 liter with distilled or deionized water.

Coated microwell strips are for one time use only.

### Precautions

For *in vitro* Diagnostic Use. All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials<sup>6</sup>.

WARNING - Sodium azide (NaN<sub>3</sub>) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

**Instructions should be followed exactly as they appear in this kit insert to ensure valid results.** Do not interchange kit components with those from other sources other than the same catalog number from Immco Diagnostics Inc. Follow good laboratory practices to minimize microbial and cross contamination of reagents when handling. Do not use beyond expiration date on the label.

### Materials provided

Immulin™ TPO Antibody ELISA **REF** 1132

Kits contain sufficient reagents to perform 96 determinations each.

12 x 8	<b>MICROPLATE</b> <b>TPO</b>	<b>Microplate</b> with individual breakaway microwells coated with TPO antigen.
1 x 1.5 ml	<b>CALIBRATOR A</b> <b>TPO</b> *	Ready to use <b>Calibrator A</b> ( <i>green cap</i> ). Human serum containing antibodies to TPO.
1 x 1.5 ml	<b>CALIBRATOR B</b> <b>TPO</b> *	Ready to use <b>Calibrator B</b> ( <i>violet cap</i> ). Human serum containing antibodies to TPO.
1 x 1.5 ml	<b>CALIBRATOR C</b> <b>TPO</b> *	Ready to use <b>Calibrator C</b> ( <i>blue cap</i> ). Human serum containing antibodies to TPO.
1 x 1.5 ml	<b>CALIBRATOR D</b> <b>TPO</b> *	Ready to use <b>Calibrator D</b> ( <i>yellow cap</i> ). Human serum containing antibodies to TPO.
1 x 1.5 ml	<b>CONTROL +</b> <b>TPO</b> *	Ready to use <b>Positive Control</b> ( <i>red cap</i> ). Contains human serum positive for TPO antibodies.
1 x 1.5 ml	<b>CONTROL -</b> *	Ready to use <b>Negative Control</b> ( <i>whitecap</i> ). Contains human serum.
1 x 12 ml	<b>IgG-CONJ</b> <b>ALKPHOS</b> *	Ready to use <b>anti-human Alk. Phos. Conjugate</b> . Color coded pink.
2 x 60 ml	<b>DIL</b> *	Ready to use <b>Serum Diluent</b> . Color coded blue.
1 x 12 ml	<b>SUBSTRATE</b> *	Ready to use <b>Enzyme Substrate</b> . Contains pNPP. <b>Protect from light.</b>
1 x 12 ml	<b>STOP</b>	Ready to use <b>Stop Solution</b> .
2 x	<b>BUF</b> <b>WASH</b>	<b>Powder Wash Buffer</b> . Reconstitute to one liter each.

\* Contains <0.1% NaN<sub>3</sub>

EN

**Symbols used on labels:**

-  Lot number
-  Catalog number
-  Use by
-  Storage temperature
-  Read instructions for use
-  In vitro diagnostic use
-  Manufacturer
-  Number of Tests

**Materials Required But Not Provided**

- Deionized or distilled water
- Squeeze bottle to hold diluted wash buffer
- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Timer
- Absorbent paper towels
- Microplate reader capable of reading absorbance values at 405 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm.
- Automatic microplate washer capable of dispensing 200 µl

**SPECIMEN COLLECTION AND HANDLING**

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples.

**PROCEDURE**

**Procedural Notes**


- Before starting with the assay read carefully the product insert.
- Let serum specimens and test reagents equilibrate at room temperature before starting with the test procedure. Return all unused specimens and reagents to refrigerator immediately after use.
- All dilutions of the patient samples should be prepared prior to starting with the assay.
- Good washing technique is critical. If washing is performed manually, adequate washing is accomplished by directing a forceful stream of wash buffer with a wide tip wash bottle across the entire microplate. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 wells simultaneously. This speeds the process and provides for a more uniform incubation time.
- For all steps, careful control of timing is important. The start of all incubation periods begins with the completion of reagent addition.

EN

- Addition of all samples and reagents should be performed at the same rate and in the same sequence.
- Remove required microwell strips from the pouch and carefully reseal the pouch to prevent condensation in the unused wells. Return pouch immediately to refrigerator.

### Test Method

- Step 1** Let all reagents and specimens equilibrate at room temperature.
- Step 2** Label protocol sheet to indicate sample placement in the wells. It is good laboratory practice to run samples in duplicate.
- Step 3** For a **qualitative determination** use only the Ready to Use Low Calibrator D (*vial with yellow cap*).  
**or** For a **semi-quantitative determination** use the Ready to Use Calibrators A through D as depicted in the sample layout below.

- Step 4** Prepare a **1:201** dilution of the patient samples by mixing **5 µl** of the patient sera with **1.0 ml** of Serum Diluent.
- Step 5** Remove the required microwells from pouch and return unused strips in the sealed pouch to refrigerator. Securely place the microwells into the extra provided holder .
- Step 6** Pipette **100 µl** of Ready to use Calibrators, Positive and Negative controls and diluted patient samples to the appropriate microwells as per protocol sheet.  
**Note:** Include one well which contains **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank.
- Step 7** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 8** Wash **4x** with wash buffer. For manual washing, fill each microwell with reconstituted wash buffer. Discard the fluid by inverting and tapping out the contents of each well or by aspirating the liquid from each well. To blot at the end of the last wash, invert strips and tap the wells vigorously on absorbent paper towels. For automatic washers, program the washer as per manufacturer's instructions.
- Step 9** Pipette **100 µl** of Conjugate into microwells.
- Step 10** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 11** Wash all microwells as in Step 8.
- Step 12** Pipette **100 µl** of Enzyme Substrate into each microwell in the same order and timing as for the Conjugate.
- Step 13** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 14** Pipette **100 µl** of Stop Solution into each microwell using the same order and timing as for the addition of the Enzyme Substrate. Read absorbance values within 1 hour from adding Stop Solution.
- Step 15** Read absorbance of each microwell at **405 nm** using a single or 405/630nm dual wavelength microplate reader against the reagent blank set at zero absorbance.

**Quality Control**

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be <0.3. The Calibrator A should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be <20 IU/ml. If the test is run in duplicate, the mean of the two readings should be taken for determining IU/ml. When performing Qualitative determinations, the optical density of the Calibrator D must be greater than that of the negative control and lesser than the absorbance of the positive control. For semi-quantitative determinations, the positive control must give values in the range stated on the vial.

**RESULTS**

**Calculations**

The concentrations of the patient samples can be determined by either of two methods:


**1. QUALITATIVE DETERMINATION**

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator D}} \times \text{IU/ml of Calibrator D} = \text{IU/ml Test Sample}$$

**2. SEMI-QUANTITATIVE DETERMINATION**

Plot absorbance of Calibrator A through D against their respective concentration on a linear-linear graph paper. Plot the concentration in IU/ml on the X-axis against the absorbance on the Y-axis and draw the best fit curve. Determine the concentrations of the patient samples from the curve against its corresponding absorbance value.

Anti-TPO (IU/ml) Immulisa™ Standard Curve

**Calibrator**

The Ready to Use Calibrators are included to provide semi-quantitation and must be used with each run. Patient samples containing higher antibody levels may give absorbance values greater than that of the Calibrator A. For determining accurate semi-quantitative values such serum sample should be further diluted so they fall within the range of the calibrator curve when retested. For determining IU/ml, multiply the units obtained by the dilution factor.

**Interpretation**

The following serves only as a guide in the interpretation of the laboratory results. These values were determined by testing 80 normal blood donors. The values depicted below are the mean of the normal plus 2SD. Each laboratory must determine its own normal values.

Anti-TPO values	Interpretation
≤20 IU/ml	Negative
>20 IU/ml	Positive

**LIMITATIONS OF THE PROCEDURE**


Thyroid autoantibodies may also be present in non-thyroid disorders. Test results obtained by this assay alone are not diagnostic for thyroid disease and should be considered in conjunction with iodine uptake and other standard thyroid tests and the clinical presentation of the patient.

**EXPECTED VALUES**

The expected values in a normal population are negative. However, 5-10% of apparently healthy, asymptomatic individuals may test positive for thyroid autoantibodies. The incidence of these autoantibodies increases with age. The following table and figure depicts the incidence of thyroid autoantibodies in normals and in various disease states using Immulisa™ antigen plates. Antibodies to TPO occur primarily in *Grave's*, *Hashimoto's*, toxic nodular goiter and in association with certain other autoimmune disorders, such as SLE.

**Significance of Thyroid Antibodies**

Ab Specificity	Disease Association	Indications for use
<b>Thyroglobulin (Tg)</b>	1) thyroid autoimmune diseases 2) incorrect thyroglobulin levels	complement to thyroglobulin levels
<b>Microsomal (TPO)</b>	Thyroid autoimmune disease	1) Goiter of unknown etiology 2) Hyperthyroidism

*N=Normals (99), G=Grave's (25), H=Hashimoto's (15), T=Toxic Multi Nodular Goiter (11), S=SLE (20), P=Primary Biliary Cirrhosis (9)*

**PERFORMANCE CHARACTERISTICS**

Results obtained with the Immulisa™ anti-TPO ELISA were compared with another commercially manufactured assay.

		Immulisa™		Total
		Positive	Negative	
Other ELISA	Positive	42	14	56
	Negative	4	62	66
	Total	46	76	122

Relative Agreement: 85.2%

Relative Sensitivity: 75%

Relative Specificity: 94%

**Precision:**

Three different anti-TPO antibody positive sera were tested with the Immulisa™ anti-TPO ELISA to determine inter- and intra-assay variability based on 10 replicates. The results are listed below:

	mean value IU/ml	inter-assay %CV	intra-assay %CV
Sample 1	196.6	10.8	6.6
Sample 2	155.8	12.3	8.0
Sample 3	110.9	18.0	7.5

**Recovery:**

Samples with known anti-TPO antibody concentrations were mixed with appropriate dilutions of another positive sample with known amounts. Anti-TPO antibody levels of the mixed samples were determined and from the values obtained the percent recovery calculated. The results are as follows:

	TPO-Ab conc. added (IU/ml)	TPO-Ab conc. obtained (IU/ml)	% Recovery
Sample 1	275.0	244.3	88.9
Sample 2	183.0	175.3	95.8
Sample 3	116.0	107.0	92.2

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**REFERENCES • ΒΙΒΛΙΟΓΡΑΦΙΑ • LITERATUR • BIBLIOGRAPHIE • BIBLIOGRAFIA**

1. Mooij P, Drexhage HA. Autoimmune thyroid disease. *Clini Lab Med.* 13:683-697, 1993.
2. Champion BR, Cooke A, Reyner DC. Thyroid autoimmunity. *Curr Opinion Immunol.* 4:770-778, 1992.
3. Lindstedt G, Berg G, Jansson S, et al Clinical use of laboratory thyroid tests and investigations. *J Int Fed Clin Chem* 6:136-141, 1994.
4. Pashke R, Vogg M, Swillens S, Usadel KH. Correlation of microsomal antibodies with the intensity of the intrathyroidal autoimmune process in Grave's disease. *J Clin Endocrinol Metab.* 77:939-943, 1993.
5. Sundbeck G, Eden S, Jagenburg R, et al. Prevalence of serum anti-thyroid peroxidase antibodies in 85 year-old women and men. *Clin Chem*; 41:707-712, 1995.
6. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health [HHS Pub. No. (CDC) 93-8395], 1993.

*For technical assistance please contact:*

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The Netherlands

Tel (+31) 345 8570, Fax (+31) 346 7299

[www.emergogroup.com](http://www.emergogroup.com)

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**IMMCO Diagnostics, Inc.**
EMERGO Group, Inc.