# ImmLiSa™ Anti-*Saccharomyces cerevisiae* Antibody (ASCA) IgA and IgG ELISA

**IVD** For *in vitro* diagnostic use

## PRODUCT INSERT

Catalog No. 1156A/G

48 Determinations

### INTENDED USE

An enzyme linked immunosorbent assay (ELISA) for the detection and semi-quantitation of anti-*Saccharomyces cerevisiae* antibodies (IgG or IgA) in human serum of patients with inflammatory bowel disorder (IBD) as an aid in the diagnosis of Crohn's disease (CD).

### SUMMARY AND EXPLANATION

Inflammatory bowel diseases (UC and CD) are chronic remitting disorders with unpredictable courses and variable responses to therapy. The diagnosis of IBD is based on clinical, radiological, endoscopic and histological findings. Immunological, environmental, infectious (*Yersinia enterocolitica*) and genetic factors have been postulated to increase the risk for developing IBD<sup>1</sup>.

CD is primarily a disease of older children and young adults and is rare in infancy. Few epidemiological studies have been conducted to find the incidence of CD in the general population. In the Scandinavian population alone, the incidence of IBD is 7 per 100,000, with CD at 1.3 per 100,000<sup>2</sup>. CD can affect any part of the intestine, and the lesions are often patchy and occasionally very extensive. Histologically, CD is characterized by a transmural inflammatory infiltrate composed predominantly of lymphocytes and macrophages. Macrophage aggregates are often seen in biopsy specimens, but well developed granulomas are present in only 50% resected specimens<sup>3</sup>. The usual surgical procedures for CD include segmental resection and stricturoplasty<sup>2</sup>. CD can also manifest itself along with other disorders like Henoch-Schonlein purpura<sup>4</sup>. Vasculitis is an uncommon complication of CD, however, few cases have been reported of cutaneous polyarteritis nodosa with CD<sup>4</sup>. A milder form of perianal CD exists and can vary from mild asymptomatic disease to a severe disabling disorder<sup>5</sup>.

ASCA of IgG and IgA isotype are present in 60% of cases with diagnosed CD<sup>1</sup>. The ASCA's are directed against phosphopeptidomannans present in the yeast (*S. cerevisiae*) cell walls<sup>6</sup>. ASCA may be elicited due to molecular mimicry and priming by a high mannose-containing bacterial or viral antigen or an "auto antigenic (self) molecule"<sup>1,7</sup>. It has been observed that levels of ASCA in CD cases are independent of disease activity, duration and treatment<sup>8</sup>.

ASCA occur frequently in CD patients with positive familial history. ASCA are not considered specific markers of genetic susceptibility, but their presence indicates a genetic predisposition<sup>9</sup>. Studies in monozygotic twins and the Jewish population indicates an influence of genetic factors on the pathogenesis of CD<sup>10</sup>. The mode of inheritance seems to be complex and heterogenous, with candidate loci of IBD located on chromosomes 2, 6, 12, and 16<sup>11</sup>.

Recently, serological testing algorithms have been proposed for diagnosis of IBD and to improve discrimination between CD and UC<sup>1</sup>. In summary, all suspected cases of IBD must be tested for ASCA (IgG and IgA).

### PRINCIPALS OF THE PROCEDURE

The ASCA test is performed as a solid phase immunoassay (ELISA). Microwells are coated with *Saccharomyces cerevisiae* phosphopeptidomannan antigen followed by blocking the unreacted sites to reduce nonspecific binding. Controls, calibrators and patient serum samples are incubated in the antigen coated wells which allows ASCA present in the serum to bind. Unbound antibody and other serum proteins are removed by washing the microwells. Antibodies bound to the microwells are detected by adding enzyme labeled anti-human IgG or IgA conjugates to the wells. These enzyme conjugated antibodies bind specifically to the human immunoglobulin of the appropriate

class. Unbound enzyme conjugate is removed by washing. Specific enzyme substrate (pNPP) is then added to the wells and the presence of antibodies to *Saccharomyces cerevisiae* is detected by a color change produced by the conversion of the pNPP substrate. The reaction is stopped and the intensity of color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 405 nm. Results are expressed in enzyme units per milliliter (EU/ml).

## REAGENTS

### Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.** Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date. Reconstitute the wash buffer to 1 liter with distilled or deionized water. Coated microwell strips are for one time use only.

### Precautions

For *in vitro* Diagnostic Use. All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials<sup>17</sup>. **WARNING** - Sodium azide (NaN<sub>3</sub>) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

**Instructions should be followed exactly as they appear in this kit insert to ensure valid results.** Do not interchange kit components with those from other sources other than the same catalog number from IMMCO DIAGNOSTICS.

Use good laboratory techniques to minimize microbial and chemical contamination. Do not use after expiration date.

### Materials provided

ImmuLisa™ IgA-ASCA ELISA REF 1156A/G

Kits contain sufficient reagents to perform 48 determinations each.

12 x 8	MICROPLATE ASCA	<b>Microplate</b> with individual breakaway microwells coated with <i>S. cereviceae</i> phosphopeptidomannan.
4 x 1.5 ml	CALIBRATOR A-D ASCA IgA <sup>*A</sup>	Ready to use <b>Calibrator A-D</b> ( <i>green cap</i> ). Human serum containing antibodies to ASCA IgA.
4 x 1.5 ml	CALIBRATOR A-D ASCA IgG <sup>*G</sup>	Ready to use <b>Calibrator A-D</b> ( <i>green cap</i> ). Human serum containing antibodies to ASCA IgG.
1 x 1.5 ml	CONTROL+ ASCA IgA <sup>*A</sup>	Ready to use <b>Positive Control</b> ( <i>red cap</i> ). Contains human serum positive for ASCA IgA.
1 x 1.5 ml	CONTROL+ ASCA IgG <sup>*G</sup>	Ready to use <b>Positive Control</b> ( <i>red cap</i> ). Contains human serum positive for ASCA IgG.
1 x 1.5 ml	CONTROL - <sup>*</sup>	Ready to use <b>Negative Control</b> ( <i>white cap</i> ). Contains human serum.
1 x 12 ml	IgA-CONJ ALKPHOS <sup>*A</sup>	Ready to use <b>anti-human Alk. Phos. Conjugate</b> . Color coded pink.
1 x 12 ml	IgG-CONJ ALKPHOS <sup>*G</sup>	Ready to use <b>anti-human Alk. Phos. Conjugate</b> . Color coded pink.

1 x 60 ml **DIL**\*

Ready to use **Serum Diluent**. Color coded blue.

1 x 12 ml **SUBSTRATE** \*

Ready to use **Enzyme Substrate**. Contains pNPP. **Protect from light.**

1 x 12 ml **STOP**

Ready to use **Stop Solution**.

2 **BUF|WASH**

Powder **Wash Buffer**. Reconstitute to one liter each.

\* Contains < 0.1% NaN<sub>3</sub>

A IgA ELISA component. Use only with other IgA ELISA components.

B IgG ELISA component. Use only with other IgG ELISA components.

### Materials Required But Not Provided

- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Deionized or distilled water
- Microplate reader capable of reading absorbance values at 405 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm.
- Squeeze bottle to hold diluted wash buffer
- Timer
- Absorbent paper
- Automatic microplate washer capable of dispensing 200 µl

### SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples.



### PROCEDURE

#### Procedural Notes

- Before starting the assay read these instructions carefully.
- Bring all reagents and samples to room temperature (20-25°C) for 30 minutes. Return materials to refrigerator immediately after use.
- Prepare all dilutions of the patient specimens before starting the test.
- **Immediately return unused strips to the pouch containing desiccants and seal securely to minimize exposure to water vapor.**
- Wash step: Good technique is critical. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 wells simultaneously. This speeds the process and provides for a more uniform incubation time.
- Careful timing is important. Incubation periods begin after dispensing reagents.

#### Assay procedure

1. **ALL REAGENTS MUST BE BROUGHT TO ROOM TEMPERATURE (20-25°C) PRIOR TO BEGINNING THE ASSAY.**
2. Label protocol record to indicate specimen placement in the microplate. It is good laboratory practice to test specimens in duplicate.
3. **Qualitative determination:** use only Calibrator D.  
**Semi-quantitative determination:** use Calibrators A - D as shown in the example below.

**QUALITATIVE DETERMINATION****SEMI-QUANTITATIVE DETERMINATION**

4. Prepare a **1:51** dilution of the patient specimen by mixing **10 µl** of the patient specimen with **0.5 ml** of Serum Diluent.
5. Add **100 µl** of Calibrators, Positive and Negative controls and diluted patient specimens to the appropriate microwells indicated on the protocol record.

**Note:** Include one well with **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank. The absorbance of this well should not be greater than 0.3.

6. Incubate 30 minutes ( $\pm$  5 minutes) at room temperature on a level surface.
7. Wash step: Thoroughly aspirate the contents of each well. Add 200-300µL of the **reconstituted** wash buffer to all wells then aspirate. Repeat this sequence thrice more for a total of four washes. Invert the plate and tap it on absorbent material to remove any residual fluid after the last wash. Do not dry wells completely.
8. Add 100µL of the Conjugate to each well.
9. Incubate the wells for 30 minutes.
10. Wash step: Repeat step 7.
11. Add 100µL of Enzyme Substrate to each well.
12. Incubate for 30 minutes at room temperature.
13. Add 100µL of Stop Solution to each well. Maintain the same sequence and timing of Stop Solution addition as was used for the Enzyme Substrate. Read the absorbance (OD) of each well at 405nm within one hour of stopping the reaction.
14. Read the absorbance (OD) of each well at 405nm using a single or dual wavelength microplate reader against the reagent blank set at zero absorbance.

**Quality Control**

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be  $<0.3$ . The Calibrator A should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be  $<20$  EU/ml. If the test is run in duplicate, take the mean of the two readings to determine the EU/ml. While performing Qualitative determinations, the optical density of the Calibrator D must be greater than that of the negative control and lesser than the absorbance of the positive control. For semi-quantitative determinations, the positive control must give values in the range stated on the vial.

## Calculations

The concentrations of the patient samples can be determined by either of two methods:

### 1. QUALITATIVE DETERMINATION

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator D}} \times \text{EU/ml of Calibrator D} = \text{EU/ml Test Sample}$$

### 2. SEMI-QUANTITATIVE DETERMINATION

Plot absorbance of Calibrator A through D against their respective concentration on a linear-linear graph paper. Plot the concentration in EU/ml on the X-axis against the absorbance on the Y-axis and draw the best fit curve. Determine the concentrations of the patient samples from the curve against its corresponding absorbance value.

#### Calibrator

The ready to use calibrators are included to provide semi-quantitation and must be used with each run. Patient samples containing higher antibody levels may give absorbance values greater than that of the Calibrator A. For determining accurate semi-quantitative values such serum sample should be further diluted so they fall within the range of the calibrator curve when retested. For determining EU/ml, multiply the units obtained by the dilution factor. Sample standard curves appear at the end of this document.

#### Interpretation

The following serves only as a guide in the interpretation of the laboratory results. The values depicted below were determined by testing 64 normal blood donors and represent the mean of the normals plus 3SD. Each laboratory must determine its own normal values.

ASCA Values	EU/ml	Interpretation
IgA or IgG	< 20	Negative
IgA or IgG	20-25	Borderline
IgA or IgG	> 25	Positive

#### LIMITATIONS OF THE PROCEDURE

The Immulisa™ Test should not be performed on grossly hemolyzed, microbially contaminated, lipemic or icteric samples. The method should be used for testing human serum samples only. Testing for both isotypes of ASCA is strongly recommended. Testing for only one and not both may lead to false negative results. A negative ASCA result does not rule out the presence of Crohn's disease. A negative ASCA antibody result does not rule out the presence of ASCA antibodies, because the concentration of antibody may be below the detection limit of the assay. Furthermore, a diagnosis cannot be made on the basis of ASCA results alone. The results of other laboratory tests and clinical findings must also be considered. The presence of immune complexes or other immunoglobulin aggregates in the patient sample may cause an increased level of non-specific binding and produce false positives in this assay. ASCA also occur in patients with UC, their first degree relatives and mixed families with CD and UC<sup>3,9</sup>. They are also present to a lesser extent in some other autoimmune diseases<sup>13</sup>. Assay performance has not been established for pediatric CD and UC patients.

## EXPECTED VALUES

The expected values in a normal population are negative (<20 EU/ml). However, it has been determined that some apparently healthy, asymptomatic individuals may test positive for IgA or IgG ASCA antibodies<sup>9</sup>. Combined evaluation of ASCA and perinuclear anti-

**ImmuLisa™ ASCA-IgA**

		Positive	Negative	Total
Disease Diagnosis	CD	12	18	30
	Other	0	36	36
	Total	12	54	66

Sensitivity: 40%                      Specificity: 100%

**ImmuLisa™ ASCA-IgG**

		Positive	Negative	Total
Disease Diagnosis	CD	18	12	30
	Other	6	31	36
	Total	23	43	66

Sensitivity: 60%                      Specificity: 86%

neutrophil cytoplasmic antibody (pANCA) for all suspected cases of IBD is suggested, to increase the positive predictive value in comparison to individual assessment of CD cases by ASCA and UC cases by pANCA<sup>14, 15</sup>. See Tables I and II at the end of this document.

## PERFORMANCE CHARACTERISTICS

Samples from patients with clinical diagnosis of CD (30 specimens) and UC (30 specimens) were tested along with normal human sera (6 specimens). Results obtained with the ImmuLisa™ anti-ASCA ELISA assays follow:

The same sample set was used in a comparative evaluation of the ASCA ImmuLisa™ (IgG and IgA) against another commercially-available ASCA ELISA. Results are depicted in table III at the end of this document.

### Cross Reactivity Study:

Sera from patients suffering from various conditions and individuals positive for certain autoantibodies were tested for ASCA with the ImmuLisa assay. Results appear in table IV at the end of this document.

### Precision:

In ten replicates, three different ASCA positive sera were tested with the ImmuLisa to determine the intra-assay variation, and three positive sera were tested to determine the inter-assay variation. The results appear in table V at the end of this document.

### Reportable Range:

Four ASCA positive samples were tested with the ImmuLisa to determine precision at even intervals through the reportable range of the assays. The results appear in table VI at the end of this document.

### Recovery:

Samples with known ASCA concentrations were mixed with appropriate dilutions of another positive sample with known amounts of ASCA. ASCA levels of the mixed samples were determined and from the values obtained the percent recovery calculated. The results appear in table VII at the end of this document.

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
## **REFERENCES•REFERENCIAS•LITERATUR•RIFERIMENTI**

1. Peeters M, Joossens S, Vermeire S et. al. Diagnostic value of anti-Saccharomyces cerevisiae and antineutrophil cytoplasmic autoantibodies in inflammatory bowel disease. *Am J Gastroenterol.* 96: 730-734. 2001.
2. Lindberg E, Lindquist B, Hildebrand H et. al. Inflammatory bowel disease in children and adolescents in Sweden, 1984-1995. *J. Pediatr. Gastroenterol Nutr.* 30:259-264, 2000
3. Seibold F, Stich O, Hufnagl R, Kamil S, Scheurlen M. Anti-Saccharomyces cerevisiae antibodies in inflammatory bowel disease: A family study. *Scand J Gastroenterol;* 36: 196-20,2001.
4. Saulsbury F.T., Hart. M. Crohn's disease presenting with Henoch-Schonlein purpura. *J. Pediatr. Gastroenterol. Nutr.* 31: 173-175, 2000.
5. Pashankar D, Schreiber R.A., Israel. D. M. Perianal Crohn's disease in Infancy. *J. Pediatr. Gastroenterol. Nutr.* 31: 80-82, 2000.
6. Sutton C.L, Yang H, Li Z, et. al. Familial expression of anti-Saccharomyces cerevisiae mannan antibodies in affected and unaffected relatives of patients with Crohn's disease. *Gut.* 46:58-63, 2000.
7. McKay D.M. Bacterial superantigens: Provocateurs of gut dysfunction and inflammation. *Trends Immunol.* 22:497-501, 2001.
8. Hoffenberg, E.J. Fidanza S., Sauaia A. Serologic testing for inflammatory bowel disease. *J. Pediatr.* 134: 447-452, 1999.
9. Annese V, Andreoli A, Andriulli A, et al. Familial expression of anti-Saccharomyces cerevisiae Mannan antibodies in Crohn's disease and ulcerative colitis: a GISC study. *Am J Gastroenterol Aug;*96:2407-12, 2001.
10. Yang H, McElree C, Roth M-P, et al. Familial empiric risk for inflammatory bowel disease: differences between Jews and non Jews. *Gut;* 34: 517-24, 1993.
11. Curran ME, Lau KF, Hampe J, et al. Genetic analysis of inflammatory bowel disease in a large European cohort supports linkage to chromosomes 12 and 16. *Gastroenterology;* 115: 1066-71, 1998.
12. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health; (HHS Pub. No. [CDC] 93-8395), 1993.
13. Reddy R.K. Colombel J-F, Poulin D et. al. Anti-Saccharomyces cerevisiae antibodies in autoimmune liver disease. *Am. J. Gastroenterol.* 96: 252-253, 2001.
14. Quinton. JF, Sendid B, Reumaux D. et. al. Anti-Saccharomyces cerevisiae mannan antibodies combined with antineutrophil cytoplasmic antibodies in inflammatory bowel disease: prevalence and diagnostic role. *Gut.* 42:788-79, 1998.
15. Kim BG, Kim YS, Kim JS et al. Diagnostic Role of Anti-Saccharomyces cerevisiae mannan antibodies combined with antineutrophil cytoplasmic antibodies in patients with inflammatory bowel disease. *Dis Colon Rectum,* 45: 1062-1069, 2002.
16. Vermeire. S, Joossens. S, Peeters. M. et. al. Comparative study of ASCA (Anti-Saccharomyces cerevisiae antibody) assays in inflammatory bowel disease. *Gastroenterology.* 120: 827-833, 2001.
17. Ruummele FM, Targan SR, Levy G, et al. Diagnostic accuracy of serological assays in pediatric inflammatory bowel disease. *Gastroenterology.* 115(4):822-9, 1998.

**Figure I: ASCA Standard Curves**

**Table I: Diagnostic significance of ASCA and pANCA in IBD<sup>14</sup>**

	pANCA + PPV %	ASCA + PPV %	ASCA - pANCA + PPV %	pANCA - ASCA + PPV %
<b>Crohn's Disease</b>		89		96
<b>Ulcerative Colitis</b>	74		92.5	

*pANCA: atypical pANCA NPV: Negative predictive value PPV: Positive predictive value*

**Table II: Prevalence of ASCA in Crohn's Disease Versus Other GI Disorders**

Study	CD (n)	ASCA Pos. (%)	UC (n)	ASCA Pos. (%)	Control (n)	ASCA Pos (%)
Quinton e al <sup>14</sup>	100	61	101	12	163 (healthy) 27 (GI)	0.6 11
Hoffenberg et al <sup>8</sup>	20	60	25	12	74 (GI & liver)	5
Ruemmele et al <sup>17</sup>	130	55	35	6	78 (GI)	5
Peeters et al <sup>1</sup>	407	60	147	14	157 (healthy) 74(non-IBD)	3 11
Kim et al <sup>15</sup>	85	49	77	20	20 (healthy)	10

**Table III: IMMCO ASCA ELISA vs. Other ELISA**

	ASCA IgA		ASCA IgG	
	Sensitivity (%)	Specificity (%)	Sensitivity (%)	Specificity (%)
<b>ImmuLisa ASCA</b>	40	100	60	86
<b>Other ELISA I</b>	52	86	62	82

**Table IV: IMMCO ASCA ELISA vs. Other ELISA**

Specimen Type	n	IgA Pos	n	IgG Pos
Disease controls	12	3	13	1
ds DNA Antibodies	23	3	23	1
Smooth Muscle Antibodies	22	3	24	3
Immune Complex	21	3	24	1
Rheumatoid Factor	22	3	22	0

**Table V: Precision**

Positive Sera	IgA		Positive Sera	IgG	
	Intra-assay %CV	Inter-assay %CV		Intra-assay %CV	Inter-assay %CV
1. 167 EU/ml	3.6%	8.2%	1. 170 EU/ml	5.5%	11.5%
2. 129 EU/ml	5.1%	10.5%	2. 140 EU/ml	3.5%	4.6%
3. 108 EU/ml	3.4%	10.3%	3. 116 EU/ml	5.2%	7.4%

**Table VI: Reportable Range**

	IgA		IgG	
	Avg EU/ml	%CV	Avg EU/ml	%CV
Sample 1	172	6.8%	160	4.0%
Sample 2	94	6.0%	72	5.1%
Sample 3	37	7.0%	37	5.6%
Sample 4	10	10.3%	13	10.7%

**Table VII: Recovery**

	ASCA-IgA		% Recovery
	Ab. conc. added (EU/ml)	Ab. conc. obtained (EU/ml)	
Sample 1	129.5	130.2	101
Sample 2	59.8	60.6	101
Sample 3	21.0	22.2	106

	ASCA-IgG		% Recovery
	Ab. conc. added (EU/ml)	Ab. conc. obtained (EU/ml)	
Sample 1	130.4	128.4	98
Sample 2	56.5	57.7	102
Sample 3	24.6	26.3	107

**Table VIII: Clinical Studies**

Clinical Group	n	IgG		IgA		IgG or IgA		IgG and IgA	
		Pos	%	Pos	%	Pos	%	Pos	%
Crohn's Disease	39	24	62%	31	79%	36	92%	19	49%
Ulcerative Colitis	38	1	3%	9	24%	9	24%	1	3%
Controls	12	1	8%	3	25%	3	25%	1	8%

*For technical assistance please contact:*

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The Netherlands

Tel (+31) 345 8570, Fax (+31) 346 7299

[www.emergogroup.com](http://www.emergogroup.com)

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**IMMCO Diagnostics, Inc.**
EMERGO Group, Inc.