EN

## Anti-Tissue Transglutaminase Antibody (tTG) IgA ELISA

IVD

### PRODUCT INSERT

CLIA Complexity: High  
CDC Analyte Identification Code: 0546  
CDC Test System Identification Code: 28569

REF 1157 IgA-tTG ELISA 96 Determinations

### INTENDED USE

An enzyme linked immunosorbent assay (ELISA) for the detection and semi-quantitation of anti-Tissue Transglutaminase IgA antibodies in human serum to aid in the diagnosis of celiac disease.

### SUMMARY AND EXPLANATION

Celiac disease (CD) is a chronic digestive disorder caused by an autoimmune reaction against gluten, a protein found in wheat. It is characterized by damage to the absorptive villi and hyperplasia of the crypts in the small intestine. Symptoms include diarrhea, weight loss and anemia from malabsorption of iron and folic acid. Active CD increases risks of gastrointestinal malignancy<sup>1</sup>, infertility<sup>2</sup>, osteoporosis<sup>3</sup> and epilepsy and other neurological syndromes<sup>4</sup>. Celiac Disease is treated with a diet free of gluten.

CD is common in the European Caucasian population. As the disease is genetically mediated and the European and American Caucasian population share a common ancestry, actual prevalence of CD in the United States is suspected to be greater than reported in the past. Serological tests have been developed to aid in diagnosis of CD (i.e., anti-endomysial antibody (EMA), anti-reticulín antibody IFA's and anti-gliadin antibody ELISA). Recently, Dietrich et al<sup>5</sup> identified the antigen associated with EMA antigen to be tissue transglutaminase (tTG). This made it possible to develop an ELISA test for EMA antibodies that would be less subjective than typical immunofluorescence methods.

Transglutaminase (EC 2.3.2.13) are a diverse family of Ca<sup>2+</sup> dependent enzymes with distinct genes, structures, and biological functions<sup>6</sup>. These enzymes are ubiquitous and highly conserved across species<sup>7</sup>. Of all the human transglutaminases, tTG is most widely distributed. Although primarily an intracellular enzyme, tTG can be secreted from cells and accumulates in the extracellular matrix<sup>8</sup>. tTG may have a role in stabilizing the extracellular matrix, functioning as a highly efficient biological glue. Autoantibodies against tTG may disrupt or interfere with its normal functioning, resulting in the histological and pathological changes seen in CD. tTG is an important analyte for the diagnosis of CD<sup>5,9,10,11</sup>.

### PRINCIPLES OF PROCEDURES

Tissue Transglutaminase antigen is bound to the wells of a polystyrene microwell plate followed by blocking the unreacted sites to reduce non-specific binding. Controls, calibrators and diluted patient sera are added to separate wells, allowing any anti-tTG antibodies present to bind to the immobilized antigen. Unbound sample is washed away and an enzyme labeled anti-human IgA conjugate is added to each well. These enzyme conjugated antibodies bind specifically to the human immunoglobulin of the IgA class. After washing away any unbound conjugate, specific enzyme substrate (pNPP) is then added to the wells. After stopping the enzymatic reaction, the intensity of color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 405 nm. Results are expressed in ELISA units per milliliter (EU/ml).

### REAGENTS

#### Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.** Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date. Reconstitute the wash buffer to 1 liter with distilled or deionized water. Coated microwell strips are for one time use only.

### Precautions

For *in vitro* Diagnostic Use. All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials<sup>12</sup>.

**WARNING** - Sodium azide (NaN<sub>3</sub>) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

**Instructions should be followed exactly as they appear in this kit insert to ensure valid results.** Do not interchange kit components with those from other sources other than the same catalog number from Immco Diagnostics Inc.

Use good laboratory techniques to minimize microbial and chemical contamination. Do not use after expiration date.

### Materials Provided

Immulin<sup>TM</sup> IgA-tTG ELISA **REF** 1157


Kits contain sufficient reagents to perform 96 determinations each.

<b>12 x 8</b>	<b>MICROPLATE</b>   <b>tTG</b>	<b>Microplate</b> with individual breakaway microwells coated with tTG antigen.
<b>1 x 1.5 ml</b>	<b>CALIBRATOR A</b>   <b>tTG</b> *	Ready to use <b>Calibrator A</b> ( <i>green cap</i> ). Human serum containing antibodies to tTG.
<b>1 x 1.5 ml</b>	<b>CALIBRATOR B</b>   <b>tTG</b> *	Ready to use <b>Calibrator B</b> ( <i>violet cap</i> ). Human serum containing antibodies to tTG.
<b>1 x 1.5 ml</b>	<b>CALIBRATOR C</b>   <b>tTG</b> *	Ready to use <b>Calibrator C</b> ( <i>blue cap</i> ). Human serum containing antibodies to tTG.
<b>1 x 1.5 ml</b>	<b>CALIBRATOR D</b>   <b>tTG</b> *	Ready to use <b>Calibrator D</b> ( <i>yellow cap</i> ). Human serum containing antibodies to tTG.
<b>1 x 1.5 ml</b>	<b>CONTROL +</b>   <b>tTG</b> *	Ready to use <b>Positive Control</b> ( <i>red cap</i> ). Contains human serum positive for tTG.
<b>1 x 1.5 ml</b>	<b>CONTROL -</b> *	Ready to use <b>Negative Control</b> ( <i>whitecap</i> ). Contains human serum.
<b>1 x 12 ml</b>	<b>IgA-CONJ</b>   <b>ALKPHOS</b> *	Ready to use <b>anti-human Alk. Phos. Conjugate</b> . Color coded pink.
<b>1 x 60 ml</b>	<b>DIL</b> *	Ready to use <b>Serum Diluent</b> . Color coded blue.
<b>1 x 12 ml</b>	<b>SUBSTRATE</b> *	Ready to use <b>Enzyme Substrate</b> . Contains pNPP. <b>Protect from light.</b>
<b>1 x 12 ml</b>	<b>STOP</b>	Ready to use <b>Stop Solution</b> .
<b>2 x</b>	<b>BUF</b>   <b>WASH</b>	<b>Powder Wash Buffer</b> . Reconstitute to one liter each.

\* Contains <0.1% NaN<sub>3</sub>

EN

### Symbols used on labels:

-  Lot number
-  Catalog number
-  Use by
-  Storage temperature
-  Read instructions for use
-  In vitro diagnostic use
-  Manufacturer
-  Number of Tests

### Materials Required But Not Provided

- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Deionized or distilled water
- Microplate reader capable of reading absorbance values at 405 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm.
- Squeeze bottle to hold diluted wash buffer
- Timer
- Absorbent paper
- Automatic microplate washer capable of dispensing 200 µl

### SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples.


### PROCEDURE

#### Procedural Notes

- Before starting the assay read these instructions carefully.
- Bring all reagents and samples to room temperature (20-26°C) for 30 minutes. Return materials to refrigerator immediately after use.
- Prepare all dilutions of the patient specimens before starting the test.
- **Immediately return unused strips to the pouch containing desiccants and seal securely to minimize exposure to water vapor.**
- Wash step: Good technique is critical. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 wells simultaneously. This speeds the process and provides for a more uniform incubation time.
- Careful timing is important. Incubation periods begin after dispensing reagents.

### Assay procedure

1. **ALL REAGENTS MUST BE BROUGHT TO ROOM TEMPERATURE (20-26°C) PRIOR TO BEGINNING THE ASSAY.**
2. Label protocol record to indicate specimen placement in the microplate. It is good laboratory practice to test specimens in duplicate.
3. **Qualitative determination:** use only Calibrator D.  
**Semi-quantitative determination:** use Calibrators A - D as shown in the example below.

4. Prepare a **1:51** dilution of the patient specimen by mixing **10 µl** of the patient specimen with **0.5 ml** of Serum Diluent.
5. Add **100 µl** of Calibrators, Positive and Negative controls and diluted patient specimens to the appropriate microwells indicated on the protocol record. 3  
**Note:** Include one well with **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank. The absorbance of this well should not be greater than 0.3.
6. Incubate **30 minutes** ( $\pm$  5 minutes) at room temperature on a level surface.
7. Wash step: Thoroughly aspirate the contents of each well. Add 200-300µL of the **reconstituted** wash buffer to all wells then aspirate. Repeat this sequence thrice more for a total of four washes. Invert the plate and tap it on absorbent material to remove any residual fluid after the last wash. Do not dry wells completely.
8. Add 100µL of the Conjugate to each well.
9. Incubate the wells for 30 minutes.
10. Wash step: Repeat step 7.
11. Add 1 00µL of Enzyme Substrate to each well.
12. Incubate for 30 minutes at room temperature.
13. Add 100µL of Stop Solution to each well. Maintain the same sequence and timing of Stop Solution addition as was used for the Enzyme Substrate. Read the absorbance (OD) of each well at 405nm within one hour of stopping the reaction.
14. Read the absorbance (OD) of each well at 405nm using a single or dual wavelength microplate reader against the reagent blank set at zero absorbance.

### Quality Control

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be  $<0.3$ . Calibrator A should have an absorbance reading  $>1.0$ , otherwise the test must be repeated. The negative control must be  $<20$  EU/ml. If the

test is run in duplicate, take the mean of the two readings to determine the concentration of anti-hu tTG antibodies. While performing qualitative determinations, the optical density of Calibrator D must be greater than that of the negative control and less than the positive control. For semi-quantitative determinations, the positive control must give values in the range stated on the vial. We recommend borderline samples be tested with a fresh sample taken at a later date to ensure accuracy.

**RESULTS**

**Calculations**

The concentrations of the patient samples can be determined by either of two methods:


**1. QUALITATIVE DETERMINATION**

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator D}} \times \text{EU/ml of Calibrator D} = \text{EU/ml Test Sample}$$

**2. SEMI-QUANTITATIVE DETERMINATION**

Plot absorbance of Calibrator A through D against their respective concentration on a linear-linear graph paper. Plot the concentration in EU/ml on the X-axis against the absorbance on the Y-axis and draw the best fit curve. Determine the concentrations of the patient samples from the curve against its corresponding absorbance value.

**Anti-tTG IgA Standard Curve**

**Calibrator**

The calibrators are included to provide semi-quantitation and must be used with each test. Patient specimens containing higher antibody levels may give absorbance values greater than that of the Calibrator A. For determining accurate semi-quantitative values such specimens should be further diluted so they fall within the range of the calibrator curve when retested. For determining EU/ml, multiply the units obtained by the dilution factor.

**Interpretation**

The following information serves only as a guide in the interpretation of the laboratory results. Each laboratory must determine its own normal values.

Anti-tTG Value	Interpretation
<20 EU/ml	Neg (-)
20 – 25 EU/ml	Borderline
>25 EU/ml	Pos (+)

**LIMITATIONS OF THE PROCEDURE**

The Immulisa™ anti-tTG IgA assay should not be performed on grossly hemolyzed, microbially contaminated or lipemic samples. The method should be used for testing human serum samples only. Results obtained serve only

as an aid in the diagnosis and should not be interpreted as diagnostic in themselves. In IgA deficient CD patients anti-tTG may be negative.

### EXPECTED VALUES

The expected values in a normal population are negative (<20 EU/ml for adults and children). However, it has been determined that some apparently healthy, asymptomatic individuals may test positive for IgA anti-tTG antibodies.

The incidence of anti-tTG antibodies has been studied by various investigators and the findings are summarized in Table 1 on page 42.

The incidence and the levels of anti-tTG antibodies is dependent upon the diet status. The levels of these antibodies decrease and eventually will become negative in patients with CD who are on a gluten-free diet. Similarly, the levels of these antibodies will increase or the anti-tTG antibodies become positive when patients with CD who were on a gluten-free diet ingest a gluten-containing diet<sup>9,10</sup>.

### PERFORMANCE CHARACTERISTICS

The utility of the Immulisa™ anti-tTG IgA Antibody ELISA was determined by comparing the results with:

- another commercially available anti-tTG IgA ELISA method and
- ImmuGlo™ anti-endomysial antibody immunofluorescence method.

Normal Range: A total of 213 samples were tested. Of these, 64 were normal blood donors and 149, obtained from a reference lab, tested negative for endomysial antibodies. Of these 213 samples tested for antibodies to tTG, only 3 were tested positive, providing a 99% specificity of the assay.

#### Comparative Specificity and Sensitivity

A. Endomysial Antibodies vs. Immulisa™ tTG IgA Antibody ELISA: A total of 244 samples were tested for antibodies to endomysium using ImmuGlo™ anti-endomysial antibody kit and the results of anti-endomysial antibodies were compared with Immulisa™ anti-tTG IgA ELISA method. The results are summarized below:

		Immco™ anti-tTg IgA		
		POS (+)	NEG (-)	TOT (=)
Immco™ EMA	POS (+)	86	9	95
	NEG (-)	1	148	149
	TOT (=)	87	157	244

relative specificity: 96%  
relative sensitivity: 91%  
relative agreement: 99%

B. Immulisa™ anti-tTG ELISA vs. Another Commercial anti-tTG IgA ELISA Method: A total of 137 samples were tested on Immulisa™ anti-tTG IgA antibody kit and another commercially available anti-tTG IgA kits. The results were as follows:

		Immco™ anti-tTg IgA		
		POS (+)	NEG (-)	TOT (=)
Other ELISA	POS (+)	77	5	82
	NEG (-)	1	54	55
	TOT (=)	78	59	137

relative specificity: 95%  
relative sensitivity: 91%  
relative agreement: 98%

EN

C. Cross Reactivity: A total of 30 samples, 10 each positive for antinuclear, anti-basement membrane zone (BMZ) and anti-intercellular (IC) antibodies were tested for anti-tTG IgA antibodies. None were found positive.

**Precision:**

Samples with known concentrations of anti-tTG IgA were assayed in 10 replicates over a period of two weeks. Intra-and inter-assay coefficient of variation (CV) were as follows:

<b>Anti-tTG IgA</b>	<b>Intra-plate CV</b>	<b>Inter-plate CV</b>
1. 186 EU/ml	5.3%	4.9%
2. 118 EU/ml	2.3%	4.9%
3. 51 EU/ml	9.2%	7.4%

**Recovery**

Samples with known anti-tTG IgA concentrations were mixed with appropriate dilutions of another positive sample with known amounts of anti-tTG IgA. Anti-tTG IgA antibody levels of the mixed samples were determined and from the values obtained the percent recovery calculated. The results are as follows:

	<b>EU/ml added</b>	<b>EU/ml obtained</b>	<b>% Recovery</b>
1.	151	162	107%
2.	126	140	111%

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## REFERENCES • ΒΙΒΛΙΟΓΡΑΦΙΑ • LITERATUR • BIBLIOGRAPHIE • BIBLIOGRAFIA

1. Cooper BT, Holmes GK, Ferguson R, et al. Celiac disease and malignancy. *Medicine*; 1980, 59:249-261.
2. Sher KS, Jayanthi V, Probert CS, et al. Infertility, obstetric and gynecological problems in coeliac sprue. *Digestive Dis*; 1994, 12:186-190.
3. McFarlane XA, Bhalla AK, Reeves DE, et al. Osteoporosis in treated adult coeliac disease. *Gut*; 1995, 36:710-714.
4. Gobbi G, Bouquet F, Greco L, et al. Coeliac disease, epilepsy and cerebral calcifications. The Italian working group on coeliac disease and epilepsy. *Lancet*; 1992, 340:439-443.
5. Dietrich W, Ehnis T, Bauer M, et al. Identification of tissue transglutaminase as the autoantigen of celiac disease. *Nature Medicine*; 1997, 3:797-801.
6. Phillips MA, Stewart BE, Quin Q, et al. Primary structure of keratinocyte trans-glutaminase. *Proc Natl Acad Sci*; 1990, 87:9333-9337.
7. Aschlimann D, Mosher D, Paulsson M, et al. Tissue transglutaminase and factor XIII in cartilage and bone remodelling. *Sem Thromb and Hemo*; 1996, 22:437-443.
8. Jones RA, Nicholas B, Mian S, et al. Reduced expression of tissue transglutaminase in a human endothelial cell line leads to changes in cell spreading, cell adhesion and reduced polymerization of fibronectin. *J Cell Sci*; 1997, 110:2461-2472.
9. Sulkanen S, Haltunen T, Laurila K, et al. Tissue transglutaminase autoantibody enzyme-linked immunosorbent assay in detecting celiac disease. *Gastroenterology*; 1998, 115:1322-1328.
10. Dietrich W, Laag E, Schopper H, et al. Autoantibodies to tissue transglutaminase as predictors of celiac disease; *Gastroenterology* 1998, 115:1317-1321.
11. Molberg O, McAdam S, Korner R. Tissue transglutaminase selectively modifies gliadin peptides that are recognized by gut-derived T-cells in celiac disease. *Nature Medicine*; 1998, 4:713-717.
12. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control National Institutes of Health; 1993 (HHS Pub No [CDC] 93-8395).
13. Dietrich W, Laag E, Bruckner-Tuderman et al. Antibodies to tissue transglutaminase as serologic markers in patients with dermatitis herpetiformis. *J Invest Dermatol*; 1999, 113:133-136.


**Table 1: Anti-tTG Incidence: GSE**

Disease/Research	Disease (Sensitivity)	Control (Specificity)
Celiac Disease	104/106 (98%)	6/14 (96%)
Dietrich W et al <sup>10</sup>	129/136 (95%)	13/207 (94%)
Sulkanen S et al <sup>9</sup>		
Dermatitis Herpetiformis	43/61 (70%)	2/84 (98%)
Dietrich W et al <sup>13</sup>		

**Table 2: CD Diagnosis Algorithm**

EMA	ARA	AGA		anti-tTG	Interpretation
		IgG	IgA		
+	+	+	+	+	Positive
+	+	-	+	+	Positive
+	+	-	-	+	Positive
+	-	+	+	+	Positive
+	-	-	-	+	CD – probable
-	+	+	+	+	CD – probable
-	-	+	+	+	CD – possible
-	-	+	-	-	CD – unlikely
-	-	-	+	-	CD – unlikely
-	-	-	-	-	Negative


**Figure 1: Laboratory Testing Algorithm**

**Table 3: Incidence of Anti-tTG IgA Autoantibodies**

Disease Group	No. of Cases	# Pos	% Pos
<b>Celiac Disease (CD)</b>			
Gluten Diet (untreated)	43	39	91
Gluten free diet	63	13	21
<b>Dermatitis Herpetiformis (DH)</b>	66	45	68
<b>Insulin Dependent Diabetes</b>	136	18	13
<b>Other Disorders</b>			
Pemphigus	72	1	2
Pemphigoid	70	7	10
Thyroid Disorders	80	1	1
ANA Positive	51	0	0
AMA Positive	10	0	0
<b>Normals</b>			
Adults	91	3	5
Children	111	5	5

**Figure 2: Distribution of Anti-tTG IgA in Celiacs, Normals and Controls**

*For technical assistance please contact:*

**IMMCO Diagnostics, Inc.**

**60 Pineview Drive**

**Buffalo, NY 14228-2120**

**Telephone: (716) 691-0091**

**Fax: (716) 691-0466**

**Toll Free USA/Canada: 1-800-537-TEST**

**E-Mail: [info@immco.com](mailto:info@immco.com)**

*or your local product distributor*

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante  
Autorizzato/Representante Autorizado/Représentant Autorisé

EMERGO Group, Inc.

Molenstraat 15, 2513 BH, The Hague,

The Netherlands

Tel (+31) 345 8570, Fax (+31) 346 7299

[www.emergogroup.com](http://www.emergogroup.com)