For technical assistance please contact:

IMMCO Diagnostics, Inc.
60 Pineview Drive
Buffalo, NY 14228-2120
Telephone: (716) 691-0091
Fax: (716) 691-0466
Toll Free USA/Canada: 1-800-537-TEST
E-Mail: info@immcodiagnostics.com

or your local product distributor

EMERGO Europe
Molenstraat 15, 2513 BH, The Hague,
The Netherlands
Tel (+31) 345 8570, Fax (+31) 346 7299
www.emergogroup.com

NOV2012
Document No. PI5118

ImmuLisa Enhanced™ Cardiolipin Antibody (ACA) ELISAs

Cardiolipin/Phospholipid IgA, IgG or IgM Antibody ELISA

IVD For *in vitro* diagnostic use

PRODUCT INSERT

REF	5118A	ACA IgA ELISA	96 Determinations
REF	5118G	ACA IgG ELISA	96 Determinations
REF	5118M	ACA IgM ELISA	96 Determinations

INTENDED USE

ImmuLisa Enhanced™ Cardiolipin IgA Antibody (ACA) ELISA: Enzyme linked immunoassay (ELISA) for the qualitative or semi-quantitative detection of Cardiolipin IgA antibodies in human serum to aid in the diagnosis of antiphospholipid syndrome (APS) and APS associated with systemic lupus erythematosus (SLE) in conjunction with other laboratory tests and clinical findings.

ImmuLisa Enhanced™ Cardiolipin IgG Antibody (ACA) ELISA: Enzyme linked immunoassay (ELISA) for the qualitative or semi-quantitative detection of Cardiolipin IgG antibodies in human serum to aid in the diagnosis of antiphospholipid syndrome (APS) and APS associated with systemic lupus erythematosus (SLE) in conjunction with other laboratory tests and clinical findings.

ImmuLisa Enhanced™ Cardiolipin IgM Antibody (ACA) ELISA: Enzyme linked immunoassay (ELISA) for the qualitative or semi-quantitative detection of Cardiolipin IgM antibodies in human serum to aid in the diagnosis of antiphospholipid syndrome (APS) and APS associated with systemic lupus erythematosus (SLE) in conjunction with other laboratory tests and clinical findings.

Warning: Confirmed active or seropositive syphilis patients can have elevated anti-cardiolipin antibody (ACA) levels. To rule out syphilis, confirmatory tests should be performed.

SUMMARY AND EXPLANATION

Antiphospholipid antibodies are a heterogeneous group of autoantibodies against negatively charged phospholipids.¹ They are detected primarily by the anti-cardiolipin antibody (ACA) test, the biological false positive test for syphilis and the lupus anticoagulant test. These three tests detect related, but not necessarily identical antibodies. Thus, more than one of these tests is sometimes necessary to identify antiphospholipid antibodies. The ELISA assay format is a highly sensitive method for detection of antiphospholipid antibodies.²

The presence of anti-cardiolipin antibodies helps to identify patients at risk of venous and/or arterial thrombosis often accompanied by thrombocytopenia, a syndrome referred to as antiphospholipid syndrome.¹⁻¹² The antiphospholipid syndrome most commonly occurs in patients with systemic lupus erythematosus (SLE) or lupus-like disease where the criteria for SLE are not fulfilled.⁵⁻⁷ High levels of anti-cardiolipin antibodies occur in thrombosis, fetal loss, thrombocytopenia and several other disorders.¹⁻¹⁵ Low levels of anti-cardiolipin antibodies are found in a variety of clinical disorders which are unrelated to antiphospholipid syndrome. Therefore, low levels of these antibodies are of limited significance.

EN

IgG and IgA class anti-cardiolipin antibodies appear to be more closely associated with antiphospholipid syndrome than the IgM class antibodies. However, IgM antibodies appear to be more influenced by treatment.^{5,10} Low levels of IgM antibodies can be identified in other autoimmune diseases such as rheumatoid arthritis, primary Sjögren's Syndrome, drug induced lupus erythematosus, Lyme disease, and syphilis.^{8,10} The Calibrators and Positive Control used in the kit have been standardized against the Harris calibration standards from the Antiphospholipid Standardization Laboratory.¹⁷

PRINCIPLES OF PROCEDURE

Cardiolipin antigen is bound to the wells of a polystyrene microwell plate followed by blocking the unreacted sites to reduce non-specific binding. Controls, calibrators and diluted patient sera are added to separate wells, allowing any cardiolipin antibodies present to bind to the immobilized antigen. Unbound sample is washed away and an enzyme labeled anti-human IgG, IgA or IgM conjugate is added to each well. These enzyme conjugated antibodies bind specifically to the human immunoglobulin of the appropriate class. After washing away any unbound conjugate, specific enzyme substrate (TMB) is then added to the wells. After stopping the enzymatic reaction, the intensity of color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 450 nm. Results are expressed in units per milliliter (U/ml) according to the Harris classification. Units for IgG are expressed as GPL, for IgM as MPL and for IgA as APL.

REAGENTS

Store all reagents at 2-8°C. **Do not freeze.** Reagents are stable until the expiration date when stored and handled as directed.

Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use.

Reconstitute the wash buffer to 1 liter with distilled or deionized water. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date.

Coated microwell strips are for one time use only. Unused microwell strips should be carefully resealed in the pouch containing desiccants to prevent condensation and stored at 2-8°C.

Precautions

All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However, human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials.¹⁸

Stop Solution is a dilute sulfuric acid solution. Sulfuric acid (H₂SO₄) is poisonous and corrosive. Do not ingest and avoid contact with skin and eyes. Avoid exposure to bases, metals or other compounds that may react with acids.

TMB Enzyme Substrate contains an irritant that may be harmful if inhaled, ingested or absorbed through the skin. Do not ingest and avoid contact with skin and eyes.

Instructions should be followed exactly as they appear in this kit insert to ensure valid results. Do not interchange kit components with those from other sources. Follow good laboratory practices to minimize microbial and cross contamination of reagents when handling. Do not use kit components beyond expiration date on the labels.

Materials provided

ACA IgA ELISA
ACA IgG ELISA
ACA IgM ELISA

REF	5118A
REF	5118G
REF	5118M

Kits contains sufficient reagents to perform 96 determinations.

REFERENCES

1. Harris EN et al. Antiphospholipid antibodies - antibodies with a difference. *Ann Rev Med.* 1988. 39:261-271.
2. Harris EN et al. Use of an enzyme-linked immunosorbent assay and inhibition studies to distinguish between antibodies to cardiolipin from patients with syphilis or autoimmune disorders. *J Infect Dis.* 1988. 157:23-31.
3. Lockshin MD. Anti-cardiolipin antibody. *Arth Rheum.* 1987. 30:471-472.
4. Derksen RHHM et al. Discordant effects of prednisone on anticardiolipin antibodies and the lupus anti-coagulant. *Arth Rheum.* 1986. 29:1295-1296.
5. Alarcón-Segovia D et al. Antiphospholipid antibodies and the antiphospholipid syndrome in systemic lupus erythematosus. A perspective analysis of 500 consecutive patients. *Medicine;* 1987. 68:353-365.
6. Triplett DA. Antiphospholipid antibodies and recurrent pregnancy loss. *Am J Reprod Immunol.* 1989. 20:52-67.
7. Montecucco C et al. Hematological abnormalities associated with anticardiolipin antibodies. *Hematologica.* 1989. 74:195-204.
8. Asherson RA et al. The "primary" antiphospholipid syndrome: major clinical and serological features. *Medicine.* 1989. 68:366-374.
9. Walport MJ. Pregnancy and antibodies to phospholipids. *Ann Rheum Dis.* 1989. 48:795-797.
10. Machworth-Young CG. Antiphospholipid antibodies and disease. *Medicine.* 1989. 72:767-777.
11. Harris EN et al. Antiphospholipid antibodies. *Clin Rheumatol.* 1985. 1:591-609.
12. Lubbe WF et al. Fetal survival after prednisone suppression of maternal lupus anticoagulant. *Lancet.* 1983. 1:1361-1363.
13. Hogan MJ et al. Lupus anticoagulant, antiphospholipid antibodies and migraine. *Can J Neurological Sci.* 1988. 15:420-425.
14. Harris EN et al. Antiphospholipid antibodies: method of detection. *Am J Reproduct Immunol.* 1992. 28:208-210.
15. Silveira LH et al. IgA, IgG and IgM anticardiolipin antibodies in black patients with systemic lupus erythematosus. *Arthr Rheum 56th Annual Meeting.* 1992. S125.
16. Harris, NE. The second international anti-cardiolipin standardization workshop/The Kingston antiphospholipid antibody study (KAPS) group. *Am J Clin Path.* 1990. 94:476-484.
17. Ichikawa K et al. A chimeric antibody with the human γ 1 constant region as a putative standard for assays to detect IgG β 2-glycoprotein I-dependent anticardiolipin and anti- β 2-glycoprotein I antibodies. *Arthritis Rheum.* 1999. 42:2461-70.
18. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control National Institutes of Health. 2007; (HHS Pub No [CDC] 93-8395).
19. Miyakis S et al. International consensus statement on an update of the classification criteria for definite antiphospholipid syndrome (APS). *J Thromb Haemost.* 2006. 4: 295-306.
20. Levy RA et al. Characteristics of IgG antiphospholipid antibodies in patients with systemic lupus erythematosus and syphilis. *J Rheumatol.* 1990. 17:1036-1041.
21. Ong SG et al. IgG anti-beta 2 glycoprotein I antibodies in Malaysian patients with antiphospholipid syndrome and systemic lupus erythematosus: prevalence and clinical correlations. *Clin Rheumatol.* 2002. 21:382-5.
22. Ebrahim RA et al. Antiphospholipid syndrome among Bahraini patients. *Saudi Med J.* 2005. 26:488-90.
23. Amoroso A et al. Antibodies to anionic phospholipids and anti-beta2-GPI: association with thrombosis and thrombocytopenia in systemic lupus erythematosus. *Hum Immunol.* 2003. 64:265-73.
24. Vianna JL et al. Comparison of the primary and secondary antiphospholipid syndrome: a European Multicenter Study of 114 patients. *Am J Med.* 1994. 96:3-9.
25. Marai I et al. Anti-cardiolipin and anti-beta2-glycoprotein I (beta2GP-I) antibody assays as screening for anti-phospholipid syndrome. *Hum Antibodies.* 2003. 12:57-62.
26. Gharavi AE et al. Antiphospholipid antibody tests. *Manual of clinical laboratory immunology sixth edition.* Rose NR et al (ed). 2002. 973-980.
27. Alarcon-Segovia D. et al. Antiphospholipid antibodies and the antiphospholipid syndrome in systemic lupus erythematosus. A perspective analysis of 500 consecutive patients. *Medicine* 1989. 68: 353-365.
28. Asherson RA et al. The "primary" antiphospholipid syndrome: major clinical and serological features. *Medicine.* 1989. 68:366-374.
29. Asherson RA et al. Anticardiolipin antibodies-clinical associations. *Post Med J.* 1986. 62, 1081-1087.
30. Gezer S. Antiphospholipid Syndrome. *Dis Mon.* 2003. 49: 691-742.
31. Bick RL. Antiphospholipid Thrombosis Syndromes. *Clin Appl Thromb/Hemosta.* 1997. 7; 4: 241-258.

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12 x 8

MICROPLATE|ACA

Microplate with individual breakaway microwells. Coated with cardiolipin antigen and bovine/human B2GP1. Ready for use.

1 x 1.75 ml

CONTROL|+|ACA-A

Ready to use Positive Control (red cap) for [REF] 5118A. Contains human serum positive for cardiolipin IgA antibodies. The expected concentration range in APL is printed on the label.

1 x 1.75 ml

CONTROL|+|ACA-G

Ready to use Positive Control (red cap) for [REF] 5118G. Contains human serum positive for cardiolipin IgG antibodies. The expected concentration range in GPL is printed on the label.

1 x 1.75 ml

CONTROL|+|ACA-M

Ready to use Positive Control (red cap) for [REF] 5118M. Contains human serum positive for cardiolipin IgM antibodies. The expected concentration range in MPL is printed on the label.

1 x 1.75 ml

CONTROL|-

Ready to use **Negative Control** (white cap). Contains human serum.

5 x 1.75 ml

CALIBRATOR|A|ACA-A

Ready to use set of 5 Calibrators for [REF] 5118A. Calibrator A (green cap) 120 APL, Calibrator B (violet cap) 60 APL, Calibrator C (blue cap) 30 APL, Calibrator D (yellow cap) 15 APL, and Calibrator E (orange cap) 1 APL. Derived from human serum containing ACA IgA. Concentrations in APL are printed on the labels.

CALIBRATOR|B|ACA-A

CALIBRATOR|C|ACA-A

CALIBRATOR|D|ACA-A

CALIBRATOR|E|ACA-A

5 x 1.75 ml

CALIBRATOR|A|ACA-G

Ready to use set of 5 Calibrators for [REF] 5118G. Calibrator A (green cap) 120 GPL, Calibrator B (violet cap) 60 GPL, Calibrator C (blue cap) 30 GPL, Calibrator D (yellow cap) 15 GPL, and Calibrator E (orange cap) 1 GPL. Derived from human serum containing ACA IgG. Concentrations in GPL are printed on the labels.

CALIBRATOR|B|ACA-G

CALIBRATOR|C|ACA-G

CALIBRATOR|D|ACA-G

CALIBRATOR|E|ACA-G

5 x 1.75 ml

CALIBRATOR|A|ACA-M

Ready to use set of 5 Calibrators for [REF] 5118M. Calibrator A (green cap) 120 MPL, Calibrator B (violet cap) 60 MPL, Calibrator C (blue cap) 30 MPL, Calibrator D (yellow cap) 15 MPL, and Calibrator E (orange cap) 1 MPL. Derived from human serum containing ACA IgM. Concentrations in MPL are printed on the labels.

CALIBRATOR|B|ACA-M

CALIBRATOR|C|ACA-M

CALIBRATOR|D|ACA-M

CALIBRATOR|E|ACA-M

1 x 15 ml

IgA-CONJ|HRP

HRP goat anti-human IgA Conjugate for [REF] 5118A. Ready for use. Color coded pink.

1 x 15 ml

IgG-CONJ|HRP

HRP goat anti-human IgG Conjugate for [REF] 5118G. Ready for use. Color coded pink.

1 x 15 ml

IgM-CONJ|HRP

HRP goat anti-human IgM Conjugate for [REF] 5118M. Ready for use. Color coded pink.

1 x 60 ml

DIL

Serum Diluent. Ready for use.

1 x 15 ml

SUBSTRATE|TMB

TMB enzyme substrate. Ready for use. **Protect from light.**

1 x 15 ml

STOP|H2SO4

Stop Solution. Ready for use.

2 x vials

BUF|WASH

Powder Wash Buffer. **Reconstitute to one liter each.**

1 x

Protocol Sheets

Optional Components

1 x 60ml **BUF** **WASH**Liquid concentrated Wash Buffer. **Reconstitute to one liter.**

Symbols used on labels

	Lot number
	Catalog number
	In vitro diagnostic use
	Use by
	Storage temperature
	Read instructions before use
	Number of tests
	Manufacturer

Materials Required But Not Provided

- Deionized or distilled water
- Squeeze bottle to hold diluted wash buffer
- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Timer
- Absorbent paper towels
- Microplate reader capable of reading absorbance values at 450 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm
- Automatic microplate washer capable of dispensing 200 µl

SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples. It is recommended that frozen specimens be tested within one year.

PROCEDURE

Procedural Notes

- Carefully read the product insert before starting the assay.
- Let patient specimens and test reagents equilibrate to room temperature before starting with the test procedure. Return all unused specimens and reagents to refrigerator immediately after use.
- Remove required microwell strips from the pouch and carefully reseal the pouch to prevent condensation in the unused wells. Return pouch immediately to refrigerator.
- All dilutions of the patient samples should be prepared prior to starting with the assay.
- **Good washing technique is critical.** If washing is performed manually, adequate washing is accomplished by directing a forceful stream of wash buffer with a wide tip

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wash bottle across the entire microplate. **An automated microplate washer is recommended.**

- Use a multichannel pipette capable of delivering 8 or 12 wells simultaneously. This speeds the process and provides more uniform incubation times.
- For all steps, careful control of timing is important. The start of all incubation periods begins with the completion of reagent addition.
- Addition of all samples and reagents should be performed at the same rate and in the same sequence.

Test Method

Step 1 Let all reagents and specimens equilibrate to room temperature.

Step 2 Label protocol sheet to indicate sample placement in the wells. It is good laboratory practice to run samples in duplicate.

Step 3 For a **qualitative determination** use Calibrator D (*vial with yellow cap*) only.

or

For a **semi-quantitative determination** use Calibrators A through E as depicted in the sample layout below.

Qualitative				Semi-Quantitative			
A	Blank	S5	Etc.	A	Blank	S1	Etc.
B	-	S6		B	-	S2	
C	+	S7		C	+	S3	
D	Cal D	S8		D	Cal A	S4	
E	S1	S9		E	Cal B	S5	
F	S2	S10		F	Cal C	S6	
G	S3	S11		G	Cal D	S7	
H	S4	S12		H	Cal E	S8	
	1	2	3		1	2	3

Step 4 Prepare a **1:101** dilution of the patient samples by mixing **5 µl** of the patient sera with **500ul** of Serum Diluent.

Step 5 Remove the required microwells from pouch and return unused strips in the sealed pouch to refrigerator. Securely place the microwells into the extra provided holder.

Step 6 Pipette **100 µl** of Ready to use Calibrators, Positive and Negative controls and diluted patient samples (**1:101**) to the appropriate microwells as per protocol sheet.

Note: Include one well which contains **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank.

Step 7 Incubate **30 minutes** (± 5 min) at room temperature.

Step 8 Wash **4x** with wash buffer. For manual washing, fill each microwell with reconstituted wash buffer. Discard the fluid by inverting and tapping out the contents of each well or by aspirating the liquid from each well. To blot at the end of the last wash, invert strips and tap the wells vigorously on absorbent paper towels. For automatic washers, program the washer as per manufacturer's instructions.

Step 9 Pipette **100 µl** of Conjugate into microwells.

Step 10 Incubate **30 minutes** (± 5 min) at room temperature.

Step 11 Wash all microwells as in Step 8.

EN

- Step 12** Pipette **100 µl** of Enzyme Substrate into each microwell in the same order and timing as for the Conjugate.
- Step 13** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 14** Pipette **100 µl** of Stop Solution into each microwell using the same order and timing as for the addition of the Enzyme Substrate. Read absorbance values within 30 minutes of adding Stop Solution.
- Step 15** Read absorbance of each microwell at **450 nm** using a single or at 450/630nm using a dual wavelength microplate reader against the reagent blank set at zero absorbance.

Quality Control

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be <0.3. The Calibrator A should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be <10 U/ml (<10 APL, GPL or MPL). If the test is run in duplicate, the mean of the two readings should be taken for determining units. While performing Qualitative determinations, the optical density of the Calibrator D must be greater than that of the negative control and lesser than the absorbance of the positive control. For semi-quantitative determinations the positive control must give values in the range stated on the vial.

RESULTS

Calculations

The concentrations of the patient samples can be determined by either of two methods:

1. QUALITATIVE DETERMINATION

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator D}} \times \text{Units of Calibrator D} = \text{Units of Test Sample}$$

It is recommended that qualitative results be reported as "positive" or "negative." Sample results greater than or equal to Calibrator D are considered positive.

2. SEMI-QUANTITATIVE DETERMINATION

Plot the absorbance of Calibrator A through E against their respective concentrations on linear-log graph paper. Plot the concentrations in U/ml (APL, GPL, MPL) on the X-axis against the absorbance on the Y-axis and draw a point-to-point curve fit. Determine the concentrations of the patient samples from the curve according to corresponding absorbance values. Alternately, a four parameter curve may be used to plot the standard curve.

It is recommended that semi-quantitative results be reported as "positive," "negative," or "indeterminate" with U/ml unit values. Indeterminate/borderline results should be retested and evaluated along with other laboratory methods, as indicated in **Limitations of Procedure.**

Interpretation

Interpretation values were determined by testing 64 normal blood donors and non-antiphospholipid syndrome disease control specimens. The mean of these subjects plus an optimized SD (IgA=3.5SD, IgG=3SD, IgM=3.5SD) was established as the assay cutoff and assigned a value of 15 U/ml based on titrations of the Harris standards. The following information serves only as a guide in the interpretation of the laboratory results. Each laboratory must determine its own normal values.

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The literature suggests that low positive anti-cardiolipin antibody levels may occur in a variety of clinical disorders unrelated to antiphospholipid antibody syndrome. Hence, according to the investigators recommendations the diagnosis of antiphospholipid antibody syndrome should be made only when the test results are moderately or highly positive.¹⁴ It is suggested that >40 GPL or MPL units be employed to differentiate between low positive ACA levels unrelated to APS and moderate to high ACA levels associated with clinical manifestations of APS.¹⁹

Calibrator

The Ready to Use Calibrators are included to provide semi-quantitation and must be used with each run. Patient samples containing high antibody levels may give absorbance values greater than that of Calibrator A. For determining accurate semi-quantitative values, such specimens should be further diluted so they fall within the range of the calibrator curve when retested. For determining U/ml values, multiply the units obtained by the dilution factor.

LIMITATIONS OF PROCEDURE

The Immulisa™ ACA Test should not be performed on grossly hemolyzed, microbially contaminated or lipemic samples. This method should be used for testing human serum samples only.

Testing for all three isotypes of ACA is strongly recommended. Testing for only one and not all isotypes may lead to false negative results. Furthermore, a diagnosis cannot be made on the basis of ACA results alone. The results of other laboratory tests and clinical findings must be considered. When a negative ACA test occurs in the presence of clinical indications, a lupus anticoagulant test or other additional testing is indicated. ACA also occur transiently in a variety of infectious diseases. In these cases, patients positive for ACA should be retested following an appropriate interval. Confirmed active or seropositive syphilis patients can have elevated ACA levels. To rule out syphilis, confirmatory tests should be performed.²⁰

Anti-cardiolipin antibodies have also been associated with neurological syndromes such as transient ischemic effects and migraine.¹³ In contrast, in patients with antiphospholipid syndrome, the antibodies usually persist for longer periods and may even precede the onset of clinical symptoms.²⁻⁴

EXPECTED VALUES

The incidence of ACA in various disease conditions is summarized in the following tables:

Incidence of ACA in SLE and APS^{15, 21-25}

Antibody Isotype	SLE % Incidence	APS % Incidence
IgG	39-44	36-80
IgA	17-57	10-30
IgM	5-33	15-60
Any isotype	53	71

EN

Disease Association of ACA^{24, 27-31}

Condition	% Incidence
Recurrent Venous Thrombosis	28-71
Recurrent Fetal Loss	28-64
Transverse Myelitis	50
Hemolytic Anemia	38
Thrombocytopenia	27-33
Arterial Occlusions	25-31
Livedo Reticularis	25
Pulmonary Hypertension	20-40

Sets of clinical samples were tested on the Immulisa™ ACA ELISAs. Results demonstrating incidence in the populations for this study are provided below.

Category	IgA			IgG			IgM		
	n	pos	%	n	pos	%	n	pos	%
APS with unknown ACA status	127	11	8.7%	99	56	56.6%	108	38	35.2%
APS with SLE involvement and unknown ACA status	53	30	56.6%	95	47	49.5%	83	19	22.9%
Sera submitted for APS testing	188	14	7.4%	184	45	24.5%	188	13	6.9%
SLE	79	13	16.5%	83	20	24.1%	86	8	9.3%
Sera submitted for SLE testing	90	5	5.6%	88	13	14.8%	90	5	5.6%
Thrombocytopenia	15	1	6.7%	15	2	13.3%	15	0	0%
Pre-eclampsia	15	1	6.7%	15	3	20.0%	15	0	0%
Other AID*	214	6	2.8%	214	10	4.7%	231	15	6.5%
NHS	167	2	1.2%	164	5	3.0%	154	4	2.6%

* Other autoimmune disease includes celiac disease, myositis, rheumatoid arthritis, scleroderma, thyroiditis and vasculitis

PERFORMANCE CHARACTERISTICS

The utility of the Immulisa™ Enhanced ACA ELISAs were evaluated by testing well-characterized specimens from cardiolipin positive SLE subjects alongside disease controls. These specimens were also tested on commercially available ELISA kits. Only specimens in the linear range of the assay were included in the method comparison. These results are summarized below.

A. Immulisa™ Enhanced ACA ELISAs vs. other ACA ELISAs:

Other ACA IgA ELISA				
IMMCO ACA IgA ELISA	Positive	Positive	Negative	Total
	Positive	105	15	120
	Negative	39	462	501
	Total	144	477	621

Positive Percent Agreement: 72.9% (95% CI 64.8% - 79.8%)
 Negative Percent Agreement: 96.9% (95% CI 95.7% - 98.2%)
 Overall Percent Agreement: 91.3% (95% CI 88.7% - 93.3%)

Disease associated subjects (APS, SLE, APS with SLE, suspected APS and SLE): n=432
 Disease controls: n=189

Other ACA IgG ELISA				
IMMCO ACA IgG ELISA	Positive	Positive	Negative	Total
	Positive	155	32	187
	Negative	25	467	492
	Total	180	499	679

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EN

Positive Percent Agreement: 86.1% (95% CI 80.0% - 90.6%)
Negative Percent Agreement: 93.6% (95% CI 91.0% - 95.5%)
Overall Percent Agreement: 91.6% (95% CI 89.2% - 93.5%)

Disease associated subjects (APS, SLE, APS with SLE, suspected APS and SLE): n=464
Disease controls: n=217

Other ACA IgM ELISA				
IMMCO ACA IgM ELISA	Positive	Positive	Negative	Total
		119	2	121
	Negative	50	419	469
	Total	169	421	590

Positive Percent Agreement: 70.4% (95% CI 62.8% - 77.1%)
Negative Percent Agreement: 99.5% (95% CI 98.1% - 99.9%)
Overall Percent Agreement: 91.2% (95% CI 88.5% - 93.3%)

Disease associated subjects (APS, SLE, APS with SLE, suspected APS and SLE): n=408
Disease controls: n=182

B. Cross Reactivity: Sets of potentially cross-reactive specimens from individuals with other autoimmune disorders were tested for ACA using the Immulisa™ Enhanced™ ACA ELISAs.

Condition	IgA n	Pos	% Pos	IgG n	Pos	% Pos	IgM n	Pos	% Pos
Celiac disease	21	0	0.0%	22	0	0.0%	24	0	0.0%
MCTD	24	1	4.2%	23	1	4.3%	29	3	10.3%
Myositis	3	0	0.0%	3	1	33.3%	4	0	0.0%
Rheumatoid arthritis	89	2	2.2%	89	6	6.7%	95	7	7.4%
Sjögren's syndrome	25	2	8.0%	25	2	8.0%	27	3	11.1%
Syphilis	40	6	15.0%	40	33	82.5%	40	14	35.0%
Systemic sclerosis	44	1	2.3%	44	0	0.0%	45	2	4.4%
Thyroiditis	8	0	0.0%	8	0	0.0%	7	0	0.0%

Precision

Precision was tested with 6 specimens selected throughout the range of the assay. Runs of 6 replicates of each specimen were conducted in 13 assays. Repeatability was determined with 12 replicates of each specimen.

Sample	Mean (U/ml)	Total Imprecision		Between days		Within run (Repeatability)	
		SD (U/ml)	CV%	SD (U/ml)	CV%	SD (U/ml)	CV%
IgA 1	2.5	0.4	14.0%	0.366	14.6%	0.234	8.8%
IgA 2	12.4	1.0	8.0%	1.014	8.2%	0.915	7.4%
IgA 3	17.9	1.4	8.0%	1.470	8.2%	1.233	6.9%
IgA 4	29.1	2.4	8.2%	2.499	8.5%	1.466	5.2%
IgA 5	53.1	4.1	7.6%	4.182	7.9%	3.088	5.9%
IgA 6	103.7	6.8	6.5%	6.922	6.7%	5.847	5.5%
IgG 1	3.8	0.183	9.2%	0.177	8.6%	0.207	9.6%
IgG 2	12.0	1.117	9.3%	1.118	9.3%	1.050	9.1%
IgG 3	18.0	1.676	9.3%	1.723	9.5%	1.383	7.7%
IgG 4	31.3	2.738	8.7%	2.805	8.9%	1.651	5.5%
IgG 5	80.4	7.586	9.4%	8.049	10.0%	3.317	4.1%
IgG 6	99.0	6.806	6.9%	7.011	7.0%	4.560	4.7%
IgM 1	11.6	0.930	8.0%	0.930	8.0%	0.407	3.5%
IgM 2	13.0	1.084	8.4%	1.084	8.4%	0.904	7.0%
IgM 3	17.0	1.257	7.4%	1.257	7.4%	0.598	3.5%
IgM 4	45.7	2.292	5.0%	2.292	5.0%	1.305	2.9%
IgM 5	79.9	5.580	7.0%	5.580	7.0%	3.609	4.5%
IgM 6	114.8	6.603	5.8%	6.603	5.8%	4.811	4.2%

Reproducibility

90 replicates of samples in the low negative range, near the cutoff, in the moderate positive range, and approximately +/- 20% of the assay cutoffs for each ACA isotype were performed to determine qualitative reproducibility. Assay results for low negative and moderate positive specimens produced 100% qualitative agreement. +20% specimens produced 98% (IgA), 99% (IgG) and 98% (IgM) qualitative agreement. Approximate cutoff specimens produced 52% (IgA), 60% (IgG) and 60% (IgM) qualitative agreement. -20% specimens produced 99% (IgA), 100% (IgG) and 97% (IgM) qualitative agreement.

Limit of Detection

Based on 60 replicates of the blank and 10 replicates each of 6 low-level (NHS) samples the limits of detection (LoD) for ACA antibodies were determined to be: IgA 1.8 APL, IgG 1 GPL, and IgM 2.8 MPL.

Linearity and Recovery

Studies were performed using equidistant dilution series of positive samples with values throughout the calibrator range to determine linear range of the assay. The linear range of the assay was determined to be: IgA 1.8 to 120 APL, IgG 1 to 120 GPL and IgM 2.8 to 120 MPL. Results are summarized below. It is suggested that samples producing results greater than the top calibrator be diluted and retested for accurate semi-quantitative results.

IgA	Test Range (APL)	Slope (95% CI)	Y-intercept (95% CI)	R2	% Recovery (Obtnd/Expctd)
	2.6 to 47.6	0.999	-0.52	1	97.7% to 114.0%
		(0.95 to 1.05)	(-1.95 to 9.91)		
	18.4 to 70.4	0.997	1.57	0.98	85.1% to 101.9%
		(0.85 to 1.15)	(-5.28 to 8.43)		
	20.0 to 107.6	1.126	0.07	0.99	84.4% to 100%
IgG		(1.00 to 1.26)	(-9.49 to 9.62)		
	Test Range (GPL)	Slope (95% CI)	Y-intercept (95% CI)	R2	% Recovery (Obtnd/Expctd)
	4.5 to 27.8	0.995	0.35	1	90.8% to 101.1%
		(0.94 to 1.05)	(-0.66 to 1.35)		
	2.4 to 98.1	0.989	-0.73	0.99	95.9% to 114.6%
		(0.90 to 1.08)	(-6.11 to 4.66)		
IgM	3.1 to 107.5	1.168	-4.73	0.97	82.6% to 102.8%
		(0.94 to 1.40)	(-20.64 to 11.17)		
	Test Range (MPL)	Slope (95% CI)	Y-intercept (95% CI)	R2	% Recovery (Obtnd/Expctd)
	5.6 to 28.3	0.925	1.70	0.99	83.0% to 102.3%
		(0.83 to 1.02)	(-0.05 to 3.44)		
	5.8 to 86.8	0.97	1.16	0.99	93.4% to 106.9%
		(0.87 to 1.07)	(-4.37 to 6.69)		
	5.5 to 103.6	0.939	4.78	0.99	81.6% to 105.2%
		(0.76 to 1.11)	(-5.89 to 15.46)		

Interference

Interference was studied by mixing sera with known ACA levels for each isotype with potentially interfering serum samples and studying deviation from expected results. No significant interference was demonstrated for the following substances at the levels indicated: Hemoglobin (2 g/L), Bilirubin (342 µmol/L), Rheumatoid Factor (100 EU/ml), Triglycerides (37 mmol/L), and Cholesterol (13 mmol/L).

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