

VASCULITIS • NEPHROLOGY

- 003 ANCA (anti-neutrophil cytoplasmic antibodies)
 014 Circulating Immune Complexes (CIC)
 056 MPO (pANCA antibody)
 057 PR3 (cANCA antibody)
 270 GBM - Glomerular basement membrane antibody titer, IgG, IgA
 290 **Nephrology Profile** includes 003, 270

SPECIMEN REQUIREMENTS

Specimen Collection Kits are available free of charge from Immco Diagnostics. Call (800) 537-8378 for an immediate shipment of collection kits. **Specimens need not be refrigerated or frozen.**

Serum Studies: Collect 5-10 ml of blood in a red top or serum separator tube. If possible, separate serum from clot and place into **clear** tube provided. Do not puncture top of clear tube. If separation facilities are not available, the blood can be sent in the tube used for collection.

Genetic Studies: Send uncoagulated whole blood in 2 **lavender/purple** top tubes. Do not freeze.

Biopsy Studies: When submitting specimens for **immunofluorescence only**, place one biopsy specimen from the edge of the lesion in the **red** capped tube provided with the collection kit and place one biopsy specimen from a normal site in the **purple** tube. These tubes contain a holding solution for immunofluorescence specimens.

When submitting specimens for **routine histological studies (H&E) only**, place one incisional or excisional biopsy specimen in the **green** tube provided with the collection kit or in a tube containing 10% formalin.

When submitting specimens for **immunofluorescence and H&E studies**, take one biopsy specimen from the edge of the lesion and divide it in half. Place one half in the **green** tube provided with the collection kit or in a tube containing 10% formalin for H&E, and place the other half in the **red** tube provided with the collection kit. Take a second biopsy specimen from a normal site and place it in the **purple** tube provided with the collection kit.

For additional information, please consult the Collection & Handling section in the Immco Laboratory Services Directory.

Serum Studies: 2-3 ml serum
 Profile Tests: 2-3 ml serum
 Biopsy Studies: see above

Ship to Immco Diagnostics

Pack Specimen Collection Kit following these simple steps:

- Place tube(s) into specimen bag.
- Place test requisition and a copy of the patient's medical insurance information and any other pertinent paperwork into the outside pocket of specimen bag.
- Fold and roll specimen bag and contents into white mailer box and place in the enclosed poly courier bag.
- Place pre-paid return airbill on outside of poly courier bag and call the number on the bag for a pick up.

Questions? Call 1.800.537.TEST for Immco Customer Service or visit us online at immco.com.

www.immco.com

United States +1.716.691.0091
 Immco Diagnostics, Inc. 800.537.8378
 10 Earhart Drive +1.716.691.6955 Fax
 Suite 100 laboratoryservices@immco.com
 Williamsville, NY 14221 www.immco.com

DERMATOPATHOLOGY • OCULAR PATHOLOGY

- 510 Direct Immunofluorescence (IgG, IgA, IgM, Fibrin, C3, C5b-9, and IgG4)
 511 Light Microscopy (H&E)
 512 Light Microscopic consultation
 513 Hereditary Epidermolysis Bullosa Classification (Collagen IV, Cytokeratin, Collagen VII, laminin).

Special Staining

Differentiation of Bullous Pemphigoid from EBA: Indicated when biopsy is positive for basement membrane zone deposits

- 550 Antigen mapping for localization of Collagen IV & Laminin. Only on Lesional Biopsies.
 551 Localization of immune deposits by "Induced in vitro split." Only on Normal Biopsies.

ORAL PATHOLOGY

- 510 Direct Immunofluorescence (IgG, IgA, IgM, Fibrin, C3, C5b-9, and IgG4)
 511 Light Microscopy (H&E)
 512 Light Microscopic consultation

BIOPSY SITE SELECTION

Proper biopsy sites for various tests and conditions indicated below.

Direct Immunofluorescence

Vesiculo-bullous diseases

Pemphigus, Pemphigoid, Linear IgA Bullous Dermatitis (LABD), EBA, Herpes Gestationis (HG)

Tissue: skin, mucosa

Site: Perilesional, erythematous adjacent to active or new blister. For DH, take skin biopsy 3-5 mm away from the active lesion.

Connective tissue disorders

Lupus erythematosus, MCTD, dermatomyositis, scleroderma.

Tissue: skin

Site: Erythematous or active border of new lesion. Take biopsy for SLE from lesional and non-lesional sun protected skin (buttock).

Lichen Planus and Lichenoid reactions

Tissue: skin, mucosa

Site: take biopsy from a new lesion.

Vasculitis

Tissue: skin

Site: Take biopsy from erythematous, active border of a fresh lesion (<48 hrs old).

Porphyria/Pseudoporphyria

Tissue: skin

Site: take from a new lesion.

Histopathology (H&E)

In all cases, take biopsy from a lesional site.



©2019 Immco Diagnostics, Inc.
 ML 032 REV 5/19

CLINICAL INFORMATION

Please note that lack of complete information may delay results.

Suspected Diagnosis _____

Type of Biopsy Skin Oral Other _____

Biopsy Site: Immunofluorescence (IF)

Lesional _____

Normal _____

Perilesional _____

Biopsy Site: Histopathology (H&E)

Lesional _____

Clinical Data _____



United States +1.716.691.0091
 Immco Diagnostics, Inc. 800.537.8378
 10 Earhart Drive +1.716.691.6955 Fax
 Suite 100 laboratoryservices@immco.com
 Williamsville, NY 14221 www.immco.com

PATIENT INFORMATION

Last Name _____ First Name _____

Social Security Number _____

Sex M F Race _____

Date of Birth _____ **Specimen Date***

Street _____

City _____ State _____ ZIP _____

Telephone (_____) _____ ***Required by State and Federal law.**

REQUESTING DOCTOR

Last Name _____ First Name _____

NPI#* _____ Practice/Facility Name _____

Street _____

City _____ State _____ ZIP _____

Telephone (_____) _____ Fax (_____) _____

e-mail _____

Physician Signature _____

***(Required by CMS)**

BILLING INFORMATION

Send Bill to: Insurance Doctor Lab Patient

For insurance billing, provide complete insurance information or send photocopy of patient's insurance card. If not billing insurance, doctor or lab, please include credit card information to ensure prompt results

Medicare No. if Applicable _____

Name/Address of Insurance _____

Street _____

City _____ State _____ ZIP _____

Telephone (_____) _____ Fax (_____) _____

Insurance ID# _____

Group # _____

Subscriber Name _____

Subscriber DOB _____

ICD-10 Codes _____

Credit Card # _____ Exp. Date _____

Visa MasterCard AmEx

I accept responsibility for bills related to testing of my specimen:

Patient Signature _____

For specimen collection kits please call 800.537.TEST or indicate kits required below:

- # ___ Serum Containers
- # ___ Skin/Mucosal Biopsies and Serum (IF and H&E)
- # ___ H&E Biopsy Studies
- # ___ Otolary (Hearing Loss)

RHEUMATOLOGY

Systemic Lupus Erythematosus (SLE)

- 001 ANA titer and pattern on HEp-2 & Mouse kidney
- 002 ANA titer and pattern on HEp-2 only
- 004 dsDNA (ndNA) antibody titer
- 005 ssDNA antibody
- 007 Histone antibody
- 008 Phospholipid antibody; IgG, IgA & IgM
- 011 Rheumatoid Factor (RF); IgG, IgA & IgM
- 013 RNA antibody
- 014 Circulating Immune Complexes (CIC)
- 055 **β2 Glycoprotein (β2GP1) antibody; IgG & IgM**

Antibodies to Extractable Nuclear (ENA) & Cytoplasmic Antigens

- 040 RNP, Sm, SS-A(Ro), and SS-B(La)
- 042 RNP
- 043 Sm
- 045 SS-A(Ro)
- 046 SS-B(La)
- 047 Scl-70
- 048 Jo-1
- 050 Ku
- 051 PCNA
- 052 PM/Scl

- 070 **CT Profile** includes 001, 004, 007, 011, 040, 047, 048, 050, 051, 052
- 071 **Drug induced LE Profile** includes 001, 004, 007
- 077 **SLE Profile I** includes 001, 004
- 078 **SLE Profile II** includes 001, 004, 040
- 079 **SLE Profile III** includes 001, 004, 007, 011, 040, 047, 048, 050

Scleroderma

- 001 ANA titer and pattern on HEp-2 & mouse kidney
- 005 ssDNA antibody
- 007 Histone antibody
- 042 RNP antibody
- 047 Scl-70 antibody
- 054 Centromere antibody titer
- 062 Centromere IgG ELISA
- 063 RNA Polymerase III
- 075 **Scleroderma Profile** includes 001, 042, 047, 052, 054, 062, 063

Sjögren's Syndrome

- 001 ANA titer and pattern on HEp-2 & Mouse kidney
- 011 Rheumatoid Factor (RF); IgG, IgA & IgM
- 040 RNP, Sm, SS-A(Ro), and SS-B(La)
- 045 SS-A(Ro)
- 046 SS-B(La)
- 088 Jo-1 WB/LB
- 076 **Sjögren's Syndrome Profile** includes 001, 011, 045, 046, 094 (Salivary Protein-1 (SP-1) IgG, IgA, IgM), 095 (Carbonic Anhydrase (CA-6) IgG, IgA, IgM), 096 (Parotid Specific Protein (PSP) IgG, IgA, IgM)
- 093 **Early Sjögren's Syndrome Profile** includes 094 Salivary Protein-1 (SP-1) IgG, IgA, IgM), 095 Carbonic Anhydrase (CA-6) (IgG, IgA, IgM), 096 (Parotid Specific Protein (PSP) (IgG, IgA, IgM)

VESICULO-BULLOUS DISEASES

Pemphigus/Pemphigoid/EBA

- 105 IC and BMZ antibodies (dual substrate)
- 104 IC Paraneoplastic pemphigus antibody titer
- 106 Differentiation of BP from EBA on "split skin"
- 109 Pemphigoid IgG4 antibody
- 112 Desmoglein 1 (DSG1) antibody
- 113 Desmoglein 3 (DSG3) antibody
- 114 Bullous Pemphigoid 180 (BP 180) antibody
- 115 Bullous Pemphigoid 230 (BP 230) antibody
- 122 **Pemphigus/Pemphigoid Profile I** includes 105, 106
- 127 **Pemphigus Profile** includes 105, 112, 113
- 128 **Pemphigoid Profile** includes 105, 106, 114, 115
- 129 **Pemphigus/Pemphigoid Profile II** includes 105, 106, 112, 113, 114, 115

Dermatitis Herpetiformis (DH)

- 100 Endomysial (EMA) IgA antibody titer
- 101 Reticulin (ARA) IgA antibody titer
- 102 Deamidated Gliadin (Celiac G+) antibody; IgA & IgG
- 108 Tissue Transglutaminase (tTG) IgA antibody
- 123 **DH Profile I** includes 100, 108
- 124 **DH Profile II** includes 100, 101, 102, 108
- 125 **DH Profile III** includes 100, 102, 108

Chronic Ulcerative Stomatitis (CUS)

- 107 SES-ANA stratified epithelium specific antinuclear antibody

Dermato/Polymyositis

- 001 ANA titer and pattern on HEp-2 & Mouse kidney
- 045 SS-A(Ro) antibody
- 049 Jo-1 antibody ELISA
- 049 Mi-2 antibody for dermatomyositis
- 050 Ku for polymyositis/scleroderma overlap
- 052 PM/Scl for myositis
- 053 SRP (Signal Recognition Particle)
- 072 **Complete Myositis Profile** includes 001, 048, 049, 050, 053, 081, 082, 083, 084, 085, 086, 087
- 198 **Myositis Interstitial Lung Disease Panel** includes 050, 053, 081, 082, 083, 084, 085, 086, 087, 088, 192, 196
- 199 **Myositis Malignancy Panel I** includes 038, 039, 049, 050, 053, 081, 082, 083, 084, 085, 086, 087, 088, 191, 192, 193, 194, 195, 196, 197
 - 038 U1 snRNP* 088 Jo-1 LIA*
 - 039 U2 snRNP* 191 Anti p155/140*
 - 081 PL7* 192 Anti-p140*
 - 082 PL12* 193 Anti-NXP-2*
 - 083 EJ* 194 Anti-SAE1*
 - 084 OJ* 195 Anti-CN1A*
 - 085 Ro-052* 196 KS*
 - 086 PM/Scl 75* 197 U3 RNP*
 - 087 PM/Scl 100*

* Test only available as part of a panel

Relapsing Polychondritis

- 015 Collagen Type II antibody

Rheumatoid Arthritis

- 001 ANA titer and pattern on HEp-2 & Mouse kidney
- 009 Keratin (AKA) antibody titer
- 011 Rheumatoid Factor (RF); IgG, IgA & IgM
- 017 Cyclic Citrullinated Peptides (CCP) IgG antibody
- 147 HLA DR 4 (Erosive Rheumatoid Arthritis)•
- 149 HLA-B27 Genotyping•
- 073 **RA Profile I** includes 001, 011, 017
- 074 **RA Profile II** includes 001, 009, 011, 017
- 069 **RA Profile III** includes 001, 009, 011, 017, 147•, 149•

Vasculitis/Thrombosis

- 003 ANCA (anti-neutrophil cytoplasmic antibodies) titer
- 056 MPO (pANCA) antibody
- 057 PR3 (cANCA) antibody
- 008 Phospholipid antibody; IgG, IgA & IgM
- 055 **β2 Glycoprotein (β2GP1) antibody; IgG & IgM**
- 058 Phosphatidylserine Antibodies (IgA, IgG, IgM)
- 059 Prothrombin Antibodies (IgG)
- 060 **Antiphospholipid Syndrome (APS) Panel I** includes 008, 055, 058
- 061 **Pre-transplant Thrombotic Risk Screen/APS II** includes 008, 055, 058, 059
- 014 Circulating Immune Complexes (CIC)
- 270 GBM - Glomerular basement membrane antibody titer, IgG, IgA
- 080 **Vasculitis Profile** includes 003, 008, 014, 270

GASTROENTEROLOGY • HEPATOLOGY

Celiac Disease (CD)

- 100 Endomysial (EMA) IgA antibody titer
- 110 Endomysial (EMA) IgG antibody titer
- 101 Reticulin (ARA) IgA antibody titer
- 102 Deamidated Gliadin (Celiac G+) antibody; IgA & IgG
- 108 Tissue Transglutaminase (tTG) IgA antibody
- 111 Tissue Transglutaminase (tTG) IgG antibody
- 150 HLA DQ.α1*0501, β1*0201•
- 120 **CD Profile I** includes 100, 108
- 121 **CD Profile II** includes 100, 101, 102, 108
- 180 **CD Profile III** includes 100, 102, 108
- 130 **Complete CD Profile** includes 100, 102, 108, 150•

Inflammatory Bowel Disease (IBD)

- 003 ANCA (anti-neutrophil cytoplasmic antibodies) titer 214
- Gastric Parietal Cell antibody titer
- 188 **Autoimmune Pancreatitis Panel** includes 001, 011, 400, 407

Autoimmune Liver Diseases

- 001 ANA titer and pattern on HEp-2 & Mouse kidney
- 152 AMA mitochondrial antibody titer
- 153 ASMA smooth muscle antibody titer
- 154 AMA & ASMA antibody titer
- 156 LKM (liver-Kidney microsomal) antibody titer
- 163 SP-100 antibody titer
- 166 SP 100 IgG ELISA
- 167 GP 210 IgG ELISA
- 169 LC-1
- 172 Actin IgA (F-Actin) ELISA
- 185 **PBC Profile** includes 001, 152, 166, 167

Pernicious Anemia

- 214 Gastric Parietal Cell antibody titer

ENDOCRINOLOGY AND INFERTILITY

Type 1 Diabetes (IDDM)

- 215 Islet cell antibody titer
- 223 **Diabetes Profile** includes 143•, 144•, 145•, 215

Thyroid

- 216 Thyroglobulin (Tg) antibody
- 217 Thyroid Peroxidase (TPO) microsomal antibody
- 218 **Thyroid Profile** includes 216, 217

Adrenal

- 213 Adrenal antibody titer

Infertility

- 008 Phospholipid antibody; IgG, IgA & IgM
- 490 Ovary antibody titer
- 492 Spermatozoa antibody titer

Please Note: Test codes listed are categorized by disease association and are not necessarily in numerical order. Tests noted in blue are profiles.

• **Genetic Tests:** 2 tubes uncoagulated whole blood

Immco Diagnostics, Inc. complies with the Health Insurance Portability and Accounting Act (HIPAA). The information provided herein will remain strictly confidential. For a copy of Immco's Policy, please contact customer service.

CARDIOLOGY

- 008 Phospholipid antibody; IgG, IgA & IgM

TRANSPLANT IMMUNOLOGY

- 132 Deceased Donor Evaluation (HLA Typing: A, B Bw, DR, DQ, DRw- Low resolution/ABO/Lymphocyte crossmatch)
- 133 Living Donor for Organ Transplant (HLA antibody Screen & Typing: A, B Bw, DR, DQ, DRw -low resolution/ABO/ Lymphocyte crossmatch)
- 134 Bone marrow Transplant Candidate/Donor (HLA antibody Screen & Typing: A, B Bw, Cw, DR, DQ, DRw - High resolution)
- 135 Donor Specific Antibodies
- 136 Organ Transplant Candidate (HLA antibody Screen & Typing: A, B Bw, DR, DQ, DRw- low resolution/ABO/Pretransplant thrombotic risk assessment)

IMMUNOGENETICS

Disease Associations•

- 143 HLA DQ 2 (Diabetes)
- 144 HLA DQ 8 (Diabetes/Narcolepsy)
- 145 HLA DQ 9 (Diabetes)
- 146 HLA DR 3/DR 4 (AIH, Addison's, Graves Disease)
- 147 HLA DR 4 (Erosive Rheumatoid Arthritis)
- 148 HLA DR 15 (Narcolepsy)
- 149 HLA B27 (Ankylosing Spondylitis)
- 150 HLA DQ a1 *501•, b1 *0201•

OTOLOGY

Autoimmune Hearing Loss (SNHL)

- 001 ANA titer and pattern on HEp-2 & mouse kidney
- 003 ANCA Anti-neutrophil cytoplasmic antibody titer
- 008 Phospholipid antibody; IgG, IgA & IgM
- 011 Rheumatoid Factor (RF); IgG, IgA & IgM
- 014 Circulating Immune Complex (CIC)
- 015 Collagen Type II antibody
- 340 69kD (hsp-70) antibodies
- 350 PO antibodies
- 370 **SNHL Profile I** includes 001, 003, 008, 011, 014, 015, 340, 350
- 375 **SNHL Profile II** includes 340, 350